# SYSTEM APPLICATION

Two part implants



# HEXACONE®

IMPLANT SYSTEM FOR ENDOSSEOUS DENTAL IMPLANTATIONS

**SIMPLADENT®** 

# IMPLANT SYSTEM FOR ENDOSSEOUS DENTAL IMPLANTIONS



Secure anti-rotation through high precision internal hexagon
Apical expanded bone thread
Excellent stability in all bone qualities: double condensation
Universal application for fixed and removable prosthodontics

Abutment alignment and 100% tightness through the taper

**HC2** implants have an apical thread and feature an internal 6-edge, an internal marginal taper and a US standard internal thread.

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).

Due to technical reasons the product dimensions shown in this brochure might deviate from reality.


**HC2** is a registered trademark. **HC2** implants are patent-protected.

In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.

#### **PREPARATORY WORK** FOR TEMPLATE APPLICATION

- 1. Ask your laboratory to prepare a drill template with the determined drill holes for the pilot drills. To be on the safe side, you can ask the laboratory to insert guide sleeves (REF BFH) into the drillholes, which specify the exact drill direction. Please use a 2.0 / 2.2 mm Ø drill for the pilot drilling.
- 2. For the following drill sequences you can use drill stops, which can be attached and tightened to the drill according to the length of drilling channel. Gingival thickness and template height are taken into account as needed. Thanks to the extremely high cutting efficiency of our drills, no ascending drilling sequences will usually be required.


Recommended RPM: 2000-5000. Apply sufficient cooling and allow the cooling to reach the working blades of the drills.



**General note: HC2** implants are used as compression screws. In order to achieve a good bone condensation and implant stability, the drilling should be carried out thinner than the core diameter of the implant. The minimal diameter of the drill depends on the bone density. It is therefore not possible to advise drill-sequences which fit all bone-qualities. Typically in the soft maxillary bone only small drill-diameters are used (e.g. the usage of **DOS1** only, for **HC2** implants with 3.3 - 5.5 mm diameter), whereas in the highly mineralized lower jaw a specific drill sequence with respect to the mineralization of the bone is necessary. For insertion under pressure use the Handgrip. Due to technical reasons **HC2** 2.9 mmd is not available with expanded apical thread. **HC2** implants with diameters 2.9 and 3.3 mm as well as 3.7 mm are not for use as single tooth restauration.

# **SURGERY**

#### 1. Recommended drill sequence



#### 2. Implant packaging



Original packaging




Open the blister using the flap. Remove the label and stick it into the patients record.



The blister (secondary packaging) contains the implant in a sterile tube (primary package).

#### 3. Remove the implant from its packaging

- 1. Open the lid. The implant is connected to the lid through a breakable section.
- 2. Remove the implant without touching the inner walls of the tube.




#### 4. Handling


Attach the insertion tool to the implant by holding the top, to which the implant is secured, with your other hand.

Alternative: Firmly attach the assembled contra-angle handpiece instrument IT 2.5 M to the implant. For ratchets ITL 2.5 can be used as well.

After you have attached the insertion tool, firmly hold the lid in your hand and break the implant off the top along the break joint.

Then insert the implant into the drill hole as much as possible.






#### 5. Insertion

Using the ratchet screw the implant clockwise into the cavity.


The endosseous part of the implant must be **completely** covered by the bone.

After insertion the implant can be turned by a ¼ rotation backwards in order to relieve the bone and allow blood access to the implant site.



# **6.** Remove insertion tool from implant

Remove the insertion tool from the implant.



#### 7. Result

Result: A correctly inserted implant



#### 8. Post-operative treatment

Close the implant with the suitable surgical cover screw.




After healing: Remove the surgical cover screw.



409. Handgrip

Use of Handgrip and HC2 adapter.



Break off the holder.



Insert the implant with axial pressure while turning.



# ${\it Max.}\ insertion\ torque\ for\ diameter$

 2.9 mm
 30 Ncm

 3.3 mm
 40 Ncm

 3.7 mm
 50 Ncm

 4.5/5.5 mm
 60 Ncm

#### 10. Pick Up Impressions


Impression taking with an A-silicone such as Safeprint® by Dr. Ihde Dental. The use of open and closed impression tray is possible.

10.1 Pick-up-procedure with an individual impression tray.

Hex tool HT 1.25

Tightening of the impression post HLT

HC2 Implant



#### 10.2 Prior to the impression


For pick up impressions the tray is inserted over the impression post until the screw peaks out on the other side and becomes accessible for the HEX-tool.

The impression post HLT must not necessarily be unscrewed from the implant in order to remove the impression tray. It can be repositioned later as well

Impression tray

Impression post HLT

HC2 Implant



# **10.3** Taking the impression

Disengage HLT from the implant: HLT remains in the impression


After the impression is taken, the implant is closed with a healing abutment (Gingiva former - straight or anatomic) and the impression is sent to the laboratory.

Loosen screw with HT 1.25

Window in impression tray

HLT

HC2 Implant



**10.4** Preparation of the impression tray for model fabrication

Screw analog against the impression post.

Fasten the laboratory analog in the impression using HT 1.25

HLT

ΙΑ



#### 11. Closed tray impressions

#### 11.1 Impression taking with a closed impression tray


To take impression use an adequately large impression tray.

Impression posts TS/TSL are mounted with the help of the knurled screw

Tighten the impression post with the knurled screw

TS/TSL HC

HC2 Implant




#### 11.2 Inserting impression

The filled impression tray is positioned sufficiently deep over the impression post to also allow an impression of the mucosa.

Impression post TS/TSL HC

HC2 Implant




#### **11.3** Removing impression

When the closed-tray method is applied, the impression post TS/TSL HC remains on the implant after the impression tray is removed. After removal of the impression tray the impression post will be unscrewed and repositioned in the impression.

After the impression is taken, the implant is closed with a healing abutment (Gingiva former - straight or anatomic) and the impression is sent to the laboratory.

TS/TSL HC

HC2 Implant



#### 11.4 Mounting the lab analog

Screw analog IA or IA HC M to the transfer post TS HC. (A)


Afterwards the impression post is repositioned in the impression.  $\ensuremath{\mathbb{B}}$ 

The impression can now be casted. In IA HC M block the lower access to the lock screw out prior to casting.

Tighten the impression post onto the laboratory analog using the knurled screw

TS HC


IA or IA HC M



#### 12. Laboratory procedures

#### 12.1

The impression is poured. Then the impression posts (HLT or TS/TSL HC) are unscrewed from the laboratory analog.



# 12.2

The laboratory analog is now in the proper position and orientation in the Gypsum.





#### 12.3

Positioning of the screwable abutments TLA15 HC, thereby the optimal position and adequate angulation must be determined.

**NOTE** The hexagon must be completely inserted into the analog.



Insert screw

TLA 15 Take care to position the hexagon correctly Arrow: View from above. Turn clockwise.

IA or IA HC M



#### 12.4

Ensure proper position of the abutment when transferring into the mouth.

Tightening torque of the screw during fastening on the implant: 20 Ncm

TLA 15 HC



# 12.5

If more than one angled abutment is used, your laboratory will prepare a detachable synthetic bar (e.g. from Pattern Resin) in order to facilitate the correct positioning in the mouth.

#### TLA 15 HC

Pattern Resin®



#### **HC2** IMPLANTS

HC2 implants have a roughened endosseous surface and a machined apical thread. They feature an internal hex, an internal marginal taper and a US standard internal thread.



#### HC2 WITH AGGRESSIVE APICAL THREAD: HC2

As a result of many years of clinical observation of products, the design of the famous HC2 implant has been revised: the broadened apical thread is fully self-cutting. Thanks to the new apical thread portion, the implant is more stable even in weak bone and higher insertion torque can be reached.

If the implant is anchored in the 2<sup>nd</sup> cortical, it may be used in immediate load protocols. Especially in the upper jaw the usage of the new handgrip (REF 13-311431, with adapter IT HC REF 13-418196,) for inserting the implant is mandatory. This tool allows to apply vertical insertion forces and will enhance the anchorage. The drill sequence remains unchanged compared to the former design of the HC2 implant. And of course all abutments and tools remain the same.

Should the first cortical be unusually firm, the insertion can be achieved using the handgrip REF 13-311431 with the adapter IT HC REF13-418196.

Dimensions HC2 4.5 + 13

a) basal thread  $\emptyset$  4.3 mm

b) height of the apical thread 3.2 mm height of polished collar 2.0 mm

c) Implant neck height 2.5 mm

d) nominal Ø 3.7 mm

# **Application limitations**

HC2/HC implant with a diameter of 2.9 mm - 3.3 mm implants may not be placed in a loaded area, especially not in the molar or premolar area. Likewise these implants may not be used where diagonal loading (off-axis loading) occurs, i.e. not for upper anteriors. Under no circumstances may Hexacone 2.9 mm implants be used for work that involves unsupported occlusal surfaces (consoles). If used in immediate load protocols, the prosthetic construction must be safely inserted on the 2<sup>nd</sup> postoperative day, and it should not be removed within the first 6 months.

In general we recommend to use implants up to (and including) the diameter  $3.7 \, \text{mm}$  with care and not to use them for single tooth replacements, unless strict force control/force distribution is guaranteed.

### **HC2** IMPLANTS WIDE APICAL THREAD

Maximum insertion torque see page 7

	Description	endosseous Ø	endosseous length	REF	Price cat.
	HC2 3.3 8	3.3 mm	8 mm	13-412220	G
	HC2 3.3 10	3.3 mm	10 mm	13-412221	G
=	HC2 3.3 11.5	3.3 mm	11.5 mm	13-412222	G
	HC2 3.3 13	3.3 mm	13 mm	13-412223	G
	HC2 3.3 15	3.3 mm	15 mm	13-412224	G
	HC2 3.7 8	3.7 mm	8 mm	13-412202	G
3	HC2 3.7 10	3.7 mm	10 mm	13-412203	G
基	HC2 3.7 11.5	3.7 mm	11.5 mm	13-412210	G
-	HC2 3.7 13	3.7 mm	13 mm	13-412204	G
	HC2 3.7 15	3.7 mm	15 mm	13-412205	G
	HC2 4.5 8	4.5 mm	8 mm	13-412206	G
	HC2 4.5 10	4.5 mm	10 mm	13-412207	G
-	HC2 4.5 11.5	4.5 mm	11.5 mm	13-412208	G
m	HC2 4.5 13	4.5 mm	13 mm	13-412209	G
	HC2 5.5 8	5.5 mm	8 mm	13-412211	G
	HC2 5.5 10	5.5 mm	10 mm	13-412212	G
BED	HC2 5.5 11.5	5.5 mm	11.5 mm	13-412213	G
	HC2 5.5 13	5.5 mm	13 mm	13-412214	G



#### HC2 6+2 IMPLANTS WITH AGGRESSIVE APICAL THREAD

HC2 6+2 was especially developed for the area of the 1st and 2nd molars in the upper and lower jaw. It is possible and recommendable to use it as a compression screw implant in the upper jaw. Endosseous length 6-8 mm. The upper edge of the polished 75° reverse cone can end at bone level or slightly above it. HC2 6+2 implants have a laser-generated surface structure (no-itis® laser) in the enossal area.

#### The conical polished implant head (c) should be submerged into the bone.

	75°conus	
C		
1	b	
-		
a	—	

a) endosseous Ø:	4.5 - 5.5 mm
b) endosseous length:	6 - 8 mm
c) reverse conus:	2 mm

Description	endosseous Ø	endosseous length	REF	Price cat.
HC2 4.5 6+2	4.5 mm	6-8 mm	13-412217	G
HC2 5.5 6+2	5.5 mm	6-8 mm	13-412218	G

Delivery includes surgical screw CSTI, REF 13-418101.



# TRADITIONAL **HC 6+2** IMPLANTS

**HC** 6+2 was especially developed for the area of the 1st and 2nd molars in the upper and lower jaw. It is possible and recommendable to use it as a compression screw implant in the upper jaw. Endosseous length 6-8 mm (8 mm incl. reverse cone). The upper edge of the polished 75° reverse cone can end at bone level or slightly above it. **HC** 6+2 implants have a laser-generated surface structure (No-Itis® laser) in the enossal area.

# The conical polished implant head (a) should be submerged into the bone.




a) endosseous Ø:	4.5 - 5.5 mm
b) endosseous length:	6 - 8 mm
c) reverse conus:	2 mm

Description	endosseous Ø	endosseous length	REF	Price cat.
HC 4.5 6+2	4.5 mm	6-8 mm	13-413217	G
HC 5.5 6+2	5.5 mm	6-8 mm	13-413218	G



### TRADITIONAL **HC** IMPLANTS

Maximum insertion torque see page 7



a) endosseous nominal Ø 2.9 - 5.5 mm

b) endosseous length 8 - 15 mm

c) length micro thread / 2.1 mm polished neck

1	
3	II.
4	

Description

HC 3.38

endoss. Ø

3.3 mm

endoss. length

8 mm

REF

13-413220

Price cat.

G











HC 4.5 11.5



HC 3.3 10	3.3 mm	10 mm	13-413221	G
HC 3.3 11.5	3.3 mm	11.5 mm	13-413222	G
HC 3.3 13	3.3 mm	13 mm	13-413223	G
HC 3.3 15	3.3 mm	15 mm	13-413224	G
HC 3.78	3.7 mm	8 mm	13-413202	G
HC 3.7 10	3.7 mm	10 mm	13-413203	G
HC 3.7 11.5	3.7 mm	11.5 mm	13-413210	G
HC 3.7 13	3.7 mm	13 mm	13-413204	G
HC 3.7 15	3.7 mm	15 mm	13-413205	G
HC 4.1 8	4.1 mm	8 mm	13-413300	G
HC 4.1 10	4.1 mm	10 mm	13-413301	G
HC 4.1 11.5	4.1 mm	11.5 mm	13-413302	G
HC 4.1 13	4.1 mm	13 mm	13-413303	G
HC 4.5 8	4.5 mm	8 mm	13-413206	G
HC 4.5 10	4.5 mm	10 mm	13-413207	G

11.5 mm



4.5 mm




13-413208

G

#### **HEXACONE® PLUS MU 0° IMPLANTS**

g] d e		Description	max. nominal Ø / without apical thread	max. nominal Ø / with apical thread	endosseous length	REF	Price cat.
c I	Maximum insertion	HC Plus MU 3.3 13 0°	3.3 mm	4 mm	13 mm	13-412250	G
i 🔙 📗	torque: 50 Ncm. Material Ti6Al4V	HC Plus MU 3.3 15 0°	3.3 mm	4 mm	15 mm	13-412251	G
		HC Plus MU 3.3 17 0°	3.3 mm	4 mm	17 mm	13-412252	G
		HC Plus MU 3.3 19 0°	3.3 mm	4 mm	19 mm	13-412253	G
		HC Plus MU 3.3 21 0°	3.3 mm	4 mm	21 mm	13-412254	G
<b>₹</b>	Dimensions	HC Plus MU 3.3 23 0°	3.3 mm	4 mm	23 mm	13-412255	G
a	HC Plus MU 0° 4.1 17	HC Plus MU 4.1 10 0°	4.1 mm	4.7 mm	10 mm	13-412259	G
		HC Plus MU 4.1 13 0°	4.1 mm	4.7 mm	13 mm	13-412260	G
a) endosseous maximal Ø	3.3 / 4.1 mm	HC Plus MU 4.1 15 0°	4.1 mm	4.7 mm	15 mm	13-412261	G
b) endosseous length	11.5 - 21.5 mm	HC Plus MU 4.1 17 0°	4.1 mm	4.7 mm	17 mm	13-412262	G
c) length micro thread	1.5 mm	HC Plus MU 4.1 19 0°	4.1 mm	4.7 mm	19 mm	13-412263	G
d) platform Ø	4.8 mm	HC Plus MU 4.1 21 0°	4.1 mm	4.7 mm	21 mm	13-412264	G
e) height head	2.6 mm	HC Plus MU 4.1 23 0°	4.1 mm	4.7 mm	23 mm	13-412265	G
f) height of the apical threa	d 3.2						
		Description			Code	REF	Price cat.
		Insertion tool incl. screv	w REF 418316. For Hexacon	e Plus MU 0°.	IT HCMU	13-418315	F

#### **HEXACONE® PLUS MU 15° IMPLANTS**





Use with IT2 BCS, IT2 S BCS, AH MU, handgrip. Tool for the screw: HT 1.25


# **ACCESSORIES** FOR HEXACONE® PLUS MU

	Description	Code	REF	Price cat.
OUT	Ratchet for all hex instruments and insertion tools.	RAT2	13-425051	K
	Torque wrench 10 - 70 Ncm. It is recommended to have the torque ratchets recalibrated by us once a year.	TW2	13-425402	S

# ACCESSORIES

Description	Code	REF	Price cat.
Insertion tool for KOS MU & BCS MU Use with IT2 BCS, IT2 S BCS, AH MU	ITX MU15	13-418203	F
Titanium base Use with SF K MU (REF 418164)	T-Base MU	13-418188	В
Castable abutment for use with T-Base and SF K MU	PA2 MU	13-418189	В
Lab analogue for MU-implants	IA K MU	13-418159	A
Prosthetic screw for KOC® MU and BECES® MU	SF K MU	13-418164	В
Long screw for prosthetic use or as pick-up screw for use with TS MU (Tool: HT 1.25). Material Ti6Al4V	SFL MU	13-418168	В
Castable abutment UCLA for direct use on MU-implants. SF K MU sold separately	PA MU	13-418119	В
Transfer Coping (Temporary base) SF K MU must be ordered separately	TC MU	13-418161	D
Transfer for pick-up, straight Delivery incl. SFL MU	HLT MU	13-418162	С
Adapter for handgrip Fits ITX MU15 (REF 13-418203)	AH MU	13-900041	F
Scan abutment for MU implants, incl. screw SSA MU. Sterilisable, two-part Material Ti6Al4V	SAB MU	13-418205	D

#### **SCANBODIES**




Please go to http://simpladent-implant.com/en/stl to download the corresponding STL files.

# **SURGICAL ACCESSORIES**

 $Application \ limitations \ Hexacone \ ^\circ 2.9 \ mm \ implants \ may \ not \ be \ placed \ in \ a \ loaded \ area, especially \ not \ in \ the \ molar \ or \ premolar \ area. \ Likewise \ these \ implants \ may \ not \ be \ used \ where \ diagonal \ loading \ (off-axis \ loading) \ occurs, i.e. \ not \ for \ upper \ anteriors. \ Under \ no \ circumstances \ may \ Hexacone \ 2.9 \ mm \ implants \ be \ used \ for \ work \ that \ involves \ unsupported \ occlusal \ surfaces \ (consoles). \ If \ used \ in \ immediate \ load \ protocols, the \ prosthetic \ construction \ must \ be \ safely \ inserted \ on \ the \ 2^{nd} \ postoperative \ day, \ and \ it \ should \ not \ be \ removed \ within \ the \ first 6 \ months.$ 

In general we recommend to use implants up to (and including) the diameter 3.7 mm with care and avoid using them for single tooth replacements, unless strict force control is guaranteed.



Screwable abutments for cemented bridges, without anti-rotation protection. Trimming and grinding is possible. Tighten with HT 1.25. Recommended insertion torque 20 Ncm.



Description	Code	REF	Price cat.
Height above implantat 8.5 mm The impression is made directly on the TCA, with tool TZ HC	TCA	13-418129	В
The impression is made directly on the TCA	TCA W	13-418173	В

Superstructure with hex and screw. Straight, for cemented bridges, without anti-rotation protection. Trimming and grinding is possible. Tighten with HT 1.25. Delivery inclusive screw SF 20. Recommended insertion torque 20 Ncm.



Description	Code	REF	Price cat.
Abutment, height above implantat 8.5 mm	TLA HC	13-418133	D
Abutment, narrow, for HC 2.9	TLAS	13-418134	D
Abutment with 2 mm gingival height	TLA HC2	13-418170	D
Abutment with 4 mm gingival height	TLA HC4	13-418171	D



Description	Code	REF	Price cat.
Abutment	TLA W	13-418193	D
Anatomical abutment	ANAB	13-418276	E



Description	Code	REF	Price cat.
15° angled, 1 mm gingival height	TLA15 HC1	13-418135	F
15° angled, 2 mm gingival height	TLA15 HC2	13-418136	F
15° angled, 3 mm gingival height	TLA15 HC3	13-418137	F
25° angled, 1 mm gingival height	TLA25 HC1	13-418139	F
25°angled, 2 mm gingival height	TLA25 HC2	13-418140	F
25° angled, 3 mm gingival height	TLA25 HC3	13-418141	F



Delivery inclusive screw SF 20

Ī	

Description	Code	REF	Price cat.
Castable abutment for TLA HC2/4	PA TLA HC	13-418172	Α
Castable abutment	PA U	13-418181	A
For TLA HC and TCA			
Ŧ. (	77.116	42 440470	
Transfer post For TLA HC and TCA	TZ HC	13-418179	A

# **IMPRESSION TAKING AND LABORATORY ACCESSORIES**

	Description	Code	REF	Price cat.
<b>W</b>	Impression post Click-on No screw is needed	HLTC	13-418107	С
<b>T</b>	Short Impression post For TLA, TLA 15 and TLA 25 For Pick-up, with screw	HLT	13-418108	С
	Pickup-Screw For HLT REF 418108	SF HLT long	13-418185	В
	Impression post for HC Height 10.6 mm	TS HC	13-418109	С
	Impression post for HC Height 15.5 mm	TSL HC	13-418110	С
	Long impression post With screw	HLTS	13-418118	С
	Lab analogue for Hex	ІА НС	13-418113	В

# **DIGITAL IMPRESSION TAKING**



Scanbody for digital impression taking	
Screw SF 20 is optional and must be	
ordered separately	

Description

Material	Unit	Code	REF	Price cat
POM	Pack of 5	Scanbody HC	13-418288	В

Screwable spacer abutment for bridges and bars. Screw in with HT 1.77. Recommended insertion torque 25 Ncm.



Description	
for gingival height 3 m	m
for gingival height 4 m	m
for gingival height 5 m	m
for gingival height 6 m	m

Code	REF	Price category
TSA 3	13-418143	В
TSA 4	13-418144	В
TSA 5	13-418145	В
TSA 6	13-418146	В







Description	TSA Analog	Castable abutment 10.5 mm high Pack of 5	Prosthetic screw for PSS on BTS/TSA
Code	BTS	PSS (white)	SF
REF	13-418152		13-418151
Price cat.	В	В	В

 $Screwable\ mesostructure\ for\ bridges\ and\ bars.\ Screw\ in\ with\ HT\ 1.77\ hex\ key.\ Recommended\ insertion\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ assigned\ with\ this\ position\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ assigned\ with\ this\ position\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ assigned\ with\ this\ position\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ assigned\ with\ this\ position\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ assigned\ with\ this\ position\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ position\ torque\ 25\ Ncm.\ The\ position\ of\ the\ TCT\ hex\ is\ position\ torque\ 25\ Ncm.\ The\ position\ posit$ approach.



Description
For gingival height 0.5 mm
For gingival height 1.5 mm
For gingival height 2.5 mm

Code	REF	Price category
TCT HC 0.5	13-418130	В
TCT HC 1.5	13-418131	В
TCT HC 2.5	13-418132	В














Description	Transfer post	long Screw	TCT analog	Castable abutment 12 mm high inside circular Pack of 5	Castable abutment 12 mm high inside edged Pack of 5	Fastening screw
Code	TST	SFL	BTT	PSTR (gray)	PSTA	SF
REF	13-418147	13-420428	13-418100	13-418124	13-418123	13-418151
Price cat.	В	Α	В	В	В	В

### **TCT SET**

 $This set contains all \, necessary \, components \, for \, the \, mesios tructure. \, For \, bridges \, and \, bars. \, Screwable \, (anti-rotation).$ 



Description	Laboratory analog for HC, HC2, with internal Hex	Long transfer post for HC, HC2, anti-rotation	Short transfer post for HC, HC2, anti-rotation	Castable abutment, round, 12 mm high Pack of 5
Code	IA HC	HLTS	HLT	PSTR
REF	13-418113	13-418118	13-418108	13-418124
Price cat.	В	С	С	В

### **ABUTMENTS**

This abutment converts the internal hexagon of the HC2 implant into an external standard-hexagon. The prosthetic screw is screwed through. It tightens the prosthetic distribution of the HC2 implant into an external standard-hexagon. The prosthetic screw is screwed through. It tightens the prosthetic distribution of the HC2 implant into an external standard-hexagon. The prosthetic screw is screwed through. It tightens the prosthetic distribution of the HC2 implant into an external standard-hexagon. The prosthetic screw is screwed through. It tightens the prosthetic distribution of the HC3 implant into an external standard-hexagon. The prosthetic screw is screwed through. It tightens the prosthetic distribution of the HC3 implant into an external standard distribution of the HC3 implant into an external standard distribution of the HC4 implant into anand the abutment at the same time.

	Description	Material	Code	REF	Price category
SF 275					
	Tempbase for HRA HC	PEEK	ТРВ Е	13-418274	С
	HEX reverse abutment incl. screw SF 275 REF <b>13-418275</b>	Ti6Al4V	HRA HC	13-418273	D

# **LOCALICER®**

 $Suitable\ tool:\ HT\ 1.77.\ We\ recommend\ a\ minimum\ of\ six\ implants\ per\ jaw\ and\ the\ use\ of\ a\ single\ denture\ as\ splint\ when\ using\ LOC\ abutments\ and\ KOS\ LOC.$ 



Description	Hight	Code	REF	Price category
Localicer® for Hexacone®	2 mm	LOC HC 2	13-418116	С
Localicer® for Hexacone®	4 mm	LOC HC 4	13-418117	С

# **ACCESSORIES FOR LOCALICER®**



Analog + impression cap Set

 Code
 REF
 Price category

 AA LOC
 13-462337
 C

 NCS
 13-462338
 D



Set with 5 Caps + 1 Housing (EXTERNAL PRODUCT)

**Pull off force**Yellow 600 g, Pink 1.200 g, Clear 1.800 g, Violet 2.700 g.
Black has no retention and is designed for temporary solutions for up to one month.

# **MULTI-UNIT** ABUTMENTS

 $Insertion \ of the angled \ MU2 \ abut ments \ with \ HT\ 1.25; Insertion \ of the straight \ MU2S-Abut ments \ with \ HT\ 1.77. \ Not \ for \ use \ on \ single \ implant \ constructions.$ 



Description	Material	Code	REF	Price cat.
MU2 17 HC, angled, incl. SF 20	Ti6Al4V	MU2 17 HC	13-418281	L
MU2 35 HC, angled, incl. SF 20	Ti6Al4V	MU2 35 HC	13-418282	L
MU2S 0.5 HC, straight	Ti6Al4V	MU2S 0.5 HC	13-418283	G
MU2S 1.5 HC, straight	Ti6Al4V	MU2S 1.5 HC	13-418284	G
MU2S 2.5 HC, straight	Ti6Al4V	MU2S 2.5 HC	13-418285	G
GF MU2 Gingivaformer incl. SF MU2 Height above abutment shoulder 6 mm	Ti6Al4V	GF MU 2	13-418286	С
MU2 Localicer incl. SF MU2 Height above abutment shoulder 6.7 mm	Ti6Al4V	MU 2	13-418287	С
Prosthetic screw for MU2	Ti6Al4V	SF 20	13-420943	Α

# **ACCESSORIES** FOR MULTI-UNIT ABUTMENTS



Description	Material	Code	REF	Price cat.
Temporary base (SF MU2 sold separately)	Ti6Al4V	TC MU2	13-418290	D
Transfer straight incl. screw SFL MU2	Ti6Al4V	TS MU2	13-418291	С
Castable for Multi-Unit, incl. screw TC MU2, for UCLA on the MU2-abutment		PA MU2	13-418292	Α
Screw for TC MU2	Ti6Al4V	SF MU2	13-418293	A
Laboratory analog for Multi-Unit	Ti6Al4V	IA MU2	13-418295	В
Hex tool long		HT 1.25	13-425100	С
Hex tool for all superstructures		HT 1.77	13-425103	С
T-Base MU2 , for use on the Multi-Unit Abutment, Screw: SF MU2		T-Base MU	13-418195	С
PA 2 MU for use on all T-Base		PA2 MU	13-418189	Α

# **ACCESSORIES**



Description

Price cat.

Ball abutment for fitting prostheses Application on TSA 3-6 abutments only Head diameter 2.5 mm

В

Code SB REF

Price cat. B

13-418153





Description	Height above implant 3-6 mm				
Code	TSA 3	TSA 4	TSA 5	TSA 6	BTS
REF	13-41814	3 13-41814	4 13-418145	5 13-418146	13-418152



Description	<b>Hight</b> 0.5 mm	Code TB 0.5	REF 13-418126	Price cat.
Ball abutment head - Ø 2.5 mm Screw in with HT 1.25 For use with NC - caps	0.5 mm	TB 2	13-418127	В
	0.5 mm	TB 4	13-418128	В



Description
Nylon cap NC 2 pieces / pack (EXTERNAL PRODUCT)
Nylon cap R-NC with increased friction strength (for use with worn out Localicer*) 2 pieces / pack (EXTERNAL PRODUCT)
Sleeve for all NC (EXTERNAL PRODUCT)

Pull-off force	Code	REF	Price cat.
ca. 1200 g, transparent	NC	13-465028	A1
ca. 800 g, pink	NC 1	13-465029	A1
ca. 500 g, yellow	NC 2	13-465030	A1
grün, strong	R-NC	13-465034	A1
rosa, medium	R-NC 1	13-465033	A1
orange, soft	R-NC 2	13-465032	A1
	Н	13-465031	В

# **INSERTION TOOLS**



# **TOOLS**

	Description	Туре	Code	REF	Price cat.
	Hex tool 1.25	long, 21 mm	HT 1.25	13-425100	С
	Torx tool 1.25	long, for all screws, 21 mm	TT 1.25	13-425105	С
	Hex tool 1.25	for contra-angle handpiece, 45 mm	HTW 1.25	13-425111	В
-	Hex tool 1.25	short, 14 mm	HTS 1.25	13-425101	С
	Hex tool 1.77	for all superstructures, 19 mm	HT 1.77	13-425103	С
G	Hex tool 1.25 M	for contra-angle handpiece, long, 26.1 mm	HT 1.25 M	13-425112	В
	Hex tool 1.77 M	for contra-angle handpiece, long, 28.6 mm	HT 1.77 M	13-425113	В
	Hex tool	extra long, 45 mm	HTX 1.25	13-425102	С
	Hex tool	for contra-angle, extra long, 45 mm	HTX 1.77	13-425104	С
	Punch	for contra-angle handpiece, 4.9 mm Ø	PUW1	13-425404	С
	Punch	manual, 5.2 mm Ø	PU	13-425406	С
	Standardized probe	Scale 1 mm for X-ray measurements 22 mm	PDG	13-425400	Α
	Drill extension contra-angle	extends by 19 mm	DX2	13-500704	D
	Guide sleeve	for pilot drill, Titanium, 10 mm, 2.2 mm Ø Pack of 5	BFH	13-425401	A
666	X-ray measuring spheres	Surgical steel, 0.5 mm Ø Pack of 5	RM	13-425403	Α

# **GUIDE JACKET**

	Description	Amount	Material	REF	Price cat.
	BFH 2.0 guide jacket for pilot drill 2.0mmd	Pack of 5	Ti6Al4V	13-425410	В
	BFH 2.5 guide jacket for pilot drill 2.5mmd	Pack of 5	Ti6Al4V	13-425411	В
	BFH 3.0 guide jacket for pilot drill 3.0mmd	Pack of 5	Ti6Al4V	13-425412	В
	BFH 3.2 guide jacket for pilot drill 3.2mmd	Pack of 5	Ti6Al4V	13-425413	В
	BFH 3.5 guide jacket for pilot drill 3.5mmd	Pack of 5	Ti6Al4V	13-425414	В
С	a) length	5 mm			
Ip	b) hight of step	0.7 mm			
	c) max. Ø top	3.7 - 5 mm			
е	d) nominal Ø	3 - 4.4 mm			
d	e) $\emptyset$ of drilling in the drill template	2.05 - 3.55 m	m		



Model with residual teeth for the fabrication of a drill guide for creating cavities for fixating the later drill guide for implant cavities.



Drill guide for creating cavities for later fixation of the surgical dill guide.



Surgical drill guide for safe BCS® placement. The drill sleeves are designed for 2.0 mm Twist drills.

#### **TITAN BASE**



# Description Titan base incl. screw

### Type Abutment basefor zirkonium. Anti-rotation (anti-rotation). Material Ti6Al4V

# Code

REF МВ НС 13-418267

Price cat. D

### **HANDGRIP-TRAY**









for all contra angle handpiece instruments, compatible with Handgrip Ratchet for all hex tools and insertion tools

Torque ratchet\* heavy duty

for all Insertion-, hex- and torx tools 10 - 70 Ncm

Code REF Price cat. Adapter Wst 13-310530 C

> RAT2 13-425051 K

TW 2 13-425402 S

Adapter for Handgrip Adapter IT HC 13-418196 C

Handgrip \*\* self locking Handgriff A 13-311431 V



<sup>\*\*</sup> To clean this tool a heatable ultrasonic bath and a thermo disinfector (i.e. Miele TD-Serie) are required. If these devices are not available in the dental office the handle with REF 13-311431 should be purchased instead.



**HANDGRIP-TRAY** empty Size of closed tray: **W** 90 mm / **L** 195 / **H** 45 REF 13-60043 Price category G

3 pieces) and adapters (max. 8 pieces). Plastic, autoclaveable up to 134° C, not suitable for dry heat

### **HEATLESS® DRILLS** FOR IMPLANTS WITH CONICAL CORE

Surgical steel, color-coded, depth-coded and autoclaveable. The drill is marked with laser depth markings. Use between 3,000 and 5,000 rpm with good cooling and intermittent drill technique. Due to the extremely high cutting performance, you can work without pressure. For the implant systems HC2 and Xign\*. Drill types DFN 3.0 - DFN 4.2-4.5.



	Ø working range	max. working depth	total length	color code	Code	REF	Price cat.
	0.1 - 1.5 mm	15 mm	31.7 mm	yellow	BCD 1	13-900240	С
M. mm	0.1 - 1.5 mm	15 mm	42 mm	yellow	BCDX 1	13-900243	С
UMG 2 2 18	2.0 / 3.6 mm	13 mm	30 mm		DFN 2.9 13	13-418102	Е
-	2.0 / 3.6 mm	15 mm	32 mm		DFN 2.9 15	13-418103	E
(1823) 11 to	2.0 mm	17 mm	36.5 mm		DS 2	13-425001	D
1111	2.8 mm	17 mm	36.5 mm		DS 2.8	13-425005	D
M 00L+24	4.5 mm	25 mm	44.5 mm		DSL+ 2.8	13-425015	E
- OK 30	2.7 mm	18 mm	36 mm		DFN 3.0	13-425030	E
	3.0 mm	18 mm	36 mm		DFN 3.4	13-425031	E
	3.4 mm	18 mm	36 mm		DFN 3.7	13-425032	Е
DFN 4:1	3.5 mm	18 mm	36 mm		DFN 4.1	13-425049	E
DFN 4.2-4.5	4.05 mm	18 mm	36 mm		DFN 4.2 - 4.5	13-425033	E
H	4.4 mm	18 mm	36 mm		DFN 5.5	13-425034	Е
U118 1.1	2.7 mm	18 mm	39 mm		DFLN 3.0	13-425035	Е
	3.0 mm	18 mm	39 mm		DFLN 3.4	13-425036	Е
WILLIAM TO THE STATE OF THE STA	3.4 mm	18 mm	39 mm		DFLN 3.7	13-425037	Е
- Druk totalar	4.05 mm	18 mm	39 mm		DFLN 4.2 - 4.5	13-425038	Е
OF DELIVERY OF THE PARTY OF THE	3 mm	25 mm	43.5 mm		DFLN+ 3.4	13-425029	Е
Married Towns	3.4 mm	11.5 mm	30 mm		DFSN 3.7	13-425039	D
	3.9 mm	11.5 mm	30 mm		DFSN 4.2 - 4.5	13-425040	D
	max. 3.8 mm	max. 5 mm	27 mm		C Drill 3.7	13-425043	D
C Drill 4.1	max. 4.1 mm	2.5 mm	27 mm		C Drill 4.1	13-425050	D
OD114245	max. 4.6 mm	max. 5 mm	27 mm		C Drill 4.2 - 4.5	13-425044	D
C Drift 5.5	max. 5.5 mm	2.5 mm	27 mm		C Drill 5.5	13-425045	D

TRAY

Autoclaveable up to 134° C. Not suitable for dry heat sterilizers. Size of closed tray: **W** 175 mm **T** 145 mm **H** 65 mm

 $Please\ read\ our\ detailed\ instructions\ for\ cleaning\ and\ re-sterilization\ of\ surgical\ instruments\ on\ {\bf www.implant.com/en/downloads}$ 



**Tray** with content: REF 13-S60017-K Tray empty: REF 13-60017-K

Description	Code	REF	Description	Code	REF
Twist drill	BCD 1	13-900240	Insertion tool short	IT 2.5	13-418174
Twist drill	DS 2	13-425001	Insertion tool medium	IT 2.5 M	13-418150
Twist drill	DS 2.8	13-425005	Universal adapter	UAW	13-425107
Form drill	DFN 2.9 13	13-418102	Hex tool 1.25 long	HT 1.25	13-425100
Form drill	DFN 2.9 15	13-418103	Hex tool 1.25 short	HTS 1.25	13-425101
Form drill	DFN 3.0	13-425030	Hex tool 1.77	HT 1.77	13-425103
Form drill	DFN 3.7	13-425032	Punch	PUW 1	13-425404
Form drill	DFN 4.2 - 4.5	13-425033	Drill extension	DX 2	13-500704
Form drill	DFN 5.5	13-425034	Standardized probe	PDG	13-425400
Form drill	DFSN 3.7	13-425039	Standardized probe	PDG	13-425400
Form drill	DFSN 4.2 - 4.5	13-425040	Standardized probe	PDG	13-425400
Cortical drill	C Drill 3.7	13-425043	Twist drill	DFLN 3.0	13-425035
Cortical drill	C Drill 4.2 - 4.5	13-425044	Twist drill	DFLN 3.7	13-425037
Cortical drill	C Drill 5.5	13-425045	Twist drill	DFLN 4.2 - 4.5	13-425038
Insertion tool long	ITL 2.5	13-418175	Torque ratchet	TW 2	13-425402

# DRILLSTOP - TRAY

Not suitable for dry heat sterilizers.



Description	Code	REF	Price €
Drillstop A		13-500881	
Drillstop C		13-500883	
Drillstop D		13-500884	
Drillstop E		13-500885	
Drillstop G		13-500887	
Drillstop I		13-500889	
Drillstop J		13-500890	
Drillstop K		13-500891	
Drillstop L		13-500892	
Formdrill	DFN 3.0	13-425030	
Formdrill	DFN 3.4	13-425031	
Formdrill	DFN 3.7	13-425032	
Formdrill	DFN 4.1	13-425049	
Formdrill	DFN 4.2 - 4.5	13-425033	
Formdrill	DFN 5.5	13-425034	
Formdrill	DFLN 3.0	13-425035	
Formdrill	DFLN 3.4	13-425036	
Formdrill	DFLN 3.7	13-425037	
Formdrill	DFLN 4.2 - 4.5	13-425038	
Drillstop Tray with content		13-60031-K	739.00

# **STARTER TRAY**

This surgical kit contains all drills and tools for first works with the HC2 system. Material: Plastic autoclaveable up to 134° C Not suitable for dry heat sterilizers.



Description	Code	REF	Price €
Insertion tool	IT 2.5	13-418174	
Insertion tool	ITL 2.5	13-418175	
Insertion tool	ITM 2.5	13-418176	
Hex tool long	HT 1.25	13-425100	
Twist drill	DS 2.0	13-425001	
Twist drill	DS 2.8	13-425005	
Formdrill	DFN 3.0	13-425030	
Formdrill	DFN 3.4	13-425031	
Formdrill	DFN 3.7	13-425032	
Formdrill	DFN 4.1	13-425049	
Formdrill	DFN 4.2-4.5	13-425033	
Corticalis drill 3.7	C-Drill 3.7	13-425043	
Corticalis drill 4.1	C-Drill 4.1	13-425050	
Corticalis drill 4.2 - 4.5	C-Drill 4.2 - 4.5	13-425044	
Torque ratchet	TW2	13-425402	
Starter Tray for HC2 with content		13-S60021-K	upon request
Starter Tray for HC2 empty		13-60021-K	upon request

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).

 $\label{thm:constraints} \mbox{Due to technical reasons the product dimensions shown in this brochure might deviate from reality.}$ 

Hexacone \$ is a registered trademark.

**Hexacone**\* implants are patent-protected.

In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.

(The products of this catalogue are CE marked (class I) and CE 1936 marked (class IIa and IIb) according to 93/42/EC Directive).

Commercial products that are not monitored by our notified body are declared as third-party products.

Compilation and explanation of symbols on the packaging:



Batch No.



Sterilized by radiation



Non-sterile



Intended for use by dentists or surgeons only



Single use product



Instruction for use



Expiry date



Store in a dry place



Store tightly keep closed



Do not use if packing is damaged



Do not resterilize



Manufacturer



Production date



Catalogue number



Secure anti-rotation through high precision internal hexagon

Apical expanded bone thread

Excellent stability in all bone qualities: double condensation

Universal application for fixed and removable prosthodontics

Abutment alignment and 100% tightness through the taper

Simpladent GmbH

Dorfplatz 11 CH - 8737 Gommiswald / SG SWITZERLAND Phone +41 (0)55 293 23 70 implants@simpladent.ch www.simpladent-implants.com

PLEASE ASK YOUR LOCAL DISTRIBUTOR FOR THE VALID PRICE LIST

SIMPLADENT®