SYSTEM APPLICATION

Single-piece implants



IMMEDIATE LOADING DENTAL IMPLANT SYSTEM

SIMPLADENT®

SIMPLADENT®

CONTENT

Impression taking and laboratory accessories

KOC® instruction for application	5	Scanbodies	33
KOC® Classic implants	13	Instruments and tools	33
KOC® Classic X implants	14	Overview drilling depth and max. working length	34
KOC® B implants with small head for bridges	14	KOC® Mini implants	35
Impression taking and laboratory accessories	15	Surgery	36
KOC® PLUS implants	17	Accessories	37
Example of use of KOC® Plus implants	17	BECES® implants	40
KOC® TX implants	18	Fitting and cementing of prosthetics	41
Insertion tools	18	BECES® implants with small abutment head	42
Impression taking and laboratory accessories	18	BECES® implants with large abutment head	44
Einteilige Multi-Unit implants	19	Use of the handgrip	49
Accessories single-piece Multi-Unit implants	20	BECES® EX implants	49
Application of single-piece Multi-Unit implants	21	STRATEGIC IMPLANT® with small abutment head	50
KDS bone expanding screws	22	STRATEGIC IMPLANT® with large abutment head	51
Insertion tools	23	Accessories for BECES® and KOC® implants	53
Instruments and tools	23	Pathfinder and Twist drills	54
Handgrip	24	Insertion tools and adapter	55
Handgrip tray	24	Intra-oral welding	55
Instrument tray for KOC® and BECES®	25	Instrument tray for KOC® & BECES®	58
Drillstop tray	26	Insertion tools	59
Starter tray	26	Starter tray	59
KOC® implantations	27	Single Piece Implant Pro Kit	60
KOC® M andMicro implants	28	ZDI implants	61
Material / functionality / drilling procedure / insertion	29	Accessories for ZDI implants	61
KOC® M implants with small abutment head	30	ZDI MU implants	62
KOC® Micro implants with large abutment head	30	Accessories for ZDI MU implants	62
KOC® Micro D implants with large abutment head	31	Reprocessing of tools and drills	62

ENDOSSEOUS DENTAL IMPLANT SYSTEM KOC®



KOC® implants are used for crowns, bridges and bars. The compression screw design permits to incorporate the restoration in an immediate loading protocol (incorporation of the prosthesis within max. three days).

KOC® implants are routinely used for more than two decades in immediately loading protocols. The single- piece design reduces costs, the danger of peri-implantitis and it eliminates the hazzles of screw loosening. In extraction cases, **KOC®** and **BECES®** are combined: **BECES®** are placed into the fresh sockets, while **KOC®** are used in healed bone areas.

KOC[®] Classic implants are straight implants with small prosthetic head for crowns, bridges and bars, for cementation.

KOC® - INSTRUCTION FOR APPLICATION


PREPARATORY WORK

Get your lab to make a drilling template with the specified drill holes for the marking

For the pilot hole, use **DOS 1** or **BCD 1** (yellow) as the primary reamer. Prepare the implant bed with the form drills at full length.

Please use an intermittent drilling technique with good NaCl cooling. If necessary, the laboratory can insert guide sleeves can in the drill holes (code **BFH**) through which the precise direction of drilling can be set.


If, due to high drilling resistance in hard bone, it is difficult to reach the complete drilling depth with **DOS 1**, the correct depth can be reached with the cylinder drill **DS 2** (diameter 2 mm).





SURGERY

1. Drilling and preparation/compaction of the implant site

DRILL SEC	UENCE normal/	hard bone		DRILL SEC	UENCE soft bor	ne	
Pilot drill	Form drill	KDS	Implant	Pilot drill	Form drill	KDS	Implant
		KDS 3.0	KOC 3.0				KOC 3.0
	DOS 2	KDS 3.2	KOC 3.2			KDS 3.0	KOC 3.2
DOS 1	DOC 2 (4)	KDS 3.7	KOC 3.7	DOS 1	DOC 3	KDS 3.2	KOC 3.7
	DOS 3 (4)	KDS 4.1	KOC 4.1		DOS 2	KDS 3.7	KOC 4.1
	DOS 5	KDS 5.0	KOC 5.0		DOS 3 (4)	KDS 4.1	KOC 5.0







Direction and depth calculation; alternatively BCD 1 "Pathfinder" drill.

Pilot drill DS 2 For use in hard bone in the cortical region only.

KDS Prepare the implant bed in the maxilla stepwise using the appropriate bone-expanding screw and ratchet or motorized insertion tool. Maximum 40-45 Ncm. Remove the bone-expanding screw again.

KOC® B To create the definitive implant cavity for KOC® B implants, it is imperative to use bone-expanding screws. These screws must be screwed to their full depth. They generate the compression and ensure that sufficient space is created for the implant thread in the cortical region.

All **KOC*** implants are used as compression screws. If possible, the hole should be created substantially thinner than the core diameter of the implant, since only in this way can good bone condensation be achieved. The minimum hole diameter depends on the bone density. For this reason, it is not possible to specify drill sequences that can be used favorably for all bone qualities. As a rule, it is necessary to drill much less into the soft maxilla (e.g. the DOS1 drill only can be used for **KOC*** implants with diameter 3.0-5.0) than into the well-mineralized mandible, which requires the use of a drilling sequence adjusted to the bone density.

2. Implant packaging



Original packaging



Open the sealed cover at the lid. Remove the label and place it into the patients record.

3. Remove the implant from its packaging



The open pack contains the implant, mounted to a plastic holder.

The pack also contains the lab-set.



Remove the implant by holding onto the plastic holder

The implant is fixed to the holder by a break joint.

4. Handling

Hold the implant by the holder and place the insertion tool on the implant head. The endosseous implant surface must not be touched. Pull out the implant with the plug and then twist off the plug with the needle holder at the predetermined breaking point.

IMPLANTS WITH LARGE HEAD

KOC® / KOC® Micro

KOC® K (for ball attachment)



KOC® implant with insertion tool IT2W (for angle piece) and IT2 BECES (manual)

KOC® K implant with insertion tool IT TB K

Twisting off the bracket with the needle holder

IMPLANTS WITH SMALL HEAD

KOC® (straight) / **KOC**® **B** (flexible)




KOC* / KOC* B implants with insertion tool ITW K (for angle piece) and IT K (manual)

Twisting off the bracket with the needle holder


5. Insertion using manual tools

Insert the implant by hand until it is firmly seated in the jaw.

IMPLANTS WITH LARGE HEAD



IMPLANTS WITH SMALL HEAD



6. Definitive implant insertion

Using the ratchet, torque ratchet or contra-angle, screw the implant clockwise into the cavity. With **KOC® B**, the use of the torque ratchet is mandatory. The endosseous (blasted) part of the implant must be completely covered by bone. The polished implant neck is located in the mucosa. We recommend screwing the implant into the bone up to 1 mm deeper into the implant neck.

IMPLANTS WITH LARGE HEAD



The head of the bendable KOC® 3.0 & 3.2 and KOC® Micro (all diameters) screws can be bent into the desired position after insertion with the aid of the mounted insertion tool and ratchet.

Maximum bend: approx. 15°. Only one bending operation may be performed. In the maxilla, the motorised insertion tool should be used due to its better implant guidance during insertion.



IMPLANTS WITH SMALL HEAD



IMPORTANT NOTE

KOC B° implants have a predetermined breaking point integrated into the head. If the preparation with bone-expanding screws was not performed sufficiently, high screwing forces can cause the upper head portion to be torn off.


So that the implant can be screwed out again, an additional square is milled below the breaking point, into which the emergency tool **Tool E** can be inserted. The **Tool E** instrument may only be used to remove the implant.

7. Removing the placement aid from the implant

IMPLANTS WITH LARGE HEAD



IMPLANTS WITH SMALL HEAD



8. Result

All implant heads (except for the **KOC® K**) can be reshaped by grinding. The implants can be prosthetically supplied immediately if indicated. The definitive superstructure should be cemented within a few days. Immediate prosthetic splinting with a provisional bridge is recommended.

IMPLANTS WITH LARGE HEAD



IMPLANTS WITH SMALL HEAD



9. Impression

IMPLANTS WITH LARGE HEAD

Bridges



Attachment of the impression post TSPA~5, internally round, for **KOC®**


Individual crowns




Attachment of the impression post TSPA 5, with anti-rotation protection, for **KOC® Micro**




Pressureless impression taking e.g. with Safeprint®



Removal of the individual scoop from the implant post. The impression post is located in the impression material. The impression can be sent to the laboratory.



Pressureless impression taking e.g. with **Safeprint**®



 $Removal\ of\ the\ individual\ scoop\ from\ the\ implant\ post.$ The impression post is located in the impression material. The impression can be sent to the laboratory.


IMPLANTS WITH SMALL HEAD

Bridges / Individual crowns



Attachment of the impression post **TSPA 4**, Internally round, for **KOC®** and KOC® T

Individual crowns




Attachment of the impression post $\mbox{\bf TSKPA 4,}$ with anti-rotation protection, for KOC® and **KOC® T**

crown


LABORATORY PROCEDURES

Attachment of the impression post onto lab analogues

IMPLANTS WITH LARGE HEAD



The modeling is performed on the castable parts $\bf PO4$ (internally round; for bridges and bars) or **PO4A** (edged inside; for individual crowns).




TSPA 5 on IAK M/13-462111 Pull impression from the model. The impression post and analogue are now separated again.



IAK with nylon cap (pink) and sleeve Polymerization of the sleeve H into the prosthesis. Press NC/NC1/NC2 into the sleeve. For initial restorations, NC1 or NC2 should be used.

Fill in plaster

IMPLANTS WITH SMALL HEAD




The modeling is performed on the castable parts $\bf PO4$ (internally round; for bridges and bars) or **PO4A** (edged inside; for individual crowns).

Polymerization of the sleeve H into the prosthesis. Press NC/NC1/NC2 into the sleeve. For initial restorations, NC1 or NC2 should be used.

KOC® IMPLANTS

KOC® CLASSIC IMPLANTS

Straight implants with small head for crowns and bars.



a) endosseus length

d) abutment length e) square SW f) neck \emptyset

b) endosseus \emptyset

c) neck length

g) abutment Ø

Description
KOC 3.0 10
KOC 3.0 12
KOC 3.0 15
KOC 3.2 12
KOC 3.2 15
KOC 3.7 6
KOC 3.7 8
KOC 3.7 10
KOC 3.7 12
KOC 3.7 15
KOC 4.1 8
KOC 4.1 10
KOC 4.1 12
KOC 4.1 15
KOC 4.1 17
KOC 4.1 19
KOC 5.0 10
KOC 5.0 12
KOC 5.0 15

6 - 19 mm
3.0 - 5.0
3.5 mm
6.8 mm
1.9 mm
2.0 / 2.5 / 2.8 mm
3.35 mm

Endosseous Ø	Length	Neck Ø	REF	Price cat.
3.0 mm	10 mm	2 mm	13-455108	F
3.0 mm	12 mm	2 mm	13-455109	F
3.0 mm	15 mm	2 mm	13-455110	F
3.2 mm	12 mm	2 mm	13-455111	F
3.2 mm	15 mm	2 mm	13-455112	F
3.7 mm	6 mm	2.5 mm	13-455106	F
3.7 mm	8 mm	2.5 mm	13-455107	F
3.7 mm	10 mm	2.5 mm	13-455114	F
3.7 mm	12 mm	2.5 mm	13-455115	F
3.7 mm	15 mm	2.5 mm	13-455120	F
4.1 mm	8 mm	2.8 mm	13-455129	F
4.1 mm	10 mm	2.8 mm	13-455130	F
4.1 mm	12 mm	2.8 mm	13-455132	F
4.1 mm	15 mm	2.8 mm	13-455135	F
4.1 mm	17 mm	2.8 mm	13-455136	F
4.1 mm	19 mm	2.8 mm	13-455137	F
5.0 mm	10 mm	2.8 mm	13-455171	F
5.0 mm	12 mm	2.8 mm	13-455172	F
5.0 mm	15 mm	2.8 mm	13-455173	F

KOC® implants are delivered incl. lab-set consisting of REF 13-462111, 13-462029 and 13-462088



KOC 3.0 - 3.2 Max. insertion torque 50 Ncm **KOC 3.7 - 5.0** Max. insertion torque 80 Ncm

KOC® CLASSIC X IMPLANTS


 $\textbf{KOC} \\ \textbf{Classic X} \\ \textbf{implants are prosthetically compatible with \textbf{BCES}} \\ \textbf{implants with diameters 3.6 mm, 4.6 mm and } \\ \textbf{> 5.5 mm, and also with \textbf{BOI}} \\ \textbf{implants.} \\ \textbf{> 1.5 mm, and also with \textbf{BOI}} \\ \textbf{> 1.5 mm, and also with \textbf{AOI}} \\ \textbf{> 1.5 mm, and also with \textbf$

(Description	endosseous Ø	Length	Neck Ø	REF	Price cat.
	KOC X 3.7 10	3.7 mm	10 mm	2.5 mm	13-455720	F
d	KOC X 3.7 12	3.7 mm	12 mm	2.5 mm	13-455721	F
- f	KOC X 3.7 15	3.7 mm	15 mm	2.5 mm	13-455722	F
c P	KOC X 4.1 8	4.1 mm	8 mm	2.8 mm	13-455730	F
1	KOC X 4.1 10	4.1 mm	10 mm	2.8 mm	13-455731	F
	KOC X 4.1 12	4.1 mm	12 mm	2.8 mm	13-455732	F
	KOC X 4.1 15	4.1 mm	15 mm	2.8 mm	13-455733	F
a S	KOC X 4.1 19	4.1 mm	19 mm	2.8 mm	13-455735	F
	KOC X 5.0 10	5.0 mm	10 mm	2.8 mm	13-455740	F
	KOC X 5.0 12	5.0 mm	12 mm	2.8 mm	13-455741	F
-	KOC X 5.0 15	5.0 mm	15 mm	2.8 mm	13-455742	F
b a) endosseous length	8 - 19 mm			KOC® X imp	o-set	
b) endosseous Ø	3.7 - 5.0 mm			sting of REF 13-46 162136 and 13-46		
c) neck length	3.0 mm				ij	
d) length abutment	7.2 mm	ACCESSORIES				
e) neck Ø	2.5 / 2.8 mm	Analogue IAB	REF 13-462106	lmpr	ession post TSPA 5	REF 13-462030
f) Ø abutment	3.9 mm	Pack of 5	Price cat. B	Pack	of 5	Price cat. B
				Red rour		d red analogue are

KOC® B IMPLANTS WITH SMALL HEAD FOR BRIDGES

KOC® B implants with bendable neck (use after pre-drilling and preparation with the bone-expanding screw). Suitable for bridges in the reduced-load range (no individual tooth restorations). The bendable implant now offers dual safety:

- $1. \ \ Torque\ reduction\ by\ pre-compression\ with\ the\ congruent\ bone-expanding\ screw$
- 2. Safety head with predetermined breaking point and dual square (Pat. Pend.)



Description	Code KDS	endosseous Ø	Length	REF	Price cat.
KOC B 3.2 12	D	3.2 mm	12 mm	13-455162	F
KOC B 3.2 15	E	3.2 mm	15 mm	13-455161	F
KOC B 3.7 12	F	3.7 mm	12 mm	13-455164	F
KOC B 3.7 15	G	3.7 mm	15 mm	13-455165	F
KOC B 4.1 15	L	4.1 mm	15 mm	13-455166	F
KOC B 4.1 17	М	4.1 mm	17 mm	13-455167	F

The predetermined fracture site integrated in the abutment prevents the twisting off of the abutment head from the endosseous implant part. The implant socket has to however always be pre-compressed using the bone-expanding screw. **Max. insertion torque 45 Ncm.**

a) endosseous length	12 - 17 mm
b) endosseous Ø	3.2 / 4.1 mm
c) length abutment	6.8 mm
d) neck length	3.0 mm
e) neck Ø	1.8 mm





IMPRESSION TAKING AND LABORATORY ACCESSORIES

	Description	Unit	Material	Code	REF	Price cat.
	Impression post castable, Internally round	Pack of 5		TSPA 4	13-462029	В
	Impression post castable, internal edges	Pack of 5		TSKPA 4	13-462028	В
	Double analogue	Pack of 5	plastic	IA4/IAU	13-462111	В
or	Double analogue	Single	metal	IA4/IAU	13-462112	A
	Castable abutment and base for provisionals. 7 mm high, white, Internally round			PO4	13-462088	В
	Castable abutment and base for provisionals. 7 mm high, white, internally edged			PO4A	13-462089	В

***TSPA 4 & 5** For impressions on ground-down implant heads.

The ring transfer constitutes the lower limit of the head, so that the impression can be poured with die stone or epoxy. The use of an implant analogue together with transfer 4 is not possible with ground-down implant heads. Material: PP

SCAN ANALOGUE

Scanner analogue for large and small implant head, self-descriptive. These analogues do not need to be sprayed with spray paint. They can be pulled out of the model with anti-rotation protection. Matching impression posts: **TSPA 4** and **TSPA 5**.



escription	REF	Price cat
can analogue IAS 4	13-462019	В
can analogue IAS 5	13-462020	В



Use example for self-descriptive scanner analogue.


CEMENTABLE ANGULATION ADAPTER (TIGAL4V)

These adapters are mounted on **KOC*** implants to compensate for the insertion direction. Plastic cements are preferably used. The implant head must be roughened beforehand. The protruding head parts are then removed. The impression is taken directly on the adapter.



Price cat.

Α



CASTABLE CROWN BASE

These adapters are used by the dental technician for modeling of bridge frames. In the metal try-in, the protruding head parts are removed by the dentist.

	Description	Height	Code	REF	Price cat.
8	Adapter Reducible Pack of 5	7.5	AAL 15 KK	13-462045	E

LAB ANALOGUE




CASTABLE PART AND IMPRESSION CAP



KOC® PLUS IMPLANTS

KOC® PLUS are made in one piece and have a polished apical cutting thread for anchorage in the opposite cortical. **KOC® Plus** combines the advantages of the compression screw with the advantages of the bicortical screw: in addition to bone compression, there is also anchorage in the 2nd cortical (opposite cortical). Can be used in the maxilla and mandible. Made of titanium alloy Ti6Al4V, laserd. Tighten with **IT2 BECES**.



EXAMPLE OF USE OF KOC® PLUS IMPLANTS



Bicortical anchorage of a **KOC® Plus** implant in the atrophied distal mandible.



Bicortical anchorage of a **KOC® Plus** implant in the area of the nasal floor.



Bicortical anchorage of **KOC® Plus** implants (Ø 3.7 and 4.1) in the area of the maxillary sinus.

NOTE - KOC® Plus may only be operated/used by validly authorized users. Only polished implant parts may penetrate a maximum of 1.5 mm into the opposite cortical. For a given indication (min. three stable implants, sufficient bone quality, etc.). **KOC® Plus** are suitable for immediate loading.

KOC® TX IMPLANTS

KOC® TX is a compression screw implant with an extended polished and bendable neck (6 mm) for anchorages in the tuberosity and in regions with large mucosal thickness. The roughened thread parts must be completely submerged in the bone. **KOC® TX** is preferably used without flap preparation. The nominal diameter is reached at the widest part of the compression thread.



Description	Max. endosseous Ø	REF	Price cat.	Description	Max. endosseous Ø	REF	Price cat.
KOC TX 3.0 12	3 mm	13-455001	F	KOC TX 3.7 18	3.7 mm	13-455007	F
KOC TX 3.0 15	3 mm	13-455002	F	KOC TX 4.0 12	4 mm	13-455175	F
KOC TX 3.0 18	3 mm	13-455003	F	KOC TX 4.0 15	4 mm	13-455176	F
KOC TX 3.7 12	3.7 mm	13-455005	F	KOC TX 4.0 18	4 mm	13-455177	F
KOC TX 3.7 15	3.7 mm	13-455006	F	KOC TX 4.0 21	4 mm	13-455178	F
Max. insertion to	orque 80 Ncm.						
a) endosseous ler	igth	12 - 21 mm			C® TX implants		
b) max. endosseo	us Ø	4 mm		are deriver	ed incl. lab-set of REF 13-462111,		//
c) abutment leng	th	7.2 mm		-	and 13-462086	8	
d) neckØ		2 mm					
e) neck height		6 mm					
ACCESSORIES							
Analogue IAB	REF 13-462106		Impression pos	t TSPA 5	REF 13-462030		
Pack of 5	Price cat. B		Pack of 5		Price cat. B		
		2	Red impressio	n caps and red an	alogue are round.		

INSERTION TOOLS



IMPRESSION TAKING AND LABORATORY ACCESSORIES



SINGLE-PIECE MULTI-UNIT IMPLANTS

KOC® MU implants feature a pre-angulation of 15 degrees. **KOC® MU** may be bent additionally, using the insertion tool. In conjunction with the clinically possible rotational positions of the head, virtually all possible angulations can be realized. Material **Ti6Al4V**.



BECES® MU implants feature a pre-angulation of 15 degrees. **BECES® MU** may be bent additionally, using the insertion tool. In conjunction with the clinically possible rotational positions of the head, virtually all possible angulations can be realized. **BECES® MU** implants may be used by authorized users only. Material **Ti6Al4V**.

C		J/g	Description	REF	Price cat.	Description	REF	Price cat.
			BECES MU 3.6 8	13-900397	N	BECES MU 4.6 23	13-900385	N
Ť	†	e	BECES MU 3.6 10	13-900398	N	BECES MU 4.6 26	13-900386	N
		d d	BECES MU 3.6 12	13-900376	N	BECES MU 4.6 29	13-900387	N
			BECES MU 3.6 14	13-900330	N	BECES MU 4.6 32	13-900388	N
a			BECES MU 3.6 17	13-900331	N	BECES MU 4.6 35	13-900389	N
		THE CONTRACTOR OF THE CONTRACT	BECES MU 3.6 20	13-900332	N	BECES MU 5.5 10	13-900334	N
	•	TI.	BECES MU 3.6 23	13-900333	N	BECES MU 5.5 12	13-900335	N
	•		BECES MU 3.6 26	13-900377	N	BECES MU 5.5 14	13-900336	N
	1 <u>—</u>	b	BECES MU 3.6 29	13-900378	N	BECES MU 5.5 17	13-900357	N
			BECES MU 3.6 32	13-900399	N	BECES MU 5.5 20	13-900358	N
a) endosseous length		8 - 38 mm	BECES MU 3.6 35	13-900339	N	BECES MU 5.5 23	13-900341	N
b) endosseous \emptyset		3.6 - 7.0 mm	BECES MU 3.6 38	13-900340	N	BECES MU 5.5 26	13-900342	N
c) hight abutment		3.7 mm	BECES MU 4.6 8	13-900379	N	BECES MU 7.0 10	13-900337	N
d) shaft Ø		2 mm	BECES MU 4.6 10	13-900380	N	BECES MU 7.0 12	13-900338	N
e) platform Ø		4.8 mm	BECES MU 4.6 12	13-900381	N	BECES MU 7.0 14	13-900360	N
f) neck height		0.8 mm	BECES MU 4.6 14	13-900382	N	BECES MU 7.0 17	13-900361	N
g) height of connecting pa	irt	2 mm	BECES MU 4.6 17	13-900383	N	BECES MU 7.0 20	13-900362	N
Prosthetic screw		SFK MU	BECES MU 4.6 20	13-900384	N			

MULTI-UNIT LAB SET



Description	Code	REF	Price cat.
Titanium base Use with SF K MU	T-Base MU	13-418188	
Castable abutment Use with T-Base and SF K MU	PA2 MU	13-418189	
Prosthetic screw For KOC® MU and BECES® MU	SF K MU	13-418164	
COMPLETE SET		13-418289	D

ACCESSORIES SINGLE-PIECE MULTI-UNIT IMPLANTS


	Description		Code	REF	
	Insertion tool for KOC® MU, BECES® MU & Hexacone Plus MU 15°. Use with IT2 BECES, IT2 S BECES, AH MU. Tool: HT 1.25		ITX MU15	13-418203	
□ ○	Insertion tool medium, for large head. Use with RAT2 and TW2. Length 19 mm		IT2 BECES	13-900030	
(<u>a</u> 0	Insertion tool short, for large head. Use with RAT2 and TW2. Length 7 mm		IT2 S BECES	13-900038	
0	Adapter for handgrip Fits ITX MU15 (REF 13-418203)		AH MU	13-900041	
	Scan abutment for MU implants, incl. screw SSA MI Sterilisable, two-part. Material Ti6AI4V	U.	SAB MU	13-418205	
Î	Prosthetic screw for KOC® MU and BECES® MU		SF K MU	13-418164	
Parts for passive connection of the bridge frame	Castable abutment for use with T-Base and SF K MI	PA2 MU	13-418189	ŗ	
A	Titanium base * Use with SF K MU (REF 13-418164), for Hexacone Plus MU, BECES® MU, KOC® MU		T-Base MU	13-418188	
•	Prosthetic screw for KOC® MU and BECES® MU		SF K MU	13-418164) <u>.</u>
Parts for UCLA technique	or UCLA technique Castable abutment UCLA for direct use on MU-implants. SF K MU sold separately			13-418119	
Part for UCLA technique & passive connection	Digital lab analogue for MU-implants* For BECES® MU, KOC® MU, Hexacone® MU		IA K MU	13-418159	
	Long screw for prosthetic use or as pick-up screw for use with HLT MU (Tool: HT 1.25). Material Ti6Al4V		SFL MU	13-418168	
	Transfer for pick-up impressions, straight. Delivery incl. SFL MU	Works with all MU implants	HLT MU	13-418162	
	Temporary base SF K MU or SFL MU must be ordered separately		TC MU	13-418161	
	Description		Code	REF)
=	Hex-instrument 1.25, length 14 mm	short	HTS 1.25	13-425101	
	Hex-instrument 1.25, length 21 mm	medium	HT 1.25	13-425100	
-	Hex-instrument, length 45 mm	long	HTX 1.25	13-425102	

APPLICATION OF SINGLE-PIECE MULTI-UNIT IMPLANTS

1.

Tighten screw SFL MU with the tool HT 1.25.

Fix the transfer with the long screw, then take pick-up-impression.



4.

T-Base is sandblasted from the outside and cleaned.

The bridge frame is sandblasted from below in the area of the implants.



2.

Connect the transfer to the implant analogue (IA K MU) and pour the impression with gypsum.



5.

All T-Base are fixed to the implants with SF K MU or the long screw SFL MU. Then all T-Base are glued with adhesive cement to the bridge frame

This guarantees a passive fit. Composite excess is removed and the site is polished.



3. a

Connect PA MU with SF K MU on the analogue IA K MU. Tighten screw SFL MU with the tool HT 1.25.

Now the modulation can be created and the frame is veneered. Veneering is possible with acryl, composite and ceramics.



6.

Now the bridge may be screwed on passive with SF K MU.

Screw canals are closed with temporary filling material or composite, taking into consideration that later access must be possible.



3. b

T-Base is positioned over the analogue and screwed on with SF KMU. The cartable PA2 MU is then fitted on top of the T-Base.

Now the modulation is made. Veneering is possible with acryl, composite and ceramics.



Application of insertion tool MU

Example for insertion tool ITX MU15 on the implant BECES® MU / KOC® MU.



KDS BONE EXPANDING SCREWS

For all **KOC® B** screw implants, bone-expanding screws are available as tools to create the definitive implant cavity. Basically, for each implant prior to insertion of a **KOC® B** screw implant, a bone compression with the bone-expanding screw should be performed. In addition, with a narrow alveolar ridge, an expansion of the alveolar ridge can be performed with the bone-expanding screw. By inserting the bone-expanding screw, it can be checked whether the **KOC® B** screw implant can be inserted into the bone easily and fully. Titanium alloy Ti6Al4V Eli, machined. Tighten with **IT K**, **ITS K** or **ITX K** using the torque ratchet **TW2** (max. 45 Ncm), or alternatively **RAT2**.

Package unit: 1 piece, non-sterile.




Description	Code KDS	Endosseous Ø	Length	Neck Ø	REF	Price cat.
KDS 3.0 10	Α	3.0 mm	10 mm	2.0 mm	13-455212	F
KDS 3.0 12	В	3.0 mm	12 mm	2.0 mm	13-455213	F
KDS 3.0 15	С	3.0 mm	15 mm	2.0 mm	13-455214	F
KDS 3.2 12	D	3.2 mm	12 mm	2.5 mm	13-455223	F
KDS 3.2 15	E	3.2 mm	15 mm	2.5 mm	13-455224	F
KDS 3.7 12	F	3.7 mm	12 mm	2.8 mm	13-455233	F
KDS 3.7 15	G	3.7 mm	15 mm	2.8 mm	13-455234	F
KDS 4.1 8	Н	4.1 mm	8 mm	2.8 mm	13-455241	F
KDS 4.1 10	1	4.1 mm	10 mm	2.8 mm	13-455242	F
KDS 4.1 12	K	4.1 mm	12 mm	2.8 mm	13-455243	F
KDS 4.1 15	L	4.1 mm	15 mm	2.8 mm	13-455244	F
KDS 4.1 17	М	4.1 mm	17 mm	2.8 mm	13-455245	F
KDS 4.1 19	N	4.1 mm	19 mm	2.8 mm	13-455246	F

a) endosseous length b) endosseous Ø c) length (non-blasted) d) length abutment e) neck Ø f) neck length 8-19 mm 3.0/3.2/3.7/4.1/5 mm 10 mm 6.8 mm 2.0-2.8 mm 3 mm The bone-expanding screws can easily be screwed in using suitable insertion tools and immediately screwed out again after reaching the full insertion depth. Subsequently, the KOC® B implant is inserted. With the KOC® B (bendable), the use of bone-expanding screws is mandatory regardless of the region, so that the shear forces occurring during insertion do not fracture the implant neck.

Do not use for KOC® implants with microthread.

INSERTION TOOLS



^{*} Only for W+H contra-angles with new drive

INSTRUMENTS AND **TOOLS**



^{*} It is recommended to have the torque ratchets recalibrated by us once a year.

HARD METAL BONE CUTTER




^{**} Emergency tool for retrieving KOC® B


HANDGRIP SELF LOCKING

Please note the cleaning instructions on p. 57-58-61 of this system application brochure.

DRILLS



HANDGRIP TRAY



Size of closed tray **B** 195 mm **T** 90 mm **H** 45 mm For all autoclaves

Description	Length	REF	Description	Length	REF
BCD 1 Adapter	100 mm	13-310511	Handgrip	110 mm	13-311430
Twist Drill 2.0	110 mm	13-310512			
Adapter AHK	70 mm	13-462319	Handgrip tray empty		13-60043
Adapter AHB	70 mm	13-900037	Handgrip tray with content		13-S60043

ACCESSORIES



DescriptionREFReplacement O-ring for demountable13-311430-ORhandgrip 311430

INSTRUMENT TRAY FOR KOC® & BECES®




Size of closed tray **B** 175 mm **T** 145 mm **H** 65 mm

For all autoclaves. Autoclaveable up to 134° C, not suitable for dry heat sterilizers.

Description	System	Head	REF	Description	System	REF	Price €
IT2 BECES	KOC/BECES	large	13-900030	Twist Drill 2.0 30	BECES	13-90020	
IT2 S BECES	KOC/BECES	large	13-900038	Twist Drill 2.0 21	BECES	13-90022	
IT2 W	KOC/BECES	large	13-900039	Twist Drill 2.5 21	BECES	13-90026	
IT K	KOC/BECES	small	13-462320	BCD 1	KOC/BECES	13-900240	
ITS K	KOC/BECES	small	13-462322	BCD 2	KOC/BECES	13-900241	
ITW K	KOC/BECES	small	13-462331	BCD 3	KOC/BECES	13-900242	
ITWH K	KOC/BECES	small	13-462323	BCDX 1	KOC/BECES	13-900243	
DOS 1	KOC		13-455311	BCDX 2	KOC/BECES	13-900244	
DOS 2	KOC		13-455312	BCDX 3	KOC/BECES	13-900245	
DOS 3	KOC		13-455313	CDG	KOC/BECES	13-420329	
DOS 4	KOC		13-455314	CDG	KOC/BECES	13-420329	
DOS 5	KOC		13-455315	DX 2	KOC/BECES	13-500704	
C-Drill KM 1	KOC		13-455300	TW 2	KOC/BECES	13-425402	
C-Drill KM 2	KOC		13-455301				
C-Drill KM 3	KOC		13-455302	Instrument tray emp	oty	13-60006-K	upon request
IT LOC K	KOC		13-462333	Instrument tray with	content	13-S60006-K	upon request
DS 2	KOC		13-425001				
IT TB K	KOC		13-462327	The content for the system BECES® is optional			

DRILL-STOP TRAY



IT HAS BEEN SCIENTIFICALLY PROVEN

Heatless® drills by Dr. Ihde Dental generate 55 % less heat than traditional bone drills from other manufacturers. This makes it possible to use higher rotational speeds: between 3,000 and 5,000 rpm are recommended with good external cooling and intermittent drill technique.

STARTER TRAY

Autoclaveable up to 134° C, not suitable for dry heat sterilizers. This surgical kit contains all drills and tools for first works with the system KOC°.

Material: autoclaveable plastic.



Description	REF	Price €
IT K	13-462320	
ITS K	13-462322	
C-Drill KM 1	13-455300	
C-Drill KM 2	13-455301	
C-Drill KM 3	13-455302	
IT 2 BECES	13-900030	
IT 2 S BECES	13-900038	
DOS 1	13-455311	
DOS 2	13-455312	
DOS 3	13-455313	
BCDX 1	13-900243	
Torque wrench TW2	13-425402	
HT 1.25	13-425100	nt nt
ITX MU 15	13-418203	optional
Starter tray empty	13-60041-K	upon request
Starter tray with content	13-S60041-K	upon request

KOC® IMPLANTATIONS

Minimally invasive immediate load implantology with $\textbf{KOC}\ensuremath{\text{\circ}}$ implants

- 1. Insertion of 21 **KOC®** implants in three hours.
- 2. Maxilla restoration with new motorized screwing technology.
- 3. Minimally invasive approach for beginners.

Authors

Dr. Werner Mander, (IMF) Dr. Thomas Fabritius



Description	REF	Price cat.
DVD	13-6668	Α

INDICATIONS KOC® II KOC® MICRO

- Anchorage of crowns, bridges and bars, with the presence of adequate bone supply in terms of bone quality, bone width and bone height.
- Anchorage of prostheses via bar and button anchorage systems.
- Not for use in combination with simultaneous bone augmentations

RESTRICTIONS FOR KOC® B APPLICATION

- These two implant types may only be used as support implants in the reduced-load area.
- Splinting of at least three and possibly several implants for cross arch stabilisation.
- At least one KOC® or KOC® Micro implant must be involved in the construction.
- The prosthetic restoration must be securely fixed (with definitive cements).
- Not to be used for segmented bridges without the involvement of at least two \mathbf{KOC}^{\otimes} screws.
- $\bullet \quad \text{If in doubt, angulation adapters on KOC$^{\circ}$ screws are preferable to the KOC$^{\circ}$ B implant.}$
- $\bullet \quad \text{Not to be used for additional abutments in combination with natural teeth.} \\$
- Not to be used under off-axis load as well as in deep-bite cases in the maxillary and mandibular anterior region.
- Max. width of occlusal surface 5 mm.
- Not to be used as terminal abutments.
- Bendable up to 13 degrees.

NOTES ON THE CARE OF SURGICAL STEEL INSTRUMENTS

Surgical steel instruments can quickly become damaged if inadequately or improperly cared for. Only the special solvents for cleaning surgical steel should be used; in case of doubt, consult **Dr. Ihde Dental GmbH / AG**.


The following are not recommended:

- Disinfectants/cleaners with a high chlorine content.
- Disinfectants/cleaners with a high oxalic acid content.

For instruments with colour coding, the following are $\underline{\mathsf{NOT}}$ recommended:

- Excessively high solvent concentrations, disinfectants/cleaners with the components mentioned above.
- Excessive temperatures during cleaning and sterilization (no dry heat sterilization).

KOC® M & MICRO IMPLANTS



CORTICAL MILLING FOR KOC® MICRO

Description		Code	REF	Price cat.
C-Drill KM1 3.0 - 3.2	Cortical milling	C-Drill KM1	13-455300	Е
C-Drill KM2 3.7 - 4.1	Cortical milling	C-Drill KM2	13-455301	Е
C-Drill KM3 5.0	Cortical milling	C-Drill KM3	13-455302	E

MATERIAL

Ti6AL4V, also known as "Grade 5", is the high-purity version of the conventional 6/4 Ti alloy, which is used for more than 50% of all metallic human implants. This material is the first choice for all applications which require high stability, corrosion resistance and mechanical strength. This is why today's most modern dental implant designs are made of this material. This titanium alloy is superior to the alternatively used pure titanium in terms of stability by more than 25%. Also regarding biocompatibility and the support of bone cell growth, this titanium alloy shows advantages compared to pure titanium.

FUNCTIONALITY

The one-piece KOC® M / KOC® Micro dental implant is preferably used in immediate loading. Unlike other compression screws, the polished neck has a cylindrical shape. Thus, the cortical is sealed, good retention is achieved in the cortical and the bone is protected against infections. At the same time, the endosseous implant part compresses the cancellous bone areas.

NOTE

The smooth microthread must be completely submerged below the bone level. The cylindrical neck must extend into the bone at least 1 mm deep. Therefore, the implant must be selected so that at least 1.5 mm more usable vertical bone is present than the nominal length of the implant.

Example For KOC Micro 3.7 15, 17 mm of usable vertical bone must be present. If in doubt, a shorter implant should be selected so as to ensure a sufficient insertion depth.

DRILLING PROCEDURE

The pilot hole is made with the drills of the KOC* system. Except in very dense mandibular bone, the pilot hole is usually sufficient with BCD1 or DOS1.

INSERTION

The implant can be inserted most easily with the handgrip (REF 311431) and the adapter (REF 13-900 037). When using the ratchet RAT2, small or medium insertion tools are used. Max. torque is 80 Ncm.

THE IMPLANTS ARE SUPPLIED WITH TWO DIFFERENT HEAD SIZES

- KOC® M implants are supplied with a small head; they also fit in small individual tooth gaps.
- KOC® Micro implants are supplied with a large head. This head permits easy and speedy prosthetic restoration.

KOC® M IMPLANTS WITH SMALL ABUTMENT HEAD

 KOC^*M with small head for the anterior tooth region and tight gaps



Description	Endosseous Ø	Length	Neck Ø	Drill *	REF	Price cat.
KOC M 3.0 10	3.0 mm	10 mm	2 mm		13-456108	F
KOC M 3.0 12	3.0 mm	12 mm	2 mm	DOS 1 or BCD 1	13-456109	F
KOC M 3.0 15	3.0 mm	15 mm	2 mm		13-456110	F
KOC M 3.2 12	3.2 mm	12 mm	2 mm		13-456111	F
KOC M 3.2 15	3.2 mm	15 mm	2 mm		13-456112	F
KOC M 3.7 6	3.7 mm	6 mm	2 mm		13-456106	F
KOC M 3.7 8	3.7 mm	8 mm	2 mm		13-456107	F
KOC M 3.7 10	3.7 mm	10 mm	2 mm	DOS 2 or BCD 2	13-456114	F
KOC M 3.7 12	3.7 mm	12 mm	2 mm	DOSZOFBEDZ	13-456115	F
KOC M 3.7 15	3.7 mm	15 mm	2 mm		13-456120	F

^{*} In very hard bone, it may be additionally necessary to make a cylindrical hole with a twist Drill 2.5 mmd to a depth of 2.5 mm.


a) endosseous length	6 - 15 mm
b) max endosseous Ø small head	3.7 mm
c) neck Ø	2.0 mm
d) Height of head	6.8 mm
Material	Ti6Al4V

KOC® M implants are delivered incl. lab-set consisting of REF 13-462111, 13-462029 and 13-462088



KOC® MICRO IMPLANTS WITH LARGE ABUTMENT HEAD

KOC® Micro with large head for all applications



Description	Endosseous Ø	Length	Neck Ø	Drill*	REF	Price cat.
KOC Micro 3.7 6	3.7 mm	6 mm	2.0 mm	DOS 2 / BCD 2	13-455910	F
KOC Micro 3.7 8	3.7 mm	8 mm	2.0 mm	DOS 2 / BCD 2	13-455911	F
KOC Micro 3.7 10	3.7 mm	10 mm	2.0 mm	DOS 2 / BCD 2	13-455912	F
KOC Micro 3.7 12	3.7 mm	12 mm	2.0 mm	DOS 2 / BCD 2	13-455913	F
KOC Micro 3.7 15	3.7 mm	15 mm	2.0 mm	DOS 2 / BCD 2	13-455914	F
KOC Micro 4.1 8	4.1 mm	8 mm	2.0 mm	DOS 3 / BCD 3	13-455920	F
KOC Micro 4.1 10	4.1 mm	10 mm	2.0 mm	DOS 3 / BCD 3	13-455921	F
KOC Micro 4.1 12	4.1 mm	12 mm	2.0 mm	DOS 3 / BCD 3	13-455922	F
KOC Micro 4.1 15	4.1 mm	15 mm	2.0 mm	DOS 3 / BCD 3	13-455923	F
KOC Micro 5 10	5.0 mm	10 mm	2.0 mm	DOS 5	13-455925	F
KOC Micro 5 12	5.0 mm	12 mm	2.0 mm	DOS 5	13-455926	F

^{*} In very hard bone, it may be additionally necessary to make a cylindrical hole with a twist Drill 2.5 mmd to a depth of 2.5 mm.

a) endosseous length
b) max endoss. Ø large head
c) neck Ø
d) height of head
e) height of neck
Material

6 - 15 mm

3.9 mm

2.0 mm

ACCESSORIES

7.2 mm

Analogue IAB

REF 13-462106

3.5 mm

Pack of 5

Price cat. B

Ti6Al4V





Impression post **TSPA 5**Pack of 5

KOC® Micro implants

REF 13-462030 Price cat. B

Red impression caps and red analogue are round.

KOC® MICRO D IMPLANTS WITH LARGE ABUTMENT HEAD



Descr	iption	Endosseous Ø	Length	Neck Ø	Drill *	REF	Price cat.
KOC N	1icro D 3.7 8	3.7 mm	8 mm	2.0 mm	DOS 2 / BCD 2	13-455951	G
KOC N	Micro D 3.7 15	3.7 mm	15 mm	2.0 mm	DOS 2 / BCD 2	13-455954	G
KOC N	Nicro D 4.1 8	4.1 mm	8 mm	2.0 mm	DOS 3 / BCD 3	13-455957	G
KOC N	Micro D 4.1 10	4.1 mm	10 mm	2.0 mm	DOS 3 / BCD 3	13-455958	G
KOC N	1icro D 4.1 15	4.1 mm	15 mm	2.0 mm	DOS 3 / BCD 3	13-455960	G
KOC N	1icro D 5.0 12	5 mm	12 mm	2.0 mm	DOS 5	13-455964	G

* In very hard bone, it may be additionally necessary to make a cylindrical hole with a twist Drill 2.5 mmd to a depth of 2.5 mm.

a) endosseous length
b) max endoss. Ø large head
c) neck Ø
d) Height of head
Material

6 - 15 mm 3.9 mm 2.0 mm 7.2 mm Ti6Al4V KOC® Micro D implants are delivered incl. lab-set consisting of REF 13-462111, 13-462136 and 13-462086



IMPRESSION TAKING AND LABORATORY ACCESSORIES

		Description	Unit	Material	Code	REF	Price cat.
FOR SMALL HEAD		Impression cap	Pack of 5	POM	TSPA 4	13-462029	В
		Castable abutment	Pack of 5		PO4	13-462088	В
		Double analogue for small and large head	Pack of 5	plastic	IA4/IAU	13-462111	В
		Double analogue		metal	IA4/IAU	13-462112	A
FOR LARGE HEAD	PA X will be sold in the color	Impression post castable Internally edged	Pack of 5		PA X	13-462136	A
	orange, starting 01-2018	Castable abutment	Pack of 5		РОВ	13-462086	В



TITANIUM CROWN BASE



Description	Material	Code	REF	Price cat.
Titanium cap, radio opaque	Ti6Al4V, weldable	MA4	13-462090	В
Titanium cap, radio opaque	Ti6Al4V, weldable	MA5	13-462093	В

SCANBODIES MATERIAL PEEK/POM



 Description
 Scanbody-4 cylyndrical

 Systems
 KOC, BECES

Price cat. B (Pack of 5)

13-462054

REF



DescriptionScanbody-5 cylyndricalSystemsKOC, BECESREF13-462055

Price cat. B (Pack of 5)



Description Scanbody-MU cylyndrical

Systems BECES® MU, KOC® MU, Hexacone® MU

REF 13-462056

Price cat. B (Pack of 5)



Description Flag-Scanbody **SCB4** for small head. For intra-oral scan.

Systems KOC, BECES

REF 13-462071

Price cat. C (Pack of 5)



Description Flag-Scanbody **SCB5** for large head. For intra-oral scan.

Systems KOC, BECES

REF 13-462072

Price cat. C (Pack of 5)



Description Flag-Scanbody **SCB MU** incl. screw SFK MU (13-418164). For intra-oral scan.

Systems BECES® MU, KOC® MU, Hexacone® MU

REF 13-462073

Price cat. B (Pack of 1)

Please go to http://simpladent-implant.com/en/stl to download the corresponding STL files.

HEATLESS® DRILLS **DOS** FOR IMPLANTS WITH CONICAL CORE

Surgical steel, color-coded, depth-coded and autoclaveable. The drill is marked with laser depth markings. Use between 3,000 and 5,000 rpm with good cooling and intermittent drill technique. Due to the extremely high cutting performance, you can work without pressure.



	Description	Colour	Max. working length	REF	Price cat.
↑ DOS 1	DOS 1	yellow	17 mm	13-455311	D
I DOS 2	DOS 2	black	17 mm	13-455312	D
1000 1	DOS 3	red	17 mm	13-455313	D
0001	DOS 4	blue	21 mm	13-455314	D
	DOS 5	green	17 mm	13-455315	D
	DOS 6	transparent	15 mm	13-455316	D


DOS 6 This drill is 2 mm shorter at the tip. It can therefore drill up to 2 mm deeper into hard bone than nominally indicated on the drill. Therefore, the conical bone cavity is only circularly extended in the crestal area without increasing the drilling depth.

INSTRUMENTS AND **TOOLS**

	Description	Length	Code	REF	Price cat.
	Insertion tool short, for large head. Use with RAT2 and TW2.	7 mm	IT2 S BECES	13-900038	E
	Insertion tool medium, for large head. Use with RAT2 and TW2.	19 mm	IT2 BECES	13-900030	E
•	Insertion tool for large head. Use with angle piece.	23 mm	IT2W	13-900039	С
	Insertion tool long, for small head. Use with RAT2 and TW2.	20 mm	IT K	13-462320	С
	Insertion tool short, for small head. Use with RAT2 and TW2.	7 mm	ITS K	13-462322	С
	Insertion tool for small head. Use with angle piece.	23 mm	ITW K	13-462331	С
	Torque wrench 10 - 70 Ncm. It is recommended to have the torque ratchets recalibrated by us once a year.	;	TW2	13-425402	S
0	Adapter for large head, use with handgrip.	70 mm	АНВ	13-900037	F
0	Adapter for small head, use with handgrip.	70 mm	АНК	13-462319	D
	Handgrip. For machine processing in the ultrasonic bath.	110 mm		13-311431	V



KOC® MINI IMPLANTS



Description	Endosseous Ø	Endosseous length	Height above bone	Material	Max. insertion torque
KOC Mini 1.8 10	1.8 mm	10 mm	5.4 mm	Ti6AL4V	25 Ncm
KOC Mini 1.8 12	1.8 mm	12 mm	5.4 mm	Ti6AL4V	25 Ncm
KOC Mini 1.8 15	1.8 mm	15 mm	5.4 mm	Ti6AL4V	25 Ncm
KOC Mini 2.4 12	2.4 mm	12 mm	5.4 mm	Ti6AL4V	40 Ncm
KOC Mini 2.4 15	2.4 mm	15 mm	5.4 mm	Ti6AL4V	40 Ncm
KOC Mini 2.4 18	2.4 mm	18 mm	5.4 mm	Ti6AL4V	40 Ncm
KOC Mini 2.8 12	2.8 mm	12 mm	5.4 mm	Ti6AL4V	60 Ncm
KOC Mini 2.8 15	2.8 mm	15 mm	5.4 mm	Ti6AL4V	60 Ncm
KOC Mini 2.8 18	2.8 mm	18 mm	5.4 mm	Ti6AL4V	60 Ncm


SURGERY


1. Drilling and drill depth

KOC* Mini implants are inserted into the bone with self-tapping and with bone compression. As a rule, an undersized pilot hole is sufficient.

Recommended drill-speed: min. 3.000 Rpm, with sufficient cooling and intermittent drilling technique.

Note regarding hard bone: The torque must not exceed 45 Ncm during insertion. Should this be the case, the implant must be rotated out and the drilling depth extended to 2/3 of the implant length.





Drilling in hard bone: 1/2 of the endosseous implant length

SAME PROCEDURE WITH REF13-462374 / REF13-462376!

2. Removing the implant from the package














3. insertion and placement





KOC® I	MINI	SET 1	.8 10
--------	------	-------	-------

KOC® Mini 1.8 10 (pack of 4) Drill for KOC® Mini DS, ø 1.3 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455050 0 KOC® MINI SET 2.4 18

KOC® Mini 2.4 18 (pack of 4) Drill for KOC® Mini DS, ø 1.9 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455055 0

KOC® MINI SET 1.8 12

KOC® Mini 1.8 12 (pack of 4) Drill for KOC® Mini DS, ø 1.3 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455051

KOC® MINI SET 2.8 12

KOC® Mini 2.8 12 (pack of 4) Drill for KOC® Mini DS, ø 2.3 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455056

KOC® MINI SET 1.8 15

KOC® Mini 1.8 15 (pack of 4) Drill for KOC® Mini DS, ø 1.3 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455052 Q KOC® MINI SET 2.8 15

KOC® Mini 2.8 15 (pack of 4) Drill for KOC® Mini DS, ø 2.3 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455057 Q

KOC® MINI SET 2.4 12

KOC® Mini 2.4 12 (pack of 4) Drill for KOC® Mini DS, ø 1.9 mm Sleeve and O-ring (set of 4)

REF Price cat.

13-455053 Q KOC® MINI SET 2.8 18

KOC® Mini 2.8 18 (pack of 4) Drill for KOC® Mini DS, ø 2.3 mm Sleeve and O-ring (set of 4)

REF Price cat.


13-455058 Q

KOC® MINI SET 2.4 15

KOC® Mini 2.4 15 (pack of 4) Drill for KOC® Mini DS, ø 1.9 mm Sleeve and O-ring (set of 4)


REF Price cat.

13-455054











ACCESSORIES

	Description Insertion tool short, 12 mm	Working length	Code IT K KOC® Mini	REF 13-462376	Price cat.
	Insertion tool long, 20 mm		IT K KOC® Mini	13-462374	С
	Insertion tool (2-piece) for manual insertion		IT H KOC® Mini	13-462375	D
C THE	Torque wrench, 10 - 70 Ncm		TW2	13-425402	S
1	Impression coping KOC® Mini Pack of 4			13-462117	A
	Analogue KOC® Mini Pack of 4			13-462116	A
	Sleeve KOC® Mini for polymerizing into the prosthesis, delivery with inserted O-ring			13-462113	В
	O-ring for sleeve REF 13-462113 for replacement. Pack of 4			13-462114	В
	Depth-measurement-probe for KOC® Mini			13-462115	А
	Drill for KOC® Mini, ø 1.3 mm	15 mm		13-455318	С
	Drill for KOC® Mini, ø 1.9 mm	18 mm		13-455319	С
	Drill for KOC® Mini, ø 2.3 mm	18.4 mm		13-455320	С

APPLICATIONS OF THE STRATEGIC IMPLANT® FOR ANCHORAGE IN THE UPPER AND LOWER JAW



BECES® implants can be used immediately in extraction sockets if the basal support is sufficient. The anti-rotation protection ensures immediate stability against unintentional unscrewing before prosthetic loading. The prosthesis should be inserted before the 3rd post-operative day. BECES® implants are made of strong, biocompatible titanium alloy Ti6Al4V ELI. BECES® implants are used typically for segments and circular bridges in an immediate splinting protocol. Their use is permitted only for authorized users. BECES® EX implants provide a sharp and cutting central thread part, which engages in to the corticals of the extraction sockets. The apical compression thread provides perfect stability both in compressed spongious and in cortical bone. The coronal micro-thread seals the 1st cortical if the implant is used in healed bone areas. BECES® EX implants may be used both in extraction sockets and in healed bone areas. They are used for circular bridges and segments. Under adequate loading conditions and if enough implants are splinted, the treatment can be performed in an immediate load protocol. Due to their polished surface, BECES® EX implants are extremely resistant against bacterial colonialization and they avoid peri-implantitis.

For implant Ø from 3.5 to 4.6 mm, not less than ten implants should be used for the complete upper jaw.

For implant Ø from 3.5 to 4.6 mm, not less than eight implants should be used for the complete mandible.

The phase of adaptation and consolidation of the cortical bone in which these implants are anchored (2nd or 3rd cortical) must be considered to be 24 months after the insertion of the implants. So there is a significant difference to the healing times for implants that are used according to the method of "osseointegration", for which the healing times are assumed to be 3-6 months.

Nevertheless BECES® implants are designed to be used (solely by specially trained and authorized users) in immediate functional loading protocols.

FITTING AND CEMENTING OF PROSTHETICS

The lower border of the abutment head of the Strategic Implant* is (only) used as a margin to hold the transfer during impression-taking. Because the implant and the abutment head are both polished, the lower margin of the implant does not typically serve as a crown margin as we know it from teeth or conventional 2-stage implants. There are no medical or technical reasons why the crown margin (or the margin of the technical abutment) should reach the lower border of the abutment head.

It is important however that enough distance between the lower margin of the prosthetic workpiece and the gums (or the bone respectively) is left after cementation. We recommend to use only strong permanent cements (e.g. Fuji Plus, GC Corp.) and to have a vertical cementing surface/zone of at least 4 mm on the abutments. The abutment head may be shortened/adjusted vertically and/or laterally in order to achieve a good aesthetic result and to allow good phonetics.

Those surfaces on the abutment head which will provide retention for the cement must be roughened and cleaned before cementation. All other surfaces of the abutment head must remain fully polished.

The main aim of this step of the treatment is the incorporation of a prosthetic workpiece which is easily cleanable or which allows self-cleaning (on the lingual or palatal side) in function.




Fig. 1

The transfer cap (3) is positioned on the abutment head until the lower border of the abutment head (1) is reached. The transfer will sit firmly in this position. Then the impression is taken with silicone putty or heavy body silicone material. This allows the transfer of the implant position to the model.




Fig. 2

The implant head was placed approximately 1.8 mm above the bone level (4). The mucosa level (5) reaches higher than the lower border of the abutment head (1). The level of the crown margin (6) and the lower border of the abutment head (1) are in a distance to each other. This avoids retention of cements and debris in the submucosal area. This is a correct result. On the x-ray the crown will appear however as "too short", considering not applicable criteria from conventional dentistry.



Fig. 3

The crown margin (6) will be the same level as the lower border of the abutment head (1) if the abutment head sits on the mucosa line (5).



Fig. 4

If the abutment head is positioned on a mucosal slope, the lower border of the abutment head is on one side (1a) deeper in the mucosa than on the other side (1b). In such a case the crown margin (6) will also run oblique, in order to avoid submucosal position of parts of the crown. See the clincial example in Fig. 9. Also in this case the crown may appear as "too short" on the x-ray, considering not applicable criteria from conventional dentistry.






Fig. 5

For aesthetic reasons it may be necessary to create vestibular overhangig portions of the prosthetic workpiece (7).




Fig. 6

It is not allowed to create such prosthetic overhangs (7a, 7b) on both sides of the prosthetic workpice, because this would lead to an non-hygienic situation without the possibility of self cleaning. Food and debris will get stuck in the area of the mucosal penetration area of the implant and this will create an inflammation.




Fig. 7

If vertical height is missing, the top part of the abutment head may be shortened (region 8 is removed). At the same time it might be necessary to keep a distance between the lower margin of the abutment head (1) and the lower crown margin (6). Nevertheless the vertical cementation area (9) should be not less than 4 mm heigh.




Fig. 8

If abutment heads are used as technical abutments, they are shortened after the final cementation of the prosthetic workpiece (region 8 is removed) and after the cement has fully set. This adjustment may be done at the first control appointment. They remain «open». The height of the cementing surface (9) should be not less than 4 mm. The lower margin of the crown does not necessarily concide with the lower border of the abutment head.



Fig. 9

The implant crowns 43 and 44 have been shortened more than 3 mm on the lingual side and on the vestibular side an overhang has been modelled. The necessary height for cementation is given both on the vestibular and the lingual side on the abutment head.

CONCLUSION

The question if the prosthetic construction is properly fitted to the abutment of the Strategic Implant * depends on the spational relationship between the crown margin to the mucosa much more than on anything else. Relevant for any judgement about the length of the crown is the moment of the cementation.

 $Only for selected bridge \ materials \ and \ bridge \ designs, subgingival \ connection \ between \ implant \ abut ment \ and \ prosthetics \ is \ possible. \ In \ such \ cases \ the \ final \ connection \ between \ the \ two \ components \ requires \ an \ open \ surgical \ cementation.$



BECES® IMPLANTS 2.7 MMD WITH SMALL ABUTMENT HEAD

These implants are used for the following indications

- Supporting (additional) implants for cortical anchorages of bridges and crowns
- Creation of a three-point support for the cortical anchorage of dental prostheses



Description
BECES 2.7 10
BECES 2.7 12
BECES 2.7 14
BECES 2.7 17
BECES 2.7 20
BECES 2.7 23
BECES 2.7 26
BECES 2.7 29
BECES 2.7 32

c	d	e	g	Drill	REF	Price cat
10	4.5	2.7	2.55	Twist Drill 1.8	13-900190	G
12	4.5	2.7	2.95	Twist Drill 1.8	13-900191	G
14	5.5	2.7	2.95	Twist Drill 1.8	13-900192	G
17	5.5	2.7	2.95	Twist Drill 1.8	13-900193	G
20	5.5	2.7	2.95	Twist Drill 1.8	13-900194	G
23	5.5	2.7	2.95	Twist Drill 1.8	13-900195	G
26	5.5	2.7	2.95	Twist Drill 1.8	13-900196	G
29	5.5	2.7	2.95	Twist Drill 1.8	13-900197	G
32	5.5	2.7	2.95	Twist Drill 1.8	13-900198	G

a) Max. abutment \emptyset 3.35 mm b) Abutment height 6.8 mmh c) Nominal length 10 - 32 mm d) Length of apical thread 4.5 / 5.5 mm e) Enossal Ø max. 2.7 mm f) Neck \emptyset in bending zone 1.9 mm 2.55 - 2.95 mm g) Length of bending zone h) Square AF (across flats) 1.9 mm Tool ITK, AHK

USE LIMITATIONS BECES 2.7 must not be used as an implant for single tooth replacement, however two or more BECES 2.7 may serve as such. If **only** BECES 2.7 is used in very thin jaws, the surgeon should try to insert at least eight, but better more (up to 12 implants) for this jaw. BECES 2.7 are considered additional dental implants and they are used with other BECES implants 3.5 mm - 12 mm in order to increase the stability of the implant-prosthetic system.

TWIST DRILL



Description							
Twist Drill	1.8/23						

Ø Max. working length REF1.8 mm 23 mm 13-9

REF Price cat.

13-90024 D




BECES® IMPLANTS 3.0 MMD

BECES® IMPLANTS WITH SMALL ABUTMENT HEAD

These implants are used for the following indications

- Supporting (additional) implants for cortical anchorages of bridges and crowns
- Creation of a three-point support for the cortical anchorage of dental prostheses



Description	a	b	e	g	Drill	REF	Price cat.
BECES 3.0 10	3.0	10	2.55	4.5	Twist Drill 1.8	13-900480	G
BECES 3.0 12	3.0	12	2.95	4.5	Twist Drill 1.8	13-900481	G
BECES 3.0 14	3.0	14	2.95	5.5	Twist Drill 1.8	13-900482	G
BECES 3.0 17	3.0	17	2.95	5.5	Twist Drill 1.8	13-900483	G
BECES 3.0 20	3.0	20	2.95	5.5	Twist Drill 1.8	13-900484	G
BECES 3.0 23	3.0	23	2.95	5.5	Twist Drill 1.8	13-900485	G
BECES 3.0 26	3.0	26	2.95	5.5	Twist Drill 1.8	13-900486	G
BECES 3.0 29	3.0	29	2.95	5.5	Twist Drill 1.8	13-900487	G
BECES 3.0 32	3.0	32	2.95	5.5	Twist Drill 1.8	13-900488	G

USE LIMITATIONS BECES 3.0 must not be used as an implant for single tooth replacement, however two or more BECES 3.0 may serve as such. If **only** BECES 3.0 is used in very thin jaws, the surgeon should try to insert at least eight, but better more (up to 12 implants) for this jaw. BECES 3.0 are considered additional dental implants and they are used with other BECES implants 3.5 mm - 12 mm in order to increase the stability of the implant-prosthetic system.

a) thread \varnothing max. 3.0 mm b) endosseous length 10 - 32 mm c) abutment Height 6.8 mm d) max. abutment \varnothing 3.35 mm e) length of bending zone 3.0 mm f) shaft diameter in the bending zone 1.9 mm

g) length of the apical thread 4.5 / 5.5 mm (depending on the endosseous implant length)

Tool IT K, AHK

TWIST DRILL



 Description
 Ø
 Max. working length
 REF
 Price cat.

 Twist Drill 1.8/23
 1.8 mm
 23 mm
 13-90024
 D

PATHFINDER DRILLS

Conical 3-edge drill as initial drill, ideally suited for all crestal implant systems. The drill also passes between narrow cortical areas without pressure.



Description	Colour
BCDX 1	yellow

15 mm

13-900243 C

BECES® IMPLANTS WITH SMALL ABUTMENT HEAD

BECES® implants 3.5 mmd - 12 mmd

For anchorage in the 1^{st} , 2^{nd} and if necessary 3^{rd} cortical, for the cortical anchorage of dental prostheses. **BECES®** implants can be used in sockets for a given indication immediate 1^{st} , 1^{st diately after extraction and loaded immediately in many cases. Mechanically smoothed surface in all areas. The abutment head is identical to the head of KOC* implants and the TSD4 a butment. Self-tapping thread with endosseous anti-rotation protection. Conditionally suitable for individual tooth prostheses.

		Description	a	b	f	REF	Price cat.
BECES 3.5 17	d	BECES 3.5 10	3.5	10	5.5	13-900208	G
		BECES 3.5 12	3.5	12	5.5	13-900226	G
	c	BECES 3.5 14	3.5	14	7.5	13-900210	G
	-	BECES 3.5 17	3.5	17	7.5	13-900211	G
		BECES 3.5 20	3.5	20	7.5	13-900212	G
	b	BECES 3.5 23	3.5	23	7.5	13-900213	G
f	7	BECES 3.5 26	3.5	26	7.5	13-900214	G
	Ŧ	BECES 3.5 29	3.5	29	7.5	13-900215	G
1	T	BECES 3.5 32	3.5	32	7.5	13-900216	G
	a	BECES 3.5 35	3.5	35	7.5	13-900217	G
		BECES 3.5 38	3.5	38	7.5	13-900218	G
		BECES 4.5 10	4.5	10	7.5	13-900238	G
		BECES 4.5 12	4.5	12	7.5	13-900239	G
		BECES 4.5 14	4.5	14	7.5	13-900220	G
		BECES 4.5 17	4.5	17	7.5	13-900221	G
		BECES 4.5 20	4.5	20	7.5	13-900222	G
		BECES 4.5 23	4.5	23	7.5	13-900223	G
		BECES 4.5 26	4.5	26	7.5	13-900224	G
		BECES 4.5 29	4.5	29	7.5	13-900225	G

Insertion tools: IT KOC, ITX KOC, ITS KOC, Adapter AHK

a) max. thread \emptyset 3.5 - 4.5 mm b) nominal length 10 - 38 mm c) height abutment 6.8 mm d) max. abutment Ø 3.35 mm e) max. crestal Ø of shaft 2.0 mm f) length of thread 5.5 - 7.5 mm Square key width 1.9 mm Max. insertion torque 80 Ncm

BECES® implants are delivered incl. lab-set



FIELD OF APPLICATION

Endosseous dental implant for $2^{\rm nd}$ cortical anchorage. For full upper jaws we recommend the usage of 10 - 12 or more implants, for full lower jaws the usage of 8 - 10 implants or more. For unilateral segments we recommend to use 4 - 6 implants.

BECES® IMPLANTS WITH LARGE ABUTMENT HEAD

		Description	a	b	f	REF	Price cat.
BECES 3.6 17		BECES 3.6 10	3.6	10	5.5	13-900285	Н
		BECES 3.6 12	3.6	12	5.5	13-900284	Н
HIIIA	С	BECES 3.6 14	3.6	14	7.5	13-900286	Н
		BECES 3.6 17	3.6	17	7.5	13-900287	Н
		BECES 3.6 20	3.6	20	7.5	13-900288	Н
e e		BECES 3.6 23	3.6	23	7.5	13-900289	Н
		BECES 3.6 26	3.6	26	7.5	13-900290	Н
	b	BECES 3.6 29	3.6	29	7.5	13-900291	Н
T .							
		BECES 4.6 8	4.6	8	5.5	13-900299	Н
f H		BECES 4.6 10	4.6	10	5.5	13-900292	Н
		BECES 4.6 12	4.6	12	5.5	13-900300	Н
TH		BECES 4.6 14	4.6	14	7.5	13-900293	Н
a		BECES 4.6 17	4.6	17	7.5	13-900294	Н
u		BECES 4.6 20	4.6	20	7.5	13-900295	Н
		BECES 4.6 23	4.6	23	7.5	13-900296	Н
		BECES 4.6 26	4.6	26	7.5	13-900297	Н
		BECES 4.6 29	4.6	29	7.5	13-900298	Н
BECES 4.6 14 BECES 5	.5 14						
(1) (1)		BECES 5.5 8	5.5	8	5.5	13-900255	K
100 100		BECES 5.5 10	5.5	10	5.5	13-900281	K
	Ì	BECES 5.5 12	5.5	12	6.0	13-900250	K
		BECES 5.5 14	5.5	14	6.0	13-900251	K
		BECES 5.5 17	5.5	17	6.0	13-900252	K
社		BECES 5.5 20	5.5	20	6.0	13-900253	K
茅世		BECES 5.5 23	5.5	23	6.0	13-900265	K

Insertion tools: IT2 BECES, IT2 S BECES, Adapter AHB

BECES 5.5 26

BECES 5.5 29

a) max. thread Ø 3.6 - 5.5 mm b) nominal length 8 - 29 mm c) length abutment 7.2 mm d) abutment \emptyset 3.9 mm e) max. crestal \emptyset of shaft 2.0 mm f) length of thread 5.5 - 7.5 mm

BECES® Implants with large head are delivered incl. lab-set consisting of REF 13-462111, 13-462136 and 13-462086

6.0

26

29



13-900266

13-900267

FIELD OF APPLICATION BECES 3.6 MMD - 5.5 MMD

 $Endos seous dental implant for 2^{nd}/3^{nd} cortical \ anchorage. We recommend for full upper jaws the usage of 10-12 implants or more, for full lower jaws the usage of 8-10 implants or more. \\$ or more. For unilateral segments we recommend to use 4-6 implants.

BECES® IMPLANTS WITH LARGE ABUTMENT HEAD

Description

BECES 7.0 8

BECES 7.0 10

BECES 7.0 12

BECES 7.0 14

BECES 7.0 17

BECES 7.0 20

BECES 9.0 8

BECES 9.0 10

BECES 9.0 12

BECES 9.0 14

BECES 10.5 10

BECES 10.5 12

BECES 10.5 14

	c
f	b

BECES 7 12 BECES 9 14 BECES 10.5 14



		BECES 10.5 17	10.5	17	2.1	6.5	13-900280	М
a) thread Ø	7 - 12 mm	BECES 12.0 8	12	8	2.1	5.5	13-900279	0
b) endosseous length	8 - 20 mm	BECES 12.0 10	12	10	2.1	5.5	13-900272	0
c) length abutment	7.2 mm	BECES 12.0 12	12	12	2.1	6.5	13-900275	0
d) abutment \emptyset	3.9 mm	BECES 12.0 14	12	14	2.1	6.5	13-900273	0

Insertion tools: IT2 BECES, IT2 S BECES, Adapter AHB e) max. crestal Ø of shaft 2.0, 2.1 mm

5.5 - 6.5 mm

 ${\it Max.}\, insertion\, torque\, 80\, Ncm$

f) length of thread

BECES® Implants with large head are delivered incl. lab-set

f

5.5

5.5

5.5

5.5

5.5

5.5

5.5

5.5

5.5

5.5

6.5

6.5

6.5

2.0

2.0

2.0

2.0

2.0

2.0

2.1

2.1

2.1

2.1

2.1

2.1

b

10

12

14

17

20

8

10

12

14

10

12

7

9

9

10.5

10.5

10.5

REF

13-900258

13-900282

13-900260

13-900261

13-900262

13-900263

13-900269

13-900270

13-900274

13-900271

13-900276

13-900278

13-900277 M

Price cat.

Κ

Κ

Κ

Κ

K

Κ

M

M

M

M

M

M



BECES® EX IMPLANTS

 $\textbf{BECES} \\ \textbf{EX} \\ \text{implants provide a sharp and cutting crestal thread, which engages in to the corticals of the extraction sockets. The apical compression thread provides perfect$ stability both in compressed spongious and in cortical bone. BECES® EX implants may be used both in extraction sockets and in healed bone areals. They are used for circular bridges and segments with at least three implants. Under adequate loading conditions and if enough implants are splinted, the treatment can be performed in an immediate load protocol. Due to their polished surface, BECES® EX implants are extremely resistant against bacterial colonialization and peri-implantitis. Material TI6AL4V ELI, triple thread.

BECES EX 3.5 + 15



Description	Max. endosseous Ø mm	Endosseous length	REF	Price cat.
BECES EX 3.5 + 12	3.5	12	13-420150	K
BECES EX 3.5 + 15	3.5	15	13-420151	K
BECES EX 3.5 + 17	3.5	17	13-420152	K
BECES EX 3.5 + 19	3.5	19	13-420153	K
BECES EX 3.5 + 21	3.5	21	13-420154	K
BECES EX 3.5 + 23	3.5	23	13-420155	K
BECES EX 3.5 + 26	3.5	26	13-420156	K
BECES EX 3.5 + 29	3.5	29	13-420157	K
BECES EX 4.1 + 12	4.1	12	13-420160	K
BECES EX 4.1 + 15	4.1	15	13-420161	K
BECES EX 4.1 + 17	4.1	17	13-420162	K
BECES EX 4.1 + 19	4.1	19	13-420163	K
BECES EX 4.1 + 21	4.1	21	13-420164	K
BECES EX 4.1 + 23	4.1	23	13-420165	K

BECES EX 4.1 + 15



a) Neck height


b) Neck diameter 2 mm (bendable up to 15°)

c) Max. endosseous Ø 3.5 / 4.1 mm d) Endosseous length 12 - 23 mm 3.9 w x 7.2 h (mm) Abutment

STRATEGIC IMPLANT® WITH SMALL ABUTMENT HEAD

STRATEGIC IMPLANT® 3.5 mmd - 12 mmd

For anchorage in the 1st, 2nd and if necessary 3rd cortical, for the cortical anchorage of dental prostheses. **Strategic Implant*** can be used in sockets for a given indication immediately after extraction and loaded immediately in many cases. Mechanically smoothed surface in all areas. The abutment head is identical to the head of **KOC*** implants and the **TSD4** abutment. Self-tapping thread with endosseous anti-rotation protection. Conditionally suitable for individual tooth prostheses.



Insertion tools: IT KOC, ITX KOC, ITS KOC, Adapter AHK

 a) max. thread Ø
 3.5 - 4.5 mm

 b) nominal length
 10 - 38 mm

 c) height abutment
 6.8 mm

 d) max. abutment Ø
 3.35 mm

 e) max. crestal Ø of shaft
 2.0 mm

 f) length of thread
 5.5 - 7.5 mm

 Square key width
 1.9 mm

Max. insertion torque 80 Ncm

Strategic Implant® are delivered **incl. lab-set** consisting of REF 13-462111, 13-462029 and 13-462088



FIELD OF APPLICATION

Endosseous dental implant for 2^{nd} cortical anchorage. For full upper jaws we recommend the usage of 10 - 12 or more implants, for full lower jaws the usage of 8 - 10 implants or more. For unilateral segments we recommend to use 4 - 6 implants.

STRATEGIC IMPLANT® IMPLANTS WITH LARGE ABUTMENT HEAD

Strategic Implant®		Description	a	b	f	REF	Price cat.
3.6 17		Strategic Implant® 3.6 10	3.6	10	5.5	14-900285	Н
		Strategic Implant® 3.6 12	3.6	12	5.5	14-900284	Н
ALI III	С	Strategic Implant® 3.6 14	3.6	14	7.5	14-900286	Н
-		Strategic Implant® 3.6 17	3.6	17	7.5	14-900287	Н
		Strategic Implant® 3.6 20	3.6	20	7.5	14-900288	Н
├── ' e		Strategic Implant® 3.6 23	3.6	23	7.5	14-900289	Н
		Strategic Implant® 3.6 26	3.6	26	7.5	14-900290	Н
- 111	b	Strategic Implant® 3.6 29	3.6	29	7.5	14-900291	Н
T 1		Strategic Implant® 4.68	4.6	8	5.5	14-900299	Н
		Strategic Implant® 4.6 10	4.6	10	5.5	14-900292	Н
f		Strategic Implant® 4.6 12	4.6	12	5.5	14-900300	Н
1		Strategic Implant® 4.6 14	4.6	14	7.5	14-900293	Н
		Strategic Implant® 4.6 17	4.6	17	7.5	14-900294	Н
a		Strategic Implant® 4.6 20	4.6	20	7.5	14-900295	Н
		Strategic Implant® 4.6 23	4.6	23	7.5	14-900296	Н
		Strategic Implant® 4.6 26	4.6	26	7.5	14-900297	Н
Strategic Implant® 4.6 14 5.5 14		Strategic Implant® 4.6 29	4.6	29	7.5	14-900298	Н
4.5 I4 5.5 I4	m	Strategic Implant® 5.5 8	5.5	8	5.5	14-900255	K
1111 1	١١.	Strategic Implant® 5.5 10	5.5	10	5.5	14-900281	K
		Strategic Implant® 5.5 12	5.5	12	6.0	14-900250	K
T I	Т	Strategic Implant® 5.5 14	5.5	14	6.0	14-900251	K
	ш	Strategic Implant® 5.5 17	5.5	17	6.0	14-900252	K
5	J.	Strategic Implant® 5.5 20	5.5	20	6.0	14-900253	K
世上	-	Strategic Implant® 5.5 23	5.5	23	6.0	14-900265	K
T 3		Strategic Implant® 5.5 26	5.5	26	6.0	14-900266	K

Insertion tools: IT2 BECES, IT2 S BECES, Adapter AHB

Strategic Implant® 5.5 29

 a) max. thread Ø
 3.6 - 5.5 mm

 b) nominal length
 8 - 29 mm

 c) length abutment
 7.2 mm

 d) abutment Ø
 3.9 mm

 e) max. crestal Ø of shaft
 2.0 mm

 f) length of thread
 5.5 - 7.5 mm

Strategic Implant® with large head are delivered **incl. lab-set** consisting of REF 13-462111, 13-462136 and 13-462086



ACCESSORIES

Analogue IAB REF 13-462106
Pack of 5 Price cat. B



Impression post **TSPA 5 REF** 13-462030

14-900267

Pack of 5 Price cat. B

Red impression caps and red analogue are round.

FIELD OF APPLICATION Strategic Implant® 3.6 MMD - 5.5 MMD

Endosseous dental implant for $2^{nd}/3^{nd}$ cortical anchorage. We recommend for full upper jaws the usage of 10-12 implants or more, for full lower jaws the usage of 8-10 implants or more. For unilateral segments we recommend to use 4-6 implants.

STRATEGIC IMPLANT® IMPLANTS WITH LARGE ABUTMENT HEAD

Strategic Implant® 12 12
d c
f e b

	U	
	а	
Stra	ategic Im	plant®
7 12	9 14	10.5 14
	4	4

Description	a	b	e	f	REF	Price cat.
Strategic Implant® 7.08	7	8	2.0	5.5	14-900258	K
Strategic Implant® 7.0 10	7	10	2.0	5.5	14-900282	K
Strategic Implant® 7.0 12	7	12	2.0	5.5	14-900260	K
Strategic Implant® 7.0 14	7	14	2.0	5.5	14-900261	K
Strategic Implant® 7.0 17	7	17	2.0	5.5	14-900262	K
Strategic Implant® 7.0 20	7	20	2.0	5.5	14-900263	K
Strategic Implant® 9.0 8	9	8	2.1	5.5	14-900269	M
Strategic Implant® 9.0 10	9	10	2.1	5.5	14-900270	M
Strategic Implant® 9.0 12	9	12	2.1	5.5	14-900274	M
Strategic Implant® 9.0 14	9	14	2.1	5.5	14-900271	M
Strategic Implant® 10.5 10	10.5	10	2.1	6.5	14-900276	M
Strategic Implant® 10.5 12	10.5	12	2.1	6.5	14-900277	M
Strategic Implant® 10.5 14	10.5	14	2.1	6.5	14-900278	M
Strategic Implant® 10.5 17	10.5	17	2.1	6.5	14-900280	M
Strategic Implant® 12.0 8	12	8	2.1	5.5	14-900279	0
Strategic Implant® 12.0 10	12	10	2.1	5.5	14-900272	0
Strategic Implant® 12.0 12	12	12	2.1	6.5	14-900275	0
Strategic Implant® 12.0 14	12	14	2.1	6.5	14-900273	0

Insertion tools: IT2 BECES, IT2 S BECES, Adapter AHB

a) thread Ø
b) endosseous length
c) length abutment
d) abutment Ø

e) max. crestal \emptyset of shaft

 ${\it Max.}\, insertion\, torque\, 80\, Ncm$

f) length of thread

i) leligiii oi tillead

7 - 12 mm 8 - 20 mm 7.2 mm 3.9 mm 2.0, 2.1 mm 5.5 - 6.5 mm

ACCESSORIES

Analogue IAB REF 13-462106 Pack of 5 Price cat. B

Impres

Strategic Implant® with large head are delivered **incl. lab-set** consisting of REF 13-462111,

Impression post **TSPA 5**

REF 13-462030 **Price cat. B**


Red impression caps and red analogue are round.

IMPRESSION TAKING AND LABORATORY ACCESSORIES FOR BECES® AND KOC® IMPLANTS

ALTERNATIVE	

B 1.4				
Description	Unit	Code	REF	Price cat.
Impression post castable, POM For small head Internally round	Pack of 5	TSPA 4	13-462029	В
Impression post castable, POM For small head Internally round	Pack of 5	TSPA 4	13-462027	В
Impression post castable, POM For large head Internally round	Pack of 5	TSPA 5	13-462030	В
Impression post castable Internally edged	Pack of 5	PA X	13-462136	В
Double analogue, metal	1 piece	IA4/IAU	13-462112	Α
Double analogue, plastic	Pack of 5	IA4/IAU	13-462111	В
Castable abutment and base for provisionals For small head 7 mm high, white, internally round	Pack of 5	PO4	13-462088	В
Castable abutment For large head Internally round	Pack of 5	РОВ	13-462086	В

GUIDE JACKET



Description	Unit	Material	REF	Price cat.
BFH 2.0 guide jacket 2.0 mmd	Pack of 5	Ti6Al4V	13-425410	В
BFH 2.5 guide jacket 2.5 mmd	Pack of 5	Ti6Al4V	13-425411	В
a) Length	5 mm			
b) Height of step	0.7 mm			
c) Max. Ø top	3.7 / 4 mm			
d) Nominal Ø	3 / 3.35 mm			
e) \varnothing of drilling in the drill template	2.05 / 2.55 mm			



Model with residual teeth for the fabrication of a drill guide for creating cavities for fixating the later drill guide for implant cavities.




Drill guide for creating cavities for later fixation of the surgical drill guide.



Surgical drill guide for safe BECES® placement. The drill sleeves are designed for 2.0 mm Twist drills.

HANDGRIP SELF LOCKING


For machine reprocessing, cannot be dismantled. Clean in an ultrasonic bath at 45° with an alkaline cleaning agent. For adapter, self-locking. Please note the cleaning instructions on www.implant.com/en/downloads







HANDGRIP TRAY



Size of closed tray **W** 195 mm **D** 90 mm **H** 45 mm For all autoclaves

Description	Length	REF	Price €
BCD 1 Adapter	100 mm	13-310511	
Twist Drill 2.0	110 mm	13-310512	
Adapter AHK	70 mm	13-462319	
Adapter AHB	70 mm	13-900037	
Handgrip	110 mm	13-311431	
Handgrip tray w/o content		13-60043	upon request
Handgrip tray with content		13-S60043	upon request

Please read our detailed instructions for cleaning and re-sterilization of surgical instruments on https://implant.com/en/downloads

PATHFINDER DRILLS

 $Conical\ 3-edge\ drill\ as\ initial\ drill, ideally\ suited\ for\ all\ crestal\ implant\ systems.\ The\ drill\ also\ passes\ between\ narrow\ cortical\ areas\ without\ pressure.$

		Description	Colour	Max. working length	REF	Price cat.
	39 mm	BCD 1	yellow	15 mm	13-900240	С
		BCD 2	black	15 mm	13-900241	С
		BCD 3	red	13 mm	13-900242	С
- BGDE 1		BCDX 1	yellow	15 mm	13-900243	С
:		BCDX 2	black	15 mm	13-900244	С
-	47 mm	BCDX 3	red	15 mm	13-900245	С
	NOR T	BCD 1 Adapter pathfinder for handgrip length 100 mm			13-310511	F

TWIST DRILLS

	Description Twist Drill 2.0/21 Twist Drill 2.0/30 Twist Drill 2.0/40	Ø 2.0 mm 2.0 mm	Max. working length 21 mm 30 mm 40 mm	REF 13-90022 13-90020 13-90019	Price cat. D D
	Twist Drill 2.5/21	2.5 mm	21 mm	13-90026	D
		Description	Max. working length	REF	Price cat.
- T.	***	Twist Drill 2.0 Cylindrical drill 2.0 mm for handgrip, length 110 mm	35 mm	13-310512	F
		Pilot drill for surgical handgrip. For chuck 2.35 mmd		13-310515	F
	ABOBOBO	Twist Drill 2.0/30 for surgical hand- grip. For chuck 2.35 mmd.	110 mm	13-310516	F

HARD METAL BONE CUTTER

Description	Length	Code	KEF	
Hard metal bone cutter short	30 mm	SHMCS	13-90030	F
Hard metal bone cutter long, for FG	36 mm	SHMCL	13-90031	F

INSERTION TOOLS AND **ADAPTER**

Description	Code	REF	Price cat.
For BECES® implants with Ø 3.5 mm + 4.5 mm	ITK	13-462320	С
Insertion tool medium, for large head. Use with RAT2 and TW2.	IT2 BECES	13-900030	E
For BECES® implants with Ø 3.6, 4.6, 5.5, 7.9, 10.5, 12 mm	IT2 S BECES	13-900038	Е
For BECES® implants with Ø 3.6, 4.6, 5.5, 7.9, 10.5, 12 mm. For handgrip REF 311431	АНВ	13-900037	F
Adapter for BECES 3.5/4.5 & KOC® For handgrip REF 311431	АНК	13-462319	D

WIRES FOR INTRA-ORAL WELDING




TITANIUM CAPS FOR LASER CONNECTION

Multi-use titanium caps for:

- use in immediately lasered bridge frames, together with the bar profiles (without bar matrices)
- the radiological control of plastic modeling
- for direct Polymerization into the bridge prosthesis
- · direct veneering with titanium ceramics
- material: Ti Grade 4



INTRA-ORAL WELDING



Nanda S., Ihde S., Nanda P. Intra-oral welding-A usefull adjunct in immediate loading implantology using BECES implants. CMF Impl. Dir. Vol 9, No.2, 13-24, 2014

SCANBODIES



CEMENTABLE ANGULATION ADAPTER (TI6AL4V)

These adapters are mounted on **BECES®** implants to compensate for the insertion direction. Plastic cements are preferably used. The implant head must be roughened beforehand. The protruding head parts are then removed. The impression is taken directly on the adapter.



CASTABLE CROWN BASE

These adapters are used by the dental technician for modeling of bridge frames. In the metal try-in, the protruding head parts are removed by the dentist.



LAB ANALOGUE



Description	Code	REF	Price cat.
Abutment analogue for angulation adapter 15° and 25°	AAA	13-462049	В

REF

13-462045

Price cat.

CASTABLE ABUTMENT AND IMPRESSION TRANSFER



Description	Code	REF	Price cat.
Castable abutment and transfer for AAA (pack of 5)	PA AAA	13-462050	А

CEMENTING ABUTMENT

Replacement abutment for cementing. For BECES implants up to a shaft diameter of 2.1 mm. Larger shafts must be ground down. Allows the vertical correction of the abutment position. Mounting e.g. with Fuji Plus. With drain hole, machined surface. Material **Ti6Al4V**.



Description	Code	REF	Price cat.
Replacement abutment for BECES,	B21	13-900316	Α
internal diameter 2.15 mm			

INSTRUMENT TRAY FOR KOC® & BECES®



Size of closed tray **B** 175 mm **T** 145 mm **H** 65 mm

For all autoclaves. Autoclaveable up to 134° C, not suitable for dry heat sterilizers.

Description	System	Head	REF	Description	System	REF	Price €
IT2 BECES	KOC/BECES	large	13-900030	Twist Drill 2.0 30	BECES	13-90020	
IT2 S BECES	KOC/BECES	large	13-900038	Twist Drill 2.0 21	BECES	13-90022	
IT2 W	KOC/BECES	large	13-900039	Twist Drill 2.5 21	BECES	13-90026	
IT K	KOC/BECES	small	13-462320	Twist Drill 1.8/23	BECES	13-90024	
ITS K	KOC/BECES	small	13-462322	BCD 1	KOC/BECES	13-900240	
ITW K	KOC/BECES	small	13-462331	BCD 2	KOC/BECES	13-900241	
ITWH K	KOC/BECES	small	13-462323	BCD 3	KOC/BECES	13-900242	
DOS 1	KOC		13-455311	BCDX 1	KOC/BECES	13-900243	
DOS 2	KOC		13-455312	BCDX 2	KOC/BECES	13-900244	
DOS 3	KOC		13-455313	BCDX 3	KOC/BECES	13-900245	
DOS 4	KOC		13-455314	CDG	KOC/BECES	13-420329	
DOS 5	KOC		13-455315	CDG	KOC/BECES	13-420329	
C-Drill KM 1	KOC		13-455300	DX 2	KOC/BECES	13-500704	
C-Drill KM 2	KOC		13-455301	TW 2	KOC/BECES	13-425402	
C-Drill KM 3	KOC		13-455302	Instrument tray empt	у	13-60006-K	upon request
IT LOC K	KOC		13-462333	Instrument tray with	content	13-S60006-K	upon request
DS 2	KOC		13-425001				
IT TB K	KOC		13-462327	The	ontent for the sys	tem KOC® is optio	nal

INSERTION TOOLS



STARTER TRAY

Autoclaveable up to 134° C, not suitable for dry heat sterilizers.

This surgical kit contains all drills and tools for first works with

the system BECES® and BECES® MU.

Material: autoclaveable plastic.



Description	REF	Price €
ITK	13-462320	
ITS K	13-462322	
IT 2 BECES	13-900030	
IT 2 S BECES	13-900038	
BCD1	13-900240	
Twist Drill 2.0 21	13-90022	
Twist Drill 2.0 30	13-90020	
Twist Drill 2.5 21	13-90026	
Twist Drill 1.8/23	13-90024	
BCDX 1	13-900243	
HT 1.25	13-425100	a s
ITX MU 15	13-418203	optional
Torque wrench TW2	13-425402	0 0
Starter tray empty	13-60040-K	upon request
Starter tray with content	13-S60040-K	upon request

SINGLE PIECE IMPLANT PRO KIT


All trays are delivered WITHOUT CONTENT. The tray offers a quick overview of the different lengths and diameters at hand, as well as the available amount of the corresponding implants.




ZDI ZYGOMA SCREW IMPLANTS

ZDI implants are placed trans-sinualy or below/inside the crista zygomatico-alveolaris, and in the alveolar crest of the maxilla. For cortical anchorage in the zygomatic bone. Depending on the anatomic situation the smooth parts of the implants are positioned below the Schneiderian mebrane or below the oral muscosa. These implants are to be used only by trained surgeons. ZDI implants provide a bending area below the cementing abutment. Therfore they can be aligned with other abutments in the center of the crest after the implant has been placed. An additional vertical osteotomy may be necessary (see the drawing). This implant can be used in combination with tubero-pterygoid screw implants (BECES®). The treatment protocol requires immediate splinting. **Material** Titanium alloy Ti6Al4V ELI.

a) max. diameter of thread	4.6 mm
b) enossal length	35 - 55 mm
c) abutment hight	7 mm
d) max. abutment-Ø	3.9 mm
f) upper Ø of shaft	2.0 mm
g) Ø of shaft	2.2 mm
h) length of thread	10 mm



Description	enossal Ø	Length	REF	Price cat.
Allfit ZDI 4.6 35	4.6	35	13-900100	F
Allfit ZDI 4.6 37.5	4.6	37.5	13-900101	F
Allfit ZDI 4.6 40	4.6	40	13-900102	F
Allfit ZDI 4.6 42.5	4.6	42.5	13-900103	F
Allfit ZDI 4.6 45	4.6	45	13-900104	F
Allfit ZDI 4.6 47.5	4.6	47.5	13-900105	F
Allfit ZDI 4.6 50	4.6	50	13-900106	F
Allfit ZDI 4.6 52.5	4.6	52.5	13-900107	F
Allfit ZDI 4.6 55	4.6	55	13-900108	F



 $ZDI \, Implants \, may \, be \, used \, in \, a \, trans-sinusal \, or \, sub-mucos al \, manner. \, The \, abut ment \, head \, is \, a ligned \, with \, the \, tooth \, arch \, through \, bending.$

TWIST DRILLS



Description	REF	Price cat.
Twist Drill 2.2 / 50 for Zygoma implants, SS	13-90021	D
Twist Drill 2.2 / 55 for Zygoma implants, SS	13-90023	D
Twist Drill 2.2 for handgrip to Zygoma implants. Length 100 mm	13-310514	F

ACCESSORIES FOR ZSI




Description	Einheit	Code	REF	Price cat.
Impression post castable, POM Internally round	Pack of 5	TSPA 5	13-462030	В
Double analogue metal		IA4/IAU	13-462112	Α
Double analogue plastic	Pack of 5	IA4/IAU	13-462111	В
Castable part	Pack of 5	POB	13-462086	В

ZDI MU ZYGOMA SCREW IMPLANTS

ZDI MU implants feature a pre-angulation of 15 degrees. **ZDI MU** may be bent additionally, using the insertion tool. In conjunction with the clinically possible rotational positions of the head, virtually all possible angulations can be realized. **ZDI MU** implants may be used by authorized users only. Material **Ti6Al4V**.

a) max. Ø of thread	4.6 mm
$\max.theoreticalthread\varnothing$	5.6 mm
b) enossal length	35 - 55 m
c) abutment hight	3.7 mm
d) platform \emptyset	4.8 mm
e) height of connecting part	2 mm
f) upper Ø of shaft	2.0 mm
g) Ø of shaft	2.2 mm
h) length of thread	10 mm
i) length of bending zone	10 mm



Description	enossal Ø	Length	REF	Price cat.
Allfit ZDI 4.6 35 MU	4.6	35	13-900110	N
Allfit ZDI 4.6 37.5 MU	4.6	37.5	13-900111	N
Allfit ZDI 4.6 40 MU	4.6	40	13-900112	N
Allfit ZDI 4.6 42.5 MU	4.6	42.5	13-900113	N
Allfit ZDI 4.6 45 MU	4.6	45	13-900114	N
Allfit ZDI 4.6 47.5 MU	4.6	47.5	13-900115	N
Allfit ZDI 4.6 50 MU	4.6	50	13-900116	N
Allfit ZDI 4.6 52.5 MU	4.6	52.5	13-900117	N
Allfit ZDI 4.6 55 MU	4.6	55	13-900118	N

ACCESSORIES FOR ZDI MU



REPROCESSING OF TOOLS AND DRILLS



MANUFACTURER'S INFORMATION regarding the prepagration of resterilisable medical devices complies with PIS 17664 Please read carefully! Medical devices which may be re-processed are Medical devices which may be re-processed are

Please read carefully!

Medical devices which may be re-processed are 1 tools for abutments and screws 1 tools for abutments for preparing endosseous bone cavilties (drills, cutters) 1 tools sleeves 1 tools for abutments for preparing endosseous bone cavilties (drills, cutters) 1 tools sleeves 1 tools for abutments for preparing endosseous bone cavilties (drills, cutters) 1 tools sleeves 1 tools for abutments and screws 2 tools for abutments and screws 2 tools for abutments and papintments and are not used on other patients. They should be stored by the operator between the treatment appointments, e.g. together with the and bone preparation.

Re-usability 1 frequent re-processing has influence on the product and the products are basically suitable for the cleaning and disinfection of instruments and bone preparation.

Re-usability 1 frequent re-processing has influence on the products are basically suitable for the cleaning and disinfection of instruments and bone preparation.

Re-usability 1 frequent re-processing has influence on the products are basically suitable for the cleaning and disinfection agent with proven efficacy (e.g. DCHM of FDA approved and CE Mark) 1 that the chemical sused are compatible with the instruments and reaction times given by the manufacturer is work and in the case of the reaction times given by the manufacturer is voice, if these restictions are not regarded.

Legal bases

The following legal bases, regulations and recommended to the products are basically suitable for ultrasonic cleaning of the instruments showing signs of wear must be discarded. Liability of the manufacturer is voice, if these restictions are not regarded.

Legal bases

The following legal bases, regulations and recommended to the products are the produ

walld for advising patients before and after the placement of the implants.

General principles

Altreusable products must be cleaned, disinfected and
sterilised before each use, this also applies to the Initial
sterilised before each use, this also applies to the Initial
sterilised before each use, this also applies to the Initial
sterilised in the Initial sterilised in the Initial
sterilised in Special cleaning fisterilised in Instructions
should be obtained from the instructions for use. The
operating instructions of the practice units must also be
operating instructions of the practice units must also be
operating instructions of the practice units must also be
operating instructions of the practice units must also be
operating instructions of the practice units must also be
operating instructions of the practice units must also
portating instructions for use. The
operating instructions for use, the
operating instruction of the unit and
product are constantly maintained during each cycle.
Please also observe all valid legal and hygiene regulations of the dental practice and dental hospital.

Instruments observe all valid legal and hygiene regulaloss of the dental practice and dental hospital.

Instruments made from different materials should
never be distincted, cleaned or sterilised.

During mechanical cleaning, instruments should be
arranged so that they cannot come into contact, as
otherwise there is the risk of damage.

Multi-part instruments such as raticests, trephine
their component parts and these should be individually disinfected, cleaned or sterilised.

These instruments should also be stored disassembled until the next use.

Care instructions of surgical steel instruments

Care instructions of surgical steel instruments
Surgical steel instruments can quickly become damaged with inadequate or incorrect care. Only commercally available solvents should be used for surgical
steel if in doubt contact Simpladent GmbH.
The following are not recommended:

- Disinfection/cleaning agent with a high chlorine
content

- content of the conten

• loo high temperatures with mechanical cleaning and sterilisation, never higher than 135°C.
Conditioning:
Corrie imputities must be removed from the products immediately after use (within 1:2 hrs maximum). Surgical residue (blood, secretion, itsue residue) should stand the product of the product of the standard stand

SIMPLADENT®

- repeat the cycle or clean manually.

 Manual cleaning

 1. Thoroughly clean disinfaction/cleaning agent
 1. Thoroughly clean disinfaction/cleaning agent
 1. Thoroughly clean disinfaction them with water
 and. I required, with the aid of a soft nylon brush

 Ultrosonic cleaner. Place the components in a basket, avoid acoustic shadows. Add an enzymatic
 cleaning agent to the water and clean the components of a temperature of 40 50° C in the ultrasonic
 components are immersed completely in the water
 without bubbles.

 2. Then remove the instruments from the cleaning solution and rinse them thoroughly (minimum 1 min.)
 under running water. Use fully desclinded water for
 1. Then day the instruments with compressed air
 4. Check the instruments would and repeat the cleaning stage, if necessary.
 5. Pack the instruments would van drepeat the cleaning
 stage, if necessary.
 5. Pack the instruments would van drepeat the cleaning
 and a clean focation.
 6. Document the approval.

 Mechanical cleaning

Mechanical cleaning
Cleaning, disinfection and drying in accordance with
DIN ENISO 1888-1 2006 and DIN EN 15883-2006
Per-Liceaning, Place the disassembled instruments with
DIN ENISO 1888-1 2006 and DIN EN 15883-2006
Instruments with a soft nylon brush under water to remove coarse impurities.
Mechanical Celaning, e.g., using the Miele 8535 CD
unit of 15° C for 3 minutes (programme Vario 1D) with
an enzymatic cleaner.

- Important points

 All instruments must be sterilised after cleaning.

 When sterilising multi-part instruments in an autoclass.

 When sterilising multi-part instruments in an autoclassing multi-part part of the instruments are always sterilised in a disassembled state!

 The instruments should always be checked for corrosion after sterilisation. Instruments must still be visible after sterilisation: otherwise the instruments should be replaced.

- after sterilization of the standard standard should be replaced. New instruments must be cleaned and sterilised without peckediging before using for the first firme. Preparation of all instruments with cavities is particularly critical. This applies especially to internally cooled drills, placement aids and instruments with cavities in the properties of the standard sterilization of the standard standard

Control
Check all instruments after cleaning and cleaning/dis-infection for corrosion, damaged surfaces, chipping, damage to the shape (e.g. bent and non-concentiumning instruments, damaged or blunt blades) as well as contamination and discard any damaged instruments, instruments that are still contaminated must be cleaned and dislinected again. Then check the function and comply care products (e.g. oil) to instruments and abutments or screws.

Special aspects to observe with drills and cutters Use cutting instruments for a maximum of 10 times. Thoroughly check these instruments after each use for cleanlines (including the internal cooling sections in of bone drills depends on the hardness of the bone at the site. If in doubt, drills should only be used once. There is a considerable loss of cutting performance if the tip is damaged. Io ensure care of the drills it is the-refore essential to observe the following points: During the operation drills should be placed partly because the properties of the properties of the properties of logical saline solution for longer than 1 hour to avoid corrosion. Never drop the drills directly on the tip The drills should not come into contact during ultra-sonic cleanling

Packaging
Sort out the instruments in the sterilisation tray and
then pack them in single-use sterilisation packaging
(single or double packaging) and/or sterilisation containar which.

Legend

STERILE LOT

A

Read instructions

Expiration date

Only use once Do not resterilize

non sterile

LOT Charge number Keep in a dry place

Store tightly keep closed Do not use if packing is damaged

Manufacturer

C €1936

STERILE R Y Gamma-sterilized

- iner, which complies with DIN EN 868-2ff/DIN EN ISO/ANSI AAMI ISO 11607
- ISO 11407

 Is valiable for steam starilisation (temperature resistant up to min. 137° C (279° F), adequate steam permeability)

 provides adequate protection of the instruments and stafilisation packaging against mechanical damage

 is regularly serviced according to the manufacturer's instructions (sterilisation container)

Sterilisation
Method:
Fractional pre-vacuum procedure
[according to ISO 17665 or ISO 13060]
In a unit that complies with EN 285
Forsperature:
Feel to 1302 (-2 mox, 137" C
Free to 150 (-2 mox, 137" C

Storage
After distillication, the instruments must be stored dry
After distillication, the instruments must be stored dry
After distillication and the production packaging. The instruments should also be protected against sumillight and
heat. The maximum storage period (expiry date) depends on several factors and must be determined and
validated by the user.

Information on handling multi-part instruments
Multi-part instruments must be disassembled before
sterilisation. Please note the schematic diagram be-

Multi-part instruments as sterilisation. Please note the schematic diagram beterilisation. Please note the schematic diagram beterilisation of the schematic diagram bepenerally list before the beginning of implant placement, the ratchet should be thinly lubricated using a silicence oil and reassembled. The function of the ratchet
should then be checked before beginning surgery.

Warnings
We do not know of any warnings, provided the instruc-tions for use are followed for the products to be used as well as the corresponding disinfection and cleaning agent.


Simpladent GmbH reserves the right to change the design of the products and components of their pa-ckaging, adapt instructions for use as well as renego-tiate prices and delivery conditions. Liability is limited to the use of defective products. Any further claims are excluded.

Further information about the preparation of medical products is available in the Internet at www.rki.de or www.a-k-i.org.

Date of the latest revision: 2018-10

Schematic diagram of the TW/TW2 torque wrench

After use the instrument should be disassembled into its individual parts – no tool is required for disassembly



Pre-clean the individual parts under running cold water using a soft brush. Do not allow blood residue and other adhering deposits to dry on the components.


Schematic diagram of the RAT2 ratchet

After use the instrument should be disassembled into its individual parts – no tool is required for disassembly

Pre-clean the individual parts under running cold water using a soft brush. Do not allow blood residue and other adhering deposits to dry on the components. The ratchet should be autoclaved in the disassembled immediately before use.

Schematic diagram of the handle REF 311430 (can be disassembled)

After use the instrument should be disassembled into its individual parts – no tool is required for disassembly



Pre-clean the individual parts under running cold water using a soft brush. Do not allow blood residue and other adhering deposits to dry on the components. The handle should be autoclaved in the disassembled mendiately before use.

Schematic diagram of the handle REF 311431 (cannot be disassembled)



- Pre-clean the instrument under running cold water using a soft brush. Do not allow blood residue and other adhering deposits to dry on the handle. The handle should be thoroughly cleaned manually using an ultrasonic cleaner before mechanical cleaning. Manual cleaning including ultrasonic cleaner (see above) and mechanical cleaning should be performed in sequence.



SIMPLADENT GmbH Dustral CH-8737 Gontine SWITZERLAND Tel. 411 (0)55 293 23 33 Fax +41 (0)55 293 23 30 simpladente@mplant.com www.simpladent-implants.com 1-8737 Gommiswald I-17ERLAND

Distributed by

SIMPLADENT implant Solutions Private Limited Shop. No. 110 & 111 F.F. Vardhman Geedee Plaza, Plot No. 10 Sector-12. Dwarka, New Delhi-110075 Tel. 101 457 875 email: omsaidental@yahoo.com



(The products of this catalogue are CE marked (class I) and CE 1936 marked (class IIa and IIb) according to 93/42/EC Directive).

Commercial products that are not monitored by our notified body are declared as third-party products.

Basal implants may only be used and operated by qualified persons with valid authorisation (para. 2 MedProdAnw Verordnung). We are certified according to DIN EN ISO 13485 and Annex II of Directive 93/42 EEC.

The product dimensions shown in this brochure may differ from reality for technical reasons.

 \mathbf{KOC}^{\bullet} and \mathbf{BECES}^{\bullet} are registered trademarks. Pat. Pend.

If implants are reprocessed, there is a risk of the development of infections, because no validated method for processing exists. Implants therefore may not be reprocessed.

Compilation and explanation of symbols on the packaging:



Batch No.



Sterilized by radiation



Non-sterile



Intended for use by dentists or surgeons only



Single use product



Instruction for use



Expiry date




















Simpladent GmbH

Dorfplatz 11 SWITZERLAND

Phone +41 (0)55 293 23 70 implants@simpladent.ch CH - 8737 Gommiswald / SG www.simpladent-implants.com

PLEASE ASK YOUR LOCAL DISTRIBUTOR FOR THE VALID PRICE LIST

SIMPLADENT®