SYSTEM APPLICATION



IMMEDIATE LOADING DENTAL IMPLANT SYSTEM

SIMPLADENT[®]



Simpladent GmbH

Dorfplatz 11 CH - 8737 Gommiswald / SG SWITZERLAND Phone +41 (0)55 293 23 70 implants@simpladent.ch www.simpladent-implants.com

IMPLANT SYSTEM FOR ENDOSSEOUS DENTAL IMPLANTIONS



Secure rotation protection by precision inner octagon and internal cone 8°

Reliable insertion – ease of use

Universal application for permanent and removable prosthodontics

Satisfied patients in each phase of the procedure

The Allfit[®] S implantsystem **SSO[®]**, **STO[®]** and **STW[®]** is for two-stage endosseous dental implants and consist of highly break-proof titanium alloy Ti6Al4V ELI.

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).

Due to technical reasons the product dimensions shown in this brochure might deviate from reality.


 $\mathsf{SSO}^\circ,\mathsf{STO}^\circ$ and STW° are registered trademarks.

In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.



SURGERY

1.1a. Recommended Drill Sequences for \textbf{SSO}^{\otimes} Implants




1.1b. Recommended Drill Sequences for STO® Implants



1.2 Tap for lower jaw (only for use with SSO[®])

TAP 1Size 3.3 (for Implants Ø 3.3 mm)Size 4.1 (for Implants Ø 4.1 mm)Size 4.8 (for Implants Ø 4.8 mm)



2. Implant packaging SSO* / STO $^{\circ}$



Original packaging

13-420021 Allfit STO 3.7 11	Implant Ti6AI4V
LOT D1668730411AS	30.11.2018
	STERLER
	entert .

Open the pack using the flap. Remove the label and stick it into your patient's record.



The blister (secondary packaging) contains the implant in a sterile tube (primary package).

3. Removing the implant from the sterile tube

- 1. Open the lid. The implant is connected to the lid through a breakable section.
- 2. Remove the implant without touching the inner walls of the tube.
- 3. For some types the implant comes with an insertion tool already pre-mounted.



6

4. Handling

Assembling the placement aid

Attach the insertion tool to the implant (SSO or STO) by holding the Lid, to which the implant is secured, with your other hand.



After you have attached the insertion tool, firmly hold the lid in your hand and break the implant off the lid. Insert the implant manually until it sits tightly in the jaw bone.



5. Insertion

Using the ratchet or contra-angle handpiece: screw the implant clockwise into the bone cavity. The endosseous part of the implant must be **completely** covered by the bone. The polished implant neck is located partially outside of the bone or at bone level. We recommend to screw min. 1 mm of the polished implant neck inside of the bone.



6. Remove Insertion Tool from Implant

Release of the insertion tool or the contra-angle handpiece from the implant:

Using HT 1.25 loosen the screw in the insertion tool while loading the insertion tool IT STO in insertion direction with the ratchet.

Pull insertion tool off the implant and separate from the contra-angle hand piece.

Attach IT ITV Ratchet adapter with RAT2 ratchet to the ITV insertion tool. Using the HAS flat wrench firmly hold the lower hexagon of the ITV. Remove the ITV from the implant using the ratchet ("Out" position).



6. Result

The implant is fully inserted and ready for healing.



8. Post-operative Treatment

Close the implant with a cover screw CST STI (tool: HT 1.25) HT 1.25 HT 1.25 HT 1.25 HT 1.25 HT 1.25 HT 1.25 TT CST STI CST ST



9. Impression taking

9.1 Impression taking with perforated, individual impression tray The long pick up screw must be clearly visible, when

the impression try is inserted over the impression post.

The impression material must be removed in the area of the screw access prior to hardening.

9.2 Prior to impression taking Impression taking with an A silicone® such as Safeprint®.

The use of open or closed impression trys is possible.


Hex-Instrument HT 1.25

Fasten impression post OLT STO

SSO[®]/STO[®] Implant

Arrow: View from above. Turn clockwise





9.3 Impression taking

Remove OLT STO from the implant. OLT STO remains in the impression.

Loosen screw with HT 1.25

Window in Impression tray

OLT STO

SSO[®]/STO[®] Implant

9.4

9.5


to the laboratory.

View of the impression post OLT STO in the impression (Pick up method)

After the impression is taken, the implant is closed with a healing screw (HS) and the impression is sent

Position of the Impression post OLT STO

Impression material




SSO[®] / STO[®] Implant

Close healing screw clockwise

HT 1.25

10. Further processing n the Laboratory



10.5

Ensure proper position of the abutment when transferring into the mouth.Tightening torque of the screw during fastening on the implant with: 25 Ncm.

TLA 15 STO



SIMPLADENT®

Gipsum

SSO[®] IMPLANTS

Successful implant with octagon socket and 8° inner cone. Octagonal system. Connection type: Inner cone + internal octagon / head-Ø: 4.8 mm. Made from titanium alloy ASTM F 136-13/ISO 5832- 3. For rotationally-secured supra-structures. Surface is roughened in endosseous region.

		Description SSO 3.3 9	endosseous Ø 3.3 mm	endosseous length 9 mm	REF 13-420000	Price cat. F
a		SSO 3.3 11	3.3 mm	11 mm	13-420001	F
b	Internal octagon	SSO 3.3 13	3.3 mm	13 mm	13-420002	F
		SSO 3.3 15	3.3 mm	15 mm	13-420003	F
		SSO 4.1 7	4.1 mm	7 mm	13-420004	F
c I		SSO 4.1 9	4.1 mm	9 mm	13-420005	F
a) aesthetic head	1.8 mm	SSO 4.1 11	4.1 mm	11 mm	13-420006	F
b) endosseous length	7 - 17 mm	SSO 4.1 13	4.1 mm	13 mm	13-420007	F
c) endosseous Ø	3.3 -4.8 mm	SSO 4.1 15	4.1 mm	15 mm	13-420008	F
d) Head Ø	4.8 mm	SSO 4.1 17	4.1 mm	17 mm	13-420009	F
		SSO 4.8 7	4.8 mm	7 mm	13-420010	F
		SSO 4.8 9	4.8 mm	9 mm	13-420011	F
		SSO 4.8 11	4.8 mm	11 mm	13-420012	F
		SSO 4.8 13	4.8 mm	13 mm	13-420013	F



Delivery incl. surgical screw CST STI

Note in contraindications:

Implants made from c.p. titanium are in general less force resistant compared to those made from titanium alloy. Implants having a nominal diameter of less than 3.8 mm are not indicated for single tooth replacement, independently if they are made from c.p. titanium or from titanium alloy. Implants Typ STO 3.7 and SSO 3.3 should be used as supporting implants.

STO 3.7 mm and SSO 3.3 mm implants may not be used in the molar region and not as single tooth implants with off axis load. STO 3.7 mm and SSO 3.3 mm implants are used as supporting implants, for example to increase the number of abutments in immediate load situation. Never use STO 3.7 mm and SSO 3.3 mm implants when cantilevers are involved. STO 3.7 and SSO 3.3 are not for single tooth replacement.

Note the following contra-indications:

STO 3.7, SSO 3.3: Do not use in areas where off-axis load is present. Do not use in areas with strong chewing forces. Do not use in flexion areas of the jaws.

STO® IMPLANT

Precision conical implant with octagon socket and 8° inner cone. Octagonal system. Connection type: Inner cone + internal octagon / head-Ø: 4.8 mm. Strengthened titanium alloy (Ti6AI4V ELI) according to ASTM F 136/ISO 5832-3. Compression screw thread and 8° internal cone for stable endosseous anchorage.





ACCESSORIES FOR SSO®, STO®



Description Cover screw Gingivaformer


	Code	REF	Price cat.
	CST STI	13-420300	А
for 2 mm gingival height, cylindrical	HS2 STI	13-420301	В
for 4 mm gingival height, cylindrical	HS4 STI	13-420302	В
for 5 mm gingival height, cylindrical	HS5 STI	13-420303	В

ABUTMENTS FOR CEMENTED PROSTHETICS FOR SSO®, STO®

Screw SF OLT, short, 14.5 mml for OLT STO. 5 pieces / pack

With Octagon. Anti-rotation, screw-in abutment with pre-assembled screw for cemented crowns or bridges. Fasten with TT 1.25.

	Description	Code	REF	Price cat.
	Straight	TLA STO	13-420306	D
\ \	15° angled. The angulation 15° runs a xial over the edge of the octagon.	TLA 15 STO	13-420307	E
	15° angled. The angulation of 15° runs axial over the area of the octagon (color: metallic yellow).	TLA 15F STO	13-420308	E
/	20° angled. The angulation of 20° runs $axial over the edge of the octagon$	TLA 20 STO	13-420309	E
	25° angled. The angulation of 20° runs axial over the edge of the octagon.	TLA 25 STO	13-420311	D
	25° angled. The angulation of 25° runs axial over the area of the octagon (color: metallic yellow).	TLA 25F STO	13-420312	D



There are 16 possible abutment positions.

	Description	Impression post with Pick-Up Screw	Laboratory analog (Inner octa)	Synthetic base, cauterizable, 7.4 mm high TLA STO, TLA STI	Synthetic base, cauterizable, 4 mm high TLA 15 STO	Synthetic base, cauter- izable, 4 mm high TLA 20 STO	Synthetic base, cauterizable, 4 mm high TLA 25 STO
	Code	OLT STO	IA STO	PS	PS 15	PS 20	PS 2 5
	REF	13-420350	13-420348	13-420354	13-420355	13-420356	13-420357
	Price cat.	В	В	С	A	A	А
ALTERNATIVE	Description			Code	REF	Price cat.	

SF OLT

13-420909

В

12

With cone (without anti-rotation). One piece abumtent with solid thread.

	Description	Code	REF	Price cat.
P	For cemented crowns and bridges. height above implant 6.5 mm, 8° incline. Can be shortened or trimmed. Flat on one side. Direct impression or transfer to laboratory analog IA STI. Fasten with HT 1.25	TCA STI	13-420313	D

Description	Impression post 9 mm high	Impression post 13 mm high	Laboratory analog with cone	Plastic base, cauterizable, 7.4 mm high TLA 20 STO	Plastic base, cauterizable, 4 mm high TLA 15 STO
Code	TS STI	TSX STI	IA STI	PS	PS 15
REF	13-420345	13-420411	13-420347	13-420354	13-420355
Price cat.	В	В	В	A	А

	Description	Height above implant	Color	Code	REF	Price cat.
	Flat on one side for increased rotation stability	4 mm	yellow	SA4 STI	13-420318	С
	6° incline, fasten with IT TCA	5.5 mm	gray	SA5 STI	13-420319	С
		7 mm	blue	SA7 STI	13-420320	С
116						

А

Description

Code

REF

Price cat.

Color coded transfer post Abutment analog for SA-Abutments use withTZ SA 4, 5, 7 Use with **REF 13-460230**

TZ SA 4, 5, 7

13-420406

А

13-420404, 13-420405,

AA SA 4, 5, 7 13-420324, 13-420325, 13-420326

В

Castable abutment 10 mm high, reducible

Anti-rotation (PA SA) Internally round (PA SR)

PA SA / PA SR 13-420327, 13-420328 Patented gingiva retraktor 4 pieces / pack (Use instead of retraction cord), Material: PP

13-460230

А





With cone connection (without anti-rotation). One-piece abutments with fixed thread and vertical cement escape grooves.

ABUTMENTS FOR SCREW-IN PROSTHETICS

		Description			Code	REF	Price cat.
m		With Octagon (Anti-rotation). Abutment with double threaded screw bases. Anti-rotation Anchorage of the p height above implant 1.5 mm. Fasten wi	rosthetics.	dges and	OSA STO	13-420337	E
P		Single-piece and anti-rotational octago crowns, bridges and bar connectors. height above implant: 1.5 mm. Fasten with HT 1.77 / HTX 1.77	nal screw-in abutment;	for	OSA STI	13-420338	D
		or				T	or
Description	Transfer post	Screws for anchorage of TST STI	White castable abutment 10 mm high, edged internally 5 pieces / pack	White cas 10 mm hig internally 5 pieces /	round	Screws	
Code	TST STI	SF 365 or SF 415	PAOA	PAOR		SF 350 or SF	365
REF	13-420339	13-420938 or 13-420928	13-420342	13-42034	3	13-420930	or 13-420938
Price cat.	D	В	В	В		В	

ABUTMENTS FOR REMOVABLE PROSTHETICS

Description		Code	REF	Price cat.
TCT SSO 0.5		TCT SSO 0.5	13-420420	D
TCT SSO 1.5	for screw-in	TCT SSO 1.5	13-420421	D
TCT SSO 2.5	prosthetics Fasten with HT 1.77	TCT SSO 2.5	13-420422	D
TCT SSO 4		TCT SSO 4	13-420423	D

IMPRESSION AND LABORATORY ACCESSORIES

Description	Code	REF	Price cat.
Transfer post	TST	13-418147	В
long screw	SFL	13-420428	A
TCT-Analog	BTT	13-418100	В
Castable abutment, 12 mm high, internally round. 5 pieces / pack	PSTR (gray)	13-418124	В
Castable abutment, 12 mm high, internally edged. 5 pieces / pack	PSTA	13-418123	В
Screw for fastening	SF	13-418151	В

LOCALICER®

Abutment for removable prosthetics. Fasten with ${\rm HT~1.77}$

R



Description	Height	REF	Price cat.
LOC S 2.5	2.5 mm	13-420386	D
LOC S 4	4.0 mm	13-420385	D

ACCESSORIES FOR LOCALICER°

	Description	Code	REF	Price category
	Analog + impression cap Set	AA LOC	13-462337	С
	Set with 5 Caps + 1 Housing (EXTERNAL PRODUCT)	NCS	13-462338	D
20	Pull off force Yellow 600 g, Pink 1.200 g, Clear 1.800 g, Violet 2.700 g.			

Yellow 600 g, Pink 1.200 g, Clear 1.800 g, Violet 2.700 g Black has no retention and is designed for temporary solutions for up to one month.



TITANIUM BASE FOR CAD CAM

Base for Octa system / SSO. Material Ti6Al4V ELI, anti-rotation. Incl. screw SF N62.

	Description	Code	REF	Price cat.
	Height 3.9 mm above implant shoulder	MB 4 SSO	13-463120	В
ð	Height 8 mm above implant shoulder	MB 8 SSO	13-463122	В

RETENTIVE BALL ATTACHMENTS

Code

REF

Description	Code	REF	Price cat.
Ball abutment fits nylon cap NC / NC1 / NC2 Head Ø 2.5 mm, height above implant 3 mm. Fasten with HT 1.25	TB STI	13-420361	D
Ball abutment fits nylon cap NC / NC 1 / NC 2 Head Ø 2.5 mm, height above implant 5.4 mm. Fasten with Tool E (alternatively with HAS2).	TB3 STI	13-420363	D



Description	 Laboratory analog with cone Transfer post, fasten with HT 1.25 Ball analog 	Nylon caps NC 2 pieces / pack (EXTERNAL PRODUCT)	Nylon caps R-NC with increased friction strength 2 pieces / pack (EXTERNAL PRODUCT)	Sleeve for all NC (EXTERNAL PRODUCT)
Code	IA STI, TS STI, IAB	NC, NC 1, NC 2	R-NC, R-NC 1, R-NC 2	Н
REF	13-420347, 13-420345, 13-420566	13-465028, 13-465029, 13-465030	13-465034, 13-465033, 13-465032	13-465031
Price cat.	В, В, С	A1	A1	В




STW® IMPLANTS

Octagonal system

Implant with Internal octagon and 8° internal cone. For anti-rotation Superstructures. Highly fracture resistant titanium alloy (Ti6Al4V ELI) according to ASTM F 136 / ISO 5832-3.

These implants with a **head diameter of 6.5 mm** can be used, if adequate bone is available, for posterior single crowns with the width of one premolar in dentate and partially/fully edentulous ridges as for superstructures retained by bars and ball attachments.



		Description	endosseous Ø	endoss. Length	REF	Price cat.
(STW 4.8 7	4.8 mm	7 mm	13-420015	F
		STW 4.8 9	4.8 mm	9 mm	13-420016	F
lr	iternal octagon	STW 4.8 11	4.8 mm	11 mm	13-420017	F
		STW 4.8 13	4.8 mm	13 mm	13-420018	F

a) aesthetic head	1.8 mm
b) endosseous length	7 - 13 mm
c) endosseous Ø	4.8 mm
d) Head Ø	6.5 mm

SURGICAL ACCESSORIES



Description		Code	REF	Price cat.
Surgical Screw Cover screw for S	TW® implant, non sterile	CST STW	13-420370	В
Gingivaformer	for 2 mm Gingival height cylindrical	HS2 STW	13-420371	В
	for 4 mm Gingival height cylindrical	HS4 STW	13-420372	В
,	tem / STW, Ti6Al4V ELI, ight 3.9 mm above implant shoulder. 2	MB4 STW	13-463121	В

ABUTMENTS FOR CEMENTED PROSTHETICS FOR STW®

LOCK ABUTMENT for SSTW[®] with Octagon. Anti-rotation screw- in abutment with pre-mounted screw, for cemented crowns and bridges. Fasten with **TT 1.25**

	Ų	Description Straight abutment 15° angledes abutment. The of the octagon.	15° line angle runs axially acros		Code TLA STW TLA 15 STW	REF 13-420373 13-420374	Price cat. D E
Description		Impression post with Pick-Up screw Fasten with HT 1.25	Laboratory analog	Plastic bas	se, castable FW	Plastic base, ca for TLA 15 STW	
Code		OLT STW	IA STW	PA STW		PA 15 STW	
REF		13-420351	13-420349	13-42035	9	13-420360	
Price cat.		В	В	А		А	
	ALTERNATIVE	Description Screw SF OLT, short, 14.5 mm 5 pieces / pack	ll for OLT STW		Code SF OLT	REF 13-420909	Price cat. B
		Replacement screw for TLA a	butments of all S-implants	:	SF TLA	13-420910	A

SOLID ABUTMENTS

Description	height above implant	Code	REF	Price cat.
6° anlge, fasten with HT 1.77 Flattened on one side for increased rotationa	4 mm Il	SA4 STW	13-420375	С
stability	6 mm	SA6 STW	13-420376	С



|--|

Code REF Price cat. Transfer post for SA6 abutment

TZ SA6 STW 13-420377

А

Abutment analog brown anodized, for **TZ SA6 STW**

AA SA6 STW 13-420378

А

13-420359 A

PA STW


Castable abutment

RETENTIVE BALL ATTACHMENTS

	Description	Code	REF	Price cat.
*	Ball abutment, fits nylon cap NC / NC1 / NC2 Head Ø 2.5 mm, height above implant 3 mm. Fasten with HT 1.25	TB STW	13-420365	D
1.2.3.	Ball abutment, fits nylon cap NC Head Ø 2.5 mm, height above implant 3.95 mm. Fasten with Tool E (alternatively with HAS2 möglich).	TB3 STW	13-420364	D

Description	1. Laboratory analog 2. Transfer post, fasten with HT 1.25 3. Ball analog	Nylon caps NC 2 pieces / pack (EXTERNAL PRODUCT)	Nylon caps R-NC with increased friction strength 2 pieces / pack (EXTERNAL PRODUCT)	Sleeve for all NC (EXTERNAL PRODUCT)
Code	IA STW, TS STW, IAB	NC, NC 1, NC 2	R-NC, R-NC 1, R-NC 2	н
REF	13-420349, 13-420346, 13-420566	13-465028, 13-465029, 13-465030	13-465034, 13-465033, 13-465032	13-465031
Price cat.	B, B, C	A1	A1	В

DRILL SEQUENCES CYLINDRICAL IMPLANTS



DRILL SEQUENCES CONICAL IMPLANTS

Implants	endoss.Ø	Length in mm	recommended	d drill sequences					
		9 / 11 / 13	DS 2 13-425001	DS 2.2 →	DRK 1 13-425020	DRA 3.3 CS 13-425046			
STO [®]	3.7 mm		DSL 2 13-425002 →	DSL 2.2 13-425004 +	DRK 1 13-425020	DRA 3.3 CS 13-425046			
510	5.7 11111	15	DS 2 13-425001 →	DS 2.2 +	DRK 2 13-425021	DRA 3.3 CS 13-425046			
			DSL 2 13-425002 →	DSL 2.2 13-425004	DRK 2 13-425021	DRA 3.3 CS 13-425046			
		7/9/11/13	DS 2 13-425001 →	DS 2.2 +	DRK 3	DRA 4.1 CS 13-425047			
		7/9/11/13	DSL 2 13-425002 →	DSL 2.2 13-425004	DRK 3	DRA 4.1 CS 13-425047			
STO [®]	4.1 mm	1.1 mm 15 / 17 / 19 21 / 23 / 25	DS 2 13-425001 →	DS 2.2 +	DRK 2 13-425021	DRA 4.1 CS 13-425047			
510	7.1 11111		DSL 2 13-425002 →	DSL 2.2 13-425004 +	DRK 2 13-425021	DRA 4.1 CS 13-425047			
			DS 2 13-425001 →	DS 2.2 +	DRK 6 13-425025	DRA 4.1 CS 13-425047			
			DSL 2 13-425002 →	DSL 2.2 13-425004 +	DRK 6 13-425025	DRA 4.1 CS 13-425047			
		7 / 0 / 11	DS 2 13-425001 →	DS 2.2 +	DRK 4 13-425023	DRA 4.1 CS 13-425047			
STO [®]	4.8 mm	7/9/11	DSL 2 13-425002 →	DSL 2.2 13-425004 +	DRK 4 13-425023	DRA 4.1 CS 13-425047			
	4.011111	12 / 15	DS 2 13-425001 →	DS 2.2 →	DRK 5 13-425024	DRA 4.1 CS 13-425047			
		13 / 15	DSL 2 13-425002 →	DSL 2.2 13-425004 +	DRK 5 13-425024	DRA 4.1 CS 13-425047			
STW [®]		4.1 mm all lengths	DS 2 13-425001 →	DS 2.2 →	DS 2.8 13-425005	DS 3.5 13-425009	DS 4.0 13-425011 →	(DSL 4.2) (13-425014)	DRA 4.8 CSW 13-425048
	-, I IIIII		DSL 2 13-425002 →	DSL 2.2 13-425004	DSL 2.8 13-425006	DSL 3.5 13-425010	DSL 4.0 13-425012 +	(DSL 4.2) (13-425014)	DRA 4.8 CSW 13-425048

Note that STO is by design a compression screw and a drill sequence can not be recommended for all bone qualities. In low density bone a smaller drill is recommended.

PATHFINDER DRILLS

	Description	Туре	Working length	REF	Price cat.
	BCD 1	non sterile	15 mm	13-900240	С
13036	BCDX 1	non sterile	15 mm	13-900243	с

TWIST DRILL FOR SSO®, STW®

Surgical steel, laser marked, externally irrigated. Non sterile. Maximum speed with external cooling 700 RPM. For lifetime extension and consistent cutting performance and easy removal of contaminated areas. Drills can be used without cooling at 50 RPM on green or blue contra-angle handpieces.

These drills are made of premium material using state-of-the-art precision milling machines. The benefit for you and your patients: Extremely good and durable cutting performance due to sophisticated geometry of the blades. Almost no heat development and excellent as a result of this geometry.

	18 mm	Description	Length	Working length	Drill Ø	REF	Price cat.
	 1 9 7	DS 2	36.5 mm	18 mm	2 mm	13-425001	D
08 2.2	620	DS 2.2	36.5 mm	18 mm	2.2 mm	13-425003	D
	1811	DS 2.8	36.5 mm	18 mm	2.8 mm	13-425005	D
05 3.2		DS 3.2	36.5 mm	18 mm	3.2 mm	13-425007	D
		DS 3.5	36.5 mm	18 mm	3.5 mm	13-425009	D
		DS 4.0	36.5 mm	18 mm	4.0 mm	13-425011	D
		DS 4.2	36.5 mm	18 mm	4.2 mm	13-425013	D
	mm	DSL 2	45.5 mm	27 mm	2 mm	13-425002	D
		DSL 2.2	45.5 mm	27 mm	2.2 mm	13-425004	D
		DSL 2.8	45.5 mm	27 mm	2.8 mm	13-425006	D
M		DSL 3.2	45.5 mm	27 mm	3.2 mm	13-425008	D
		DSL 3.5	45.5 mm	27 mm	3.5 mm	13-425010	D
		DSL 4.0	45.5 mm	27 mm	4.0 mm	13-425012	D
		DSL 4.2	45.5 mm	27 mm	4.2 mm	13-425014	D

HEATLESS® DRILL "DRK" FOR IMPLANTS WITH CONICAL CORE

Form drills made from surgical steel (SS), length and color coded, resterilizable. Maximum motor speed with external cooling 700 rpm. Note the information in the instructions for use for hygiene and preparation requirements with medical products intended for multiple use in implantology. Pilot drilling using DS2 (P.33), DSL2 (P.33) or Pathfinder (P.34)

Description	System	Color	Length	Implant Ø	REF	Price cat.
DRK 1	STO [®]	yellow	9, 11, 13 mm	3.3 / 3.7 mm	13-425020	D
DRK 2	STO [®]	black	15, 17, 19 mm	3.3 / 3.7 / 4.1 mm	13-425021	D
DRK 3	STO [®]	red	7, 9, 11, 13 mm	4.1 mm	13-425022	D
DRK 4	STO [®]	blue	7, 9, 11 mm	4.8 mm	13-425023	D
DRK 5	STO®	green	13, 15 mm	4.8 mm	13-425024	D
DRK 6	STO®	metallic	21, 23, 25 mm	3.7 / 4.1 mm	13-425025	D

It has been scientifically proven that heatless drills generate 55% less heat than traditional bone drills from other manufacturers. This makes it possible to use higher rotational speeds: between 3,000 and 5,000 rpm are recommended with good external cooling and intermittent drill technique.

COUNTERSINK HEADS



Description	System	Туре	Implant Ø	REF	Price cat.
DRA 3.3 CS	SSO [®] , STO [®]	non sterile	3.3 mm	13-425046	D
DRA 4.1 CS	SSO [®] , STO [®]	non sterile	4.1 mm	13-425047	D
DRA 4.8 CSW	SSw®	non sterile	4.8 mm	13-425048	D

TAP (MADE OF CRONIDUR) FOR SSO[®], STW[®]

	Description	Finish	Length	Implant Ø	REF	Price cat.
3	TAP1 3.3	9-17 mm	34 mm	3.3 mm	13-425055	G
	TAP1 4.1	9-17 mm	34 mm	4.1 mm	13-425056	G
	TAP1 4.8	9-17 mm	30 mm	4.8 mm	13-425057	G
	TAP2 3.3	9-17 mm	23 mm	3.3 mm	13-425058	G
	TAP2 4.1	9-17 mm	23 mm	4.1 mm	13-425059	G

GUIDE JACKET



Description	Amount	Material	REF	Price cat.
BFH 2.0 guide jacket for pilot drill 2.0mmd	Pack of 5	Ti6Al4V	13-425410	A
BFH 2.5 guide jacket for pilot drill 2.5mmd	Pack of 5	Ti6Al4V	13-425411	A
BFH 3.0 guide jacket for pilot drill 3.0mmd	Pack of 5	Ti6Al4V	13-425412	A
BFH 3.2 guide jacket for pilot drill 3.2mmd	Pack of 5	Ti6Al4V	13-425413	A
BFH 3.5 guide jacket for pilot drill 3.5mmd	Pack of 5	Ti6Al4V	13-425414	А



a) length	5 mm
b) hight of step	0.7 mm
c) max. Ø top	3.7 - 5 mm
d) nominal Ø	3 - 4.4 mm
e) Ø of drilling in the drill template	2.05 - 3.55 mm



Model with residual teeth for the fabrication of a drill guide for creating cavities for fixating the later drill guide for implant cavities.



Drill guide for creating cavities for later fixation of the surgical dill guide.



Surgical drill guide for safe BECES® placement. The drill sleeves are designed for 2.0 mm Twist drills.




DRILLSTOP TRAY

Not suitable for dry heat sterilizers.



Description	REF	Price €
Drillstop A	13-500881	
Drillstop C	13-500883	
Drillstop F	13-500886	
Drillstop H	13-500888	
Drillstop K	13-500891	
Drillstop L	13-500892	
Drillstop M*	13-500893	
Drillstop N*	13-500894	
Drillstop O*	13-500895	
Drill DS 2.8	13-425005	
Drill DS 3.2	13-425007	
Drill DS 3.5	13-425009	
Drill DS 4.0	13-425011	
Drill DS 4.2	13-425013	
Drill DSL 2.8	13-425006	
Drill DSL 3.2	13-425008	
Drill DSL 3.5	13-425010	
Drill DSL 4.0	13-425012	
Drill DSL 4.2	13-425014	
Tray with content	13-60034-K	779.00

INSERTION TOOLS FOR SSO®, STO®, STW®



* This instrument is with an additional 6-edge for force transmission to the contra-angle handpiece equipped. Fits to W&H contra-angle handpiece only.

INSTRUMENTS FOR SCREWS

	Description	Туре	Length	Ø	Code	REF	Pricekat.
	Hex-Instrument	long	21 mm	1.25 mm	HT 1.25	13-425100	С
ä	Hex-Instrument	long, for contra-angle	26.1 mm	1.25 mm	HT 1.25 M	13-425112	В
	Hex-Instrument	short	14 mm	1.25 mm	HTS 1.25	13-425101	С
	Hex-Instrument	extralong	45 mm	1.25 mm	HTX 1.25	13-425102	С
	Hex-Instrument	long	19 mm	1.77 mm	HT 1.77	13-425103	С
	Hex-Instrument	long, for contra-angle	28.6 mm	1.77 mm	HT 1.77 M	13-425113	В
	Hex-Instrument	extra long	45 mm	1.77 mm	HTX 1.77	13-425104	С
*	Torx-Instrument	for SSO and STO Abutments	21 mm	1.25 mm	TT 1.25	13-425105	C



TOOLS

	Description	Туре	Code	REF	Price cat.
	Guide jacket	for pilot drill, Titanium, 10 mm, 2.2 mm Ø 5 Pieces / pack	BFH	13-425401	А
680	X-ray measuring sphere	Surgical steel, 0.5 mm Ø 5 Pieces / pack	RM	13-425403	А
	Drill extension contra-angle	extends by 19 mm	DX2	13-500704	D
Батана (Punch	for contra-angle, 4.9 mm Ø	PUW1	13-425404	С
	Punch	manual, 5.2 mm Ø	PU	13-425406	С
	Standardized probe	Scale 1 mm for radiographical measuring for cylindrical implant types, 22 mm	PDG	13-425400	А
A ANDERE NAME OFFICIAL CONTINUES	HAS	Flat spanner for ITV, Vierkant, 85 mm	HAS	13-463108	н
HOLDENALT HAS 2 CHOIL NO. 4ESTICE	HAS2	Flat spanner for TB3 STI, 85 mm	HAS2	13-463109	н
	Ratchet	angled, for all Hex-Instruments and Insertion Tools	RAT 2	13-425051	К
	Torque wrench	Torque wrench 10 - 70 Ncm for all Insertion Tools, Hex- and Torx-instruments	TW 2	13-425402	S

STARTER TRAY

Autoclavable up to 134° C, not suitable for dry heat sterilizers This surgical kit contains all drills and tools for first works with the S-System. Material: autoclavable plastic



Description	REF	Price €
IT3 STO	13-425063	
ITWH SSO	13-463111	
HT 1.25	13-425100	
TT 1.25	13-425105	
BCD 1	13-900240	
DS 2.2	13-425003	
DS 2.8	13-425005	
DS 3.5	13-425009	
DS 4.2	13-425013	
DRA 4.1 CS	13-425047	
PDG	13-425400	
PDG	13-425400	
TW2	13-425402	
Starter Tray with content	13-S60042-K	upon request
Starter Tray empty	13-60042-K	upon request

IMPLANT SYSTEM FOR ENDOSSEOUS DENTAL IMPLANTIONS



Safely anti-rotational thanks to its internal precision square

Cone technology for a tight seal

Universally suitable for fixed and removable prosthodontics

The cone centers the abutment and provides 100% tightness

The surface of **Bone Level Plus**[®] implants provide a specially lasered surface with exactly defined properties. For anti-rotation an internal square connects with pressfit to the abutment. The cone in combination with the internal stare provides stability and 100% tightness. **Bone Level Plus**[®] implants are universally applicable for fixed and removable prosthetics.

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007). Product dimension described in this brochure may differ from reality for technical reasons. Bone Level Plus® implants are protected by patents. Bone Level Plus® is a registered trademark. In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.



PREPARATORY STEPS WHEN USING A DRILLING TEMPLATE

- 1. Have your laboratory produce a drilling template with the appropriate drill holes for the marker bore. To be on the safe side, the laboratory might insert guide sleeves (REF BFH) into the drill holes to ensure that the drilling angle is correct. Use a 2 mm ø drill for pilot drilling.
- For subsequent drilling sequences, drill stops can be used that are slid over the drill according to the appropriate depth of the drill hole and screwed in place. Consider the thickness of the mucosa and the height of the template as appropriate.


Thanks to the extremely high cutting efficiency of our drills, no ascending drilling sequences will usually be required.


Recommended RPM: 2000-5000

Apply sufficient cooling and allow the cooling to reach the working blades of the drills.

SURGERY

1. Recommended drill sequence





Owing to the high quality and geometry of the blades of our drills, the final preparation may be performed immediately after the pilot drilling.



Original packaging

13-900509 Bone Level P Ti6Al4V	lus 4.1 8 Implant
LOT 160661007	30.11.2018
	-30.11.2013
	STERAL
001	mototan
Mersheller / Manufachrer Simplodent Cristii - 610 Comme	-


Open the blister using the flap. Remove the label and stick it into the patients record.



The blister (secondary packaging) contains the implant in a sterile tube (primary package).

3. Remove the implant from its packaging

- 1. Open the lid. The implant is connected to the lid through a breakable section.
- 2. Remove the implant without touching the inner walls of the tube.




4. Handling


4.1 Connect

Attach the placement aid to the implant, holding the cap to which the implant is attached with the other hand.

4.2 Mounting the adapter ITV WST / contra-angle

Place the ITV Wst (angled handpiece) or IT ITV (ratchet) adapter on the ITV BLP placement aid. Mount the placement aid. Hold the cap firmly in one hand and break off the implant at the pre-determined breaking line.







4.3 Alternative to 4.2:

Place the IT ITV (ratchet) adapter on the ITV BLP placement aid.

Mount the placement aid. Hold the cap firmly in one hand and break off the implant at the pre-determined breaking line.



5. Insertion

Use the angled handpiece, ratchet or shank to screw the implant into the implant bed (clockwise).

The enossal aspect of the implant must be submerged in the bone. Upon **complete** insertion, the implant may be turned back ¼ revolution to reduce the load on the bone.


The system is suitable for deep insertion (below bone level).



6. Remove insertion tool from implant

The insertion tool may be retrieved from the implant with the help of the contra-angle instrument.

As an alternative use the ratchet RAT2 + IT ITV + HAS (flat spanner).



7. Result

Bone Level Plus® implant

Jaw




8. Aftercare

Seal the implant with a matching cover screw.




Jaw

Arrow: Top view.



9. Pick-up impressions



10. Closed tray impression taking


10.1 Impression with closed tray

Impression with custom tray.

Securing the impression post with the thumbscrew

TS B NC/RC

Bone Level Plus® implant



10.2 Before taking the impression

Take an impression in an A silicone. You can use the open-tray or the closed-tray technique.

With the closed impression technique, the TS BLP NC/RC will always remain on the implant when removing the impression. Impression post TS B NC/RC


Bone Level Plus® implant


10.3 Removing the impression

In the case of closed impressions, the TS BLP NC/RC impression post will remain on the implant after removing the impression tray.

> Impression post TS B NC/RC

Bone Level Plus® implant





10.4

Once the impression has been taken, the implant is closed with an HA NC/RC healing abutment, while the impression is sent to the laboratory.

TT 1.25

insert surgical screw CSTB NC/RC Arrow: Top view. Rotation clockwise

Bone Level Plus® implant

11. Procedures in the laboratory

11.1 Preparing the impression tray for fabricating the cast

Pick-up technique Tighten the IAB or M analogue IABM against the LTB (NC/RC) impression post.

Use the TT 1.25 to insert the lab analogue

HLTB NC/RC

IAB NC or IAB RC



Secure the IAB or M analogue or IAB NC/ RC against the TS BLP (NC or RC) (A)

Reposition the impression post inside the impression $\textcircled{\textbf{B}}$

Pour the impression.

Use the thumbscrew to tighten the impression post on the lab analogue.

TSB NC/RC



Arrow: Top view.

Rotation clockwise

IAB NC or IAB RC

11.3

Pour the impression in dental stone, then remove the impression posts from the lab analogues.

Lab analogue

Fill with gypsum



11.4

The lab analogue will now be embedded in the gypsum in the correct position.

IAB NC/RC

11.5

Positioning of the screw-retained TLA2 15 BLP RC abutment, determining its optimal position and correct angulation.

NOTE The hexagon must be inserted completely into the analogue.

TT 1.25

Insert screw

TLA2 15 BLP RC Watch out for the correct hexagon position

IAB NC/RC

Arrow: Top view. Rotation clockwise

Gypsum

11.6

The correct position of the abutment must be ensured during transfer to the mouth.

Tighten the implant screw to a torque of 20 Ncm.

TLA15 BLP RC



11.7

If multiple angled abutments are used, the laboratory will produce a removable resin splint (e.g. from pattern resin) to facilitate positioning within the mouth.

TLA2 15 BLP RC



Pattern Resin

BONE LEVEL PLUS® IMPLANTS


The surface is roughed. The implant body is made of Ti6Al4V.

a	_	Description BLP 3.3 8 NC	Enossal Ø 3.3 mm	Enossal length 8 mm	REF 13-900500	Price cat. H
			3.3 mm	10 mm	13-900501	н
	b	BLP 3.3 12 NC	3.3 mm	12 mm	13-900502	н
b		BLP 3.3 14 NC	3.3 mm	14 mm	13-900503	н
		BLP 4.1 8 RC	4.1 mm	8 mm	13-900504	Н
		BLP 4.1 10 RC	4.1 mm	10 mm	13-900505	Н
		BLP 4.1 12 RC	4.1 mm	12 mm	13-900506	Н
		BLP 4.1 14 RC	4.1 mm	14 mm	13-900507	н
a) enossal Ø	3.3 - 4.8 mm	BLP 4.8 8 RC	4.8 mm	8 mm	13-900508	н
b) enossal length	8 - 14 mm	BLP 4.8 10 RC	4.8 mm	10 mm	13-900509	н
		BLP 4.8 12 RC	4.8 mm	12 mm	13-900510	н
		BLP 4.8 14 RC	4.8 mm	14 mm	13-900511	н



NC = Narrow Collar

RC = Regular Collar



BLP* implants are delivered incl. insertion tool ITV BLP and surgical screw REF 13-900518 or 13-900519.



Safely anti-rotational thanks to its internal precision square

- Cone technology for a tight seal
- Universally suitable for fixed and removable prosthodontics
- The cone centers the abutment and provides 100% tightness



GINGIVAFORMER



Description		Code	REF	Price cat.
Gingivaformer	conical	GF NC 3.6 2	13-900590	В
Gingivaformer	conical	GF NC 3.6 3.5	13-900591	В
Gingivaformer	conical	GF NC 4.8 3.5	13-900594	В
Gingivaformer	conical	GF RC 4.5 2	13-900596	В
Gingivaformer	conical	GF RC 4.5 4	13-900597	В
Gingivaformer	conical	GF RC 4.5 6	13-900598	В
Gingivaformer	conical	GF RC 6 2	13-900599	В
Gingivaformer	bottle shape	GFB NC 3.3 3.5	13-900602	В
Gingivaformer	bottle shape	GFB NC 3.3 5	13-900603	В
Gingivaformer	bottle shape	GFB RC 4.4 4	13-900604	В
Gingivaformer	bottle shape	GFB RC 4.7 6	13-900605	В

BUR CYLINDER



Description	Code	REF	Price cat.
Bur cylinder for BLP 3.3 for telescope crowns	FZB NC	13-900524	D
Bur cylinder for BLP 4.1 and 4.8 for telescope crowns	FZB RC	13-900527	D

Recommended insertion torque 30 Ncm

ANALOGUES

	l

Description	Code	REF	Price cat.
Implant analogue for BLP 3.3	IA BLP NC	13-900525	В
Implant analogue for BLP 4.1 and 4.8	IA BLP RC	13-900526	В

STANDARD ABUTMENTS



Description	Code	REF	Price cat.
Abutment for cementing on BLP 3.3, step 1 mm high Hight above step 4 mm, incl. matching screw SFBC NC	CAB 1 NC	13-900554	E
Abutment for cementing on BLP 4.1 und 4.8, step 1 mm high Hight above step 5.5 mm, incl. matching screw SFBC RC	CAB 1 RC	13-900551	E
Abutment for cementing on BLP 3.3, step 3 mm high Hight above step 4 mm, incl. matching screw SFBC NC	CAB 3 NC	13-900555	E
Abutment for cementing on BLP 4.1 und 4.8, step 3 mm high Hight above step 5.5 mm, incl. matching screw SFBC RC	CAB 3 RC	13-900552	E
Recommended insertion torque 20 Ncm			

SCREW-RETAINED ABUTMENTS (REDUCIBLE, GRINDABLE)

	Description	Code	REF	Price cat.
	Abutment Incl. matching screw SF B	TAB BLP NC/RC	13-900521	D
ļ	Abutment for BLP 3.3 15° angle, anti-rotational Incl. matching screw SF B NC	TLA2 15 BLP NC	13-900528	F
	Abutment for BLP 4.1 und 4.8 15° angle, anti-rotational Incl. matching screw SF B NC	TLA2 15 BLP RC	13-900523	F

Recommended insertion torque 20 Ncm

ANATOMICAL ABUTMENTS

Description	Code	REF	Price cat.
Anatomical abutment for BLP 3.3 Incl. matching screw SFB NC	ANAB NC	13-900544	F
Anatomical abutment for BLP 4.1 and 4.8 Incl. matching screw SFB RC	ANAB RC	13-900543	F

Recommended insertion torque 20 Ncm



TITAN BASE FOR CAD CAM



Description	REF	Price cat.
MB BLP NC, anti-rotation	13-900560	D
MB BLP RC, anti-rotation	13-900562	D

CASTABLE ABUTMENTS



C	Description	Material	Code	REF	Price cat.
	Castable abutment incl. metal base ind screw	Ti6Al4V	PLAB BLP NC	13-900620	D
	Castable abutment incl. metal base and screw	Ti6Al4V	PLAB BLP RC	13-900622	D

PICK-UP IMPRESSION POST FOR PICK-UP IMPRESSIONS



Description	Code	REF	Price cat.
Impression post for BLP 3.3	HLT BLP NC	13-900584	С
Impression post for BLP 4.1 and 4.8	HLT BLP RC	13-900585	С

IMPRESSION POST FOR CONVENTIONAL IMPRESSIONS

Description	Code	REF	Price cat.
Impression post for BLP 3.3	TS BLP NC	13-900586	С
Impression post for BLP 4.1 and 4.8	TS BLP RC	13-900587	С
Impression post long for BLP 3.3	TSL BLP NC	13-900588	С
Impression post long for BLP 4.1 and 4.8	TSL BLP RC	13-900589	С

ABUTMENTS FOR SCREW-ON PROSTHETIC



Description TCT BLP RC 0.5	Code TCT BLP RC 0.5	REF 13-900632	Price cat. D
TCT BLP RC 1.5	TCT BLP RC 1.5	13-900633	D
TCT BLP RC 3.5	TCT BLP RC 3.5	13-900634	D
TCT BLP NC 0.5	TCT BLP NC 0.5	13-900635	D
TCT BLP NC 1.5	TCT BLP NC 1.5	13-900636	D
TCT BLP NC 3.5	TCT BLP NC 3.5	13-900637	D

Fasten with HT 1.77

IMPRESSION TAKING AND LABORATORY ACCESSORIES

In this approach the position of the TCT hex is assigned.





LOCALICER° FOR REMOVABLE PROSTHETIC

LOC abutments are mounted with the tool HT 1.77. LOC abutments are used for connection of removable prosthetics with Bone Level Plus® implants. If LOC abutments (or KOC LOC implants) are used in the upper jaw, we recommend to place at least six implants and to splint them through prosthetics in a stable manner.

		Description	Hight	Code	REF	Price cat.
		Localicer [®] for BLP 3.3	2 mm	LOC BLP NC 2	13-900539	D
		Localicer [®] for BLP 4.1 und 4.8	2 mm	LOC BLP RC 2	13-900540	D
		Localicer [®] for BLP 3.3	3 mm	LOC BLP NC 3	13-900606	D
		Localicer [*] for BLP 3.3	4 mm	LOC BLP NC 4	13-900607	D
Y		Localicer* for BLP 4.1 and 4.8	3 mm	LOC BLP RC 3	13-900608	D
	Y	Localicer [*] for BLP 4.1 and 4.8	4 mm	LOC BLP RC 4	13-900609	D

ACCESSORIES FOR LOCALICER°





MULTI-UNIT ABUTMENTS

Insertion of the angled MU2-abutments with HT 1.25. Insertion of the straight MU2S-abutments with HT 1.77.

I	h= 3 mm	mm	•
	6 mm		
ŀ		6,7 mm 	4

Description	Material	Code	REF	Price cat.
MU2 17 BLP RC, angled, incl. SFB RC	Ti6Al4V	MU2 17 BLP RC	13-900640	L
MU2 35 BLP RC, angled, incl. SFB RC	Ti6Al4V	MU2 35 BLP RC	13-900641	L
MU2S 0.5 BLP RC, straight	Ti6Al4V	MU2S 0.5 BLP RC	13-900642	G
MU2S 1.5 BLP RC, straight	Ti6Al4V	MU2S 1.5 BLP RC	13-900643	G
MU2S 2.5 BLP RC, straight	Ti6Al4V	MU2S 2.5 BLP RC	13-900644	G
GF MU2 gingivaformer incl. SF MU2 Hight above abutment shoulder 6 mm	Ti6Al4V	GF MU 2	13-418286	С
MU2 Localicer incl. SF MU2 Hight above abutment shoulder 6.7 mm	Ti6Al4V	MU 2	13-418287	С
Prosthetic screw for MU2	Ti6Al4V	SFB RC	13-900532	A

ACCESSORIES FOR MULTI-UNIT ABUTMENTS



Description	Material	Code	REF	Price cat.
Temporary base SF MU2 sold separately	Ti6Al4V	TC MU2	13-418290	D
Transfer straight incl. screw SFL MU2	Ti6Al4V	TS MU2	13-418291	С
Castable for Multi-Unit incl. screw		PA MU2	13-418292	А
Screw for TC MU2	Ti6Al4V	SF MU2	13-418293	В
Lab analogue for Multi-Unit	Ti6Al4V	IA MU2	13-418295	В
Hex instrument long, Ø 1.25 mm		HT 1.25	13-425100	С
Hex instrument extralong, 45 mm, Ø 1.25 mm		HTX 1.25	13-425102	С
Hex instrument for suprastructures, Ø 1.77 mm		HT 1.77	13-425103	С

INSTRUMENTS

Description	Code	REF	Price cat.
Pilot drill short/long 2.0 mm Ø	DS 2 / DSL 2	13-425001 / 13-425002	D
Pilot drill short/long 2.8 mm Ø	DS 2.8 / DSL 2.8	13-425005 / 13-425006	D
Form drill short 2.8 mm Ø	DBL 2.8	13-900570	E
Form drill short 3.5 mm Ø	DBL 3.5	13-900571	E
Form drill short 4.0 mm Ø	DBL 4.0	13-900572	E
Cortical countersink 3.3	CSBL 3.3	13-900576	D
Cortical countersink 4.1	CSBL 4.1	13-900577	D
Cortical countersink 4.8	CSBL 4.8	13-900578	D
Тар	TAP BLP 3.3	13-900579	D
Тар	TAP BLP 4.1	13-900580	D
Тар	TAP BLP 4.8	13-900581	D




GUIDE JACKET





Description	Amount	Material	REF	Price cat.
BFH 2.0 guide jacket for pilot drill 2.0mmd	Pack of 5	Ti6Al4V	13-425410	A
BFH 2.5 guide jacket for pilot drill 2.5mmd	Pack of 5	Ti6Al4V	13-425411	A
BFH 3.0 guide jacket for pilot drill 3.0mmd	Pack of 5	Ti6Al4V	13-425412	A
BFH 3.2 guide jacket for pilot drill 3.2mmd	Pack of 5	Ti6Al4V	13-425413	A
BFH 3.5 guide jacket for pilot drill 3.5mmd	Pack of 5	Ti6Al4V	13-425414	A

. .



a) length	5 mm
b) hight of step	0.7 mm
c) max. Ø top	3.7 - 5 mm
d) nominal Ø	3 - 4.4 mm
e) Ø of drilling in the drill template	2.05 - 3.55 mm



Model with residual teeth for the fabrication of a drill guide for creating cavities for fixating the later drill guide for implant cavities.



Drill guide for creating cavities for later fixation of the surgical dill guide.




Surgical drill guide for safe BECES® placement. The drill sleeves are designed for 2.0 mm Twist drills.

ADAPTER



Description	For	Length	Code	REF	Price cat.
Adapter short / contra-angle	ITV 500850	22 mm	ITV S WST	13-500851	С
Adapter long / contra-angle	ITV 500850	32 mm	ITV L WST	13-500852	С
Adapter medium / contra-angle	ITV 500850	27 mm	ITV M WST	13-500853	С
Ratchet adapter	Adapter zu ITV		IT ITV	13-500854	С
Drill extension, contra-angle, extends by 1	9 mm		DX2	13-500704	D

INSTRUMENTS AND **TOOLS**



* It is recommended to have the torque ratchets recalibrated by us once a year.

** for cleaning this instrument an ultrasonic cleaning device and a thermo-desinfector (i.e. Miele TD-series) are required.

If these devices are not available in the dental office the handle with REF 311430 should be purchased instead.

SIMPLADENT®





Description	REF	Price €
ITV S adapter short	13-500851	
ITV M	13-500852	
ITITV	13-500854	
TT 1.25	13-425105	
CSBL 3.3	13-900576	
CSBL 4.1	13-900577	
CSBL 4.8	13-900578	
DS 2.0	13-425001	
DBL 2.8	13-900570	
DBL 3.5	13-900571	
DBL 4.0	13-900572	
TAP BLP 3.3	13-900579	
TAP BLP 4.1	13-900580	
TAP BLP 4.8	13-900581	
TW2 torque wrench	13-425402	
Starter Tray empty	13-60045-K	upon request
Starter Tray with content	13-S60045-K	upon request

SURGICAL INSTRUMENT TRAY



Description	REF	Price €
DS 2	13-425001	
DBL 2.8	13-900570	
DBL 3.5	13-900571	
DBL 4.0	13-900572	
DLBL 2.8	13-900573	
DLBL 3.5	13-900574	
DLBL 4.0	13-900575	
PDG	13-425400	
PDG	13-425400	
PDG	13-425400	
CSBL 3.3	13-900576	
CSBL 4.1	13-900577	
CSBL 4.8	13-900578	
TAP BLP 3.3	13-900579	
TAP BLP 4.1	13-900580	
TAP BLP 4.8	13-900581	
ITITV	13-500854	
ITV S adapter short	13-500851	
ITV M adapter medium	13-500853	
ITV L adapter long	13-500852	
UAW	13-425107	
PUW 1	13-425404	
TT 1.25	13-425105	
DX 2	13-500704	
TW2 torque wrench	13-425402	
Starter Tray empty	13-60018-K	upon request
Starter Tray with content	13-S60018-K	upon request



IMPLANT SYSTEM FOR ENOSSAL DENTAL IMPLANTATIONS



The Allfit[®] implant system **Xign[®]** is designed for enossal dental implantations. **Xign[®]** implants feature two high-precision internal cylindrical guides and hexagon. Implants with diameter 3.0 and 3.4 are not indicated for single tooth replacement. The **Xign[®]** implant system is suitable for dual-step implantation protocols. **Xign[®]** accessory parts are colour coded.


We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007). Due to technical reasons the product dimensions shown in this brochure might deviate from reality. Xign[®] implants are patent-protected. Xign[®] is a registered trademark. In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.

PRELIMINARY WORK FOR TEMPLATE APPLICATION

- Ask your laboratory to prepare a drill template with boreholes for the pilot drill. To be on the safe side, you can ask the laboratory to insert guide sleeves (code BFH) into the template, which specify the exact drill direction. Please use a 2.0 Ø drill for the pilot drilling.
- 2. For the following drill sequences you can use drill stops, which can be attached and tightened to the drill according to the length of the drilling channel. Mucous membrane thickness and template height are taken into account as needed. Due to the extreme high cutting performance of our drills after the pilot drill the final drill may be used right away. In this case the drill sequence may be disregarded.

Recommended RPM: 2000-5000


Apply sufficient cooling and allow the cooling to reach the working blades of the drills.



GENERAL NOTE Xign® implants with diameter 4.5 mm and larger are used as compression screws. In order to acchieve a good bone condensation and implant stability, the drilling should be carried out thinner than the core diameter of the implant. The minimal diameter of the drill depends on the bone density. It is therefore not possible to advise drill-sequences which fit all bone-qualities. Typically in the soft maxillary bone only small drill-diameters are used (e.g. the usage of DOS1 only for larger Xign® implants), whereas in the highly mineralized lower jaw a specific drill sequence with respect to the mineralisation of the bone is necessary. For insertion ander pressure use the handgrip.



SURGERY



1. Recommended drill sequence (for contra-angle handgrip)

* For areas with hard cortical bone, enlargement using a C-drill orCortical Tap (KGS) is recommended down to a cortical depth of up to 6 mm. Form drills create the bone cavity for the respective implant. The actual diameter of the form drill is not mentioned in its order code (REF), as the drill correponds to the diameter of the implants core.



Original packaging

	(I 3.8 13 Ilfit Implant Syster iGrade4	em Xi,
	920708A 2 31.0	
Allfi	t®	


Open the blister using the flap. Remove the label and stick it into the patients record.



The blister (secondary packaging) contains the implant in a sterile tube (primary package).

3. Remove the implant from its packaging

- 1. Open the lid.
- 2. The implant is fixed to the lid by a break joint.
- 3. Remove the implant without touching the inner wall of the tube.



4. Handling

4.1 Connect

Attach the insertion tool to the implant by holding the lid, to which the implant is secured, with your other hand.


4.2

After you have attached the insertion tool, firmly hold the top in your hand and break the implant off the top along the break joint.

4.3

Using the contra-angle handpiece screw the implant clockwise into the cavity.

The enossal part of the implant must be completely covered by the bone. After insertion the implant can be turned by a 1/4 rotation backwards in order to relieve the bone and allow blood access to the implant surface.













Max. insertion torque Xign 3.0 = 20 Ncm Xign 3.4 = 30 Ncm Xign 3.8 = 40 Ncm

8. Pick-up impressions







8.2 Prior to the impression


Impression taking with an A silicone. The use of open and closed impression tray is possible.

The impression post HLT Xi must be unscrewed from the implant before the impression tray is removed form the patients mouth.





Xign® implant



9. Impression taking with a closed impression spoon




9.2 Prior to the impression

Impression taking with an A silicone The use of open and closed impression tray is possible.

Impression post TS(L) Xi

Xign[®] implant



9.3 Remove impression

When the closed impression method is applied, then the impression post TS Xi remains on the implant after the impression tray is removed.

After removal of the impression tray the transferpost TS(L) Xi is unscrewed.

9.4

After the impression is taken, the implant is closed with a cover screw (i.e. Healing screw) and the impression is sent to the laboratory.

HT 1.25


TS(L) Xi

Xign[®] implant

Fasten Gingiva shaper HS Xi



Xign® implant





10.4 Impression taking	HT 1.25	Arrow Top view.	
Positioning of the screwed abutments TLA, whereu- pon the optimal position and adequate angulation must be determined.	Tighten screw	Direction clockwise.	
Note The hexagon must be completely inserted into the analogue.	TLA 15 Xi Ensure proper position of the hexagon	Ý	
	IA Xi / IA Xi M		Gypsum

10.5

Ensure proper position of the abutment when transferring into the mouth.

Tightening torque of the screw during fastening on the implant: 15 Ncm



10.6


If more than one angled abutment is used, your laboratory will prepare a detachable synthetic bar (e.g. Pattern Resin) in order to facilitate positioning in the mouth. TLA 15 Xi

TLA 15 Xi

Pattern Resin



Xign* implants feature two high-precision internal cylindrical guides and hexagon. Implants with diameter 3.0 and 3.4 are not indicated for single tooth replacement. Material Ti6Al4V.



For Implants with 3.0 mmd the max. torque for the prosthetic screw is **30 Ncm**. For Implants with 3.4 mmd the max. torque for the prosthetic screw is **80 Ncm**. For Implants with 3.8 mmd the max. torque for the prosthetic screw is **60 Ncm**.

COLOR CODING

Implant Ø





CYLINDRICAL GINGIVAFORMER

P

P

Y

Φ

Description	Code	REF	Price cat.
Gingivaformer 3.0 mm Ø, brown for 2 mm mucous membrane height	HS Xi 3.0 2	13-422304	В
Gingivaformer 3.0 mm Ø, brown for 4 mm mucous membrane height	HS Xi 3.0 4	13-422305	В
Gingivaformer 3.4 mm Ø, silver for 2 mm mucous membrane height	HS Xi 3.4 2	13-422306	В
Gingivaformer 3.4 mm Ø, silver for 2 mm mucous membrane height	HS Xi 3.4 4	13-422307	В
Gingivaformer 3.8 mm Ø, yellow for 2 mm mucous membrane height	HS Xi 3.8 2	13-422308	В
Gingivaformer 3.8 mm Ø, yellow for 4 mm mucous membrane height	HS Xi 3.8 4	13-422309	В
Gingivaformer 4.5 mm Ø, blue for 2 mm mucous membrane height	HS Xi 4.5 2	13-422310	В
Gingivaformer 4.5 mm Ø, blue for 4 mm mucous membrane height	HS Xi 4.5 4	13-422311	В

ANATOMICAL GINGIVAFORMER

	Description	Hight/width	Code	REF	Price cat.
Y	Gingivaformer 3.0 mm Ø, brown Hight 4 mm, width at top 6 mm	h4w6	HS Xi 3.0	13-422480	В
Y	Gingivaformer 3.4 mm Ø, silver Hight 4 mm, width at top 6 mm	h4w6	HS Xi 3.4	13-422481	В
Y	Gingivaformer 3.8 mm Ø, yellow Hight 4 mm, width at top 6 mm	h4w6	HS Xi 3.8	13-422482	В
Y	Gingivaformer 4.5 mm Ø, blue Hight 4 mm, width at top 6 mm	h4w6	HS Xi 4.5	13-422483	В

SURGICAL SCREWS

Description	Code	REF	Price cat.
Surgical screw for 3.0 mm Ø brown	CST Xi 3.0	13-422300	В
Surgical screw for 3.4 mm Ø silver	CST Xi 3.4	13-422301	В
Surgical screw for 3.8 mm Ø yellow	CST Xi 3.8	13-422302	В
Surgical screw for 4.5 mm Ø blue	CST Xi 4.5	13-422303	В

ABUTMENTS FOR CEMENTED PROSTHETICS

		Screw / REF	Description	H mucous membran	Color	REF	Price cat.
		-	TLA Xi 3.0 2	2 mm	brown	13-422317	D
		SF 2305	TLA Xi 3.0 3	3 mm	brown	13-422318	D
	¥.	0	TLA Xi 3.0 5	5 mm	brown	13-422319	D
1			TLA Xi 3.4 1	1 mm	silver	13-422320	D
			TLA Xi 3.4 2	2 mm	silver	13-422321	D
	W		TLA Xi 3.4 3	3 mm	silver	13-422322	D
	. e		TLA Xi 3.4 5	5 mm	silver	13-422323	D
			TLA Xi 3.4 8	8 mm	silver	13-422324	D
Straight		SF 2301	TLA Xi 3.8 1	1 mm	yellow	13-422325	D
	1		TLA Xi 3.8 2	2 mm	yellow	13-422326	D
	П		TLA Xi 3.8 3	3 mm	yellow	13-422327	D
			TLA Xi 3.8 5	5 mm	yellow	13-422328	D
	¥.		TLA Xi 4.5 1	1 mm	blue	13-422329	D
			TLA Xi 4.5 2	2 mm	blue	13-422330	D
1	NI)		TLA Xi 4.5 3	3 mm	blue	13-422331	D
	4		TLA 15 Xi 3.4 1	1 mm	silver	13-422332	F
			TLA 15 Xi 3.4 2	2 mm	silver	13-422333	F
			TLA 15 Xi 3.4 3	3 mm	silver	13-422334	F
	All		TLA 15 Xi 3.8 1	1 mm	yellow	13-422335	F
Angled 15°	1	SF 2302	TLA 15 Xi 3.8 2	2 mm	yellow	13-422336	F
	1		TLA 15 Xi 3.8 4	4 mm	yellow	13-422337	F
	AT .		TLA 15 Xi 4.5 1	1 mm	blue	13-422338	F
			TLA 15 Xi 4.5 2	2 mm	blue	13-422339	F
	¥		TLA 15 Xi 4.5 4	4 mm	blue	13-422340	F
			TLA 25 Xi 3.4 1	1 mm	silver	13-422460	F
	III		TLA 25 Xi 3.4 2	2 mm	silver	13-422461	F
			TLA 25 Xi 3.4 3	3 mm	silver	13-422462	F
	in the second se		TLA 25 Xi 3.8 1	1 mm	yellow	13-422463	F
Angled 25°	III	SF 2409	TLA 25 Xi 3.8 2	2 mm	yellow	13-422464	F
			TLA 25 Xi 3.8 3	3 mm	yellow	13-422469	F
	E		TLA 25 Xi 3.8 4	4 mm	yellow	13-422465	F
	M		TLA 25 Xi 4.5 1	1 mm	blue	13-422466	F
			TLA 25 Xi 4.5 2	2 mm	blue	13-422467	F
			TLA 25 Xi 4.5 4	4 mm	blue	13-422468	F
			Castable abutment,	plastic	PA U	13-418181	А

Tightening torque:

For Xi 3.0 = 13 Ncm For Xi 3.4, 3.8, 4.5 = 20 Ncm

Incl. screw

ABUTMENTS FOR ZIRKONIUM

4 1	S

	Description	Material	Code	REF	Price cat.
	Titanbasis for XI 3.0, antirotation	Ti6AIA4	MX B Xi 3.0	13-422473	D
SF 2301	Titanbasis for XI 3.4, antirotation	Ti6AIA4	MX B Xi 3.4	13-422470	D
	Titanbasis for XI 3.8, antirotation	Ti6AIA4	MX B Xi 3.8	13-422471	D
	Titanbasis for XI 4.5, antirotation	Ti6AIA4	MX B Xi 4.5	13-422472	D
	Incl. screw SF 2301, insertion with HT 1.25				

IMPRESSION AND LAB ACCESSORIES FOR CEMENTED PROSTHETICS

Description	Code	REF	Price cat.
Impression post with top screw for manual operation, short	TS Xi 3.0 S	13-422345	С
	TS Xi 3.4 S	13-422346	С
	TS Xi 3.8 S	13-422347	С
	TS Xi 4.5 S	13-422348	С
Impression post with long screw	HLT Xi 3.0	13-422349	С
	HLT Xi 3.4	13-422350	С
	HLT Xi 3.8	13-422351	С
	HLT Xi 4.5	13-422352	С
Implant analogue with internal hexagon	IA Xi 3.0	13-422353	В
	IA Xi 3.4	13-422354	В
	IA Xi 3.8	13-422355	В
	IA Xi 4.5	13-422356	В
Castable abutment, plastic	PA U	13-418181	A

ABUTMENTS FOR SCREW RETAINED RESTORATIONS



OSA abutment		Code	REF	Price cat.
Abutment for Xi 3.4, silver 2 mm transmucosal height	two part with screw	OSA Xi 3.4 2	13-422361	E
Abutment for Xi 3.8, yellow, 2 mm transmucosal height	two part with screw	OSA Xi 3.8 2	13-422362	Е
Abutment for Xi 3.8, yellow, 2 mm transmucosal height	two part with screw	OSA Xi 3.8 4	13-422363	E
Abutment for Xi 4.5, blue, 2 mm transmucosal height	two part with screw	OSA Xi 4.5 2	13-422364	E

IMPRESSION AND LAB ACCESSORIES



BUR CYLINDER

I	I	
1		
P	Y	

Description	Code	REF	Price cat.
Bur cylinder for Xi 3.0 incl. screw SF 2301*	FZ Xi 3.0	13-422365	D
Bur cylinder for Xi 3.4 incl. screw SF 2301	FZ Xi 3.4	13-422366	D
Bur cylinder for Xi 3.8 incl. screw SF 2301	FZ Xi 3.8	13-422367	D
Bur cylinder for Xi 4.5 incl. screw SF 2301	FZ Xi 4.5	13-422368	D

*The bur cylinder FZ Xi 3.0 is not for removable prosthetics and especially not for telescopic restorations.

LOCALICER®



Description	Hight	Code	REF	Price cat.
Localizer [®] for Xign [®] 3.4 mm	2 mm	LOC Xi 3.4 2	13-422372	D
Localizer [®] for Xign [®] 3.4 mm	3 mm	LOC Xi 3.4 3	13-422373	D
Localizer [®] for Xign [®] 3.4 mm	4 mm	LOC Xi 3.4 4	13-422374	D
Localizer [®] for Xign [®] 3.8 mm	2 mm	LOC Xi 3.8 2	13-422375	D
Localizer [®] for Xign [®] 3.8 mm	3 mm	LOC Xi 3.8 3	13-422376	D
Localizer [®] for Xign [®] 3.8 mm	4 mm	LOC Xi 3.8 4	13-422377	D
Localizer [®] for Xign [®] 4.5 mm	2 mm	LOC Xi 4.5 2	13-422378	D
Localizer [®] for Xign [®] 4.5 mm	3 mm	LOC Xi 4.5 3	13-422379	D
Localizer [®] for Xign [®] 4.5 mm	4 mm	LOC Xi 4.5 4	13-422380	D

Required tool HT 1.77.

When using LOC-abutments and KOS LOC, we recommend to use minimally 6 implants per jaw and to use a single denture as splinting.

ACCESSORIES FOR LOCALICER°



Description Analogue + impression Set	Code AA LOC	REF 13-462337	Price cat. C
Set with 5 caps + 1 housing (EXTERNAL PRODUCT)	NCS	13-462338	D
Pull off force Yellow 600 g, Pink 1.200 g, Clear 1.800 g, Violet 2.700 g. Black has no retention and is designed for temporary solutions for up to one month.			

BALL ABUTMENTS

H	

Description	Code	REF	Price cat.
Ball abutment for Xign® 3.0 headØ 2.5mm; 2 mm Hight	TB Xi 3.0 2	13-422381	E
Ball abutment for Xign® 3.0 headØ 2.5mm; 4 mm Hight	TB Xi 3.0 4	13-422382	E
Ball abutment for Xign® 3.4 headØ 2.5mm; 2 mm Hight	TB Xi 3.4 2	13-422383	E
Ball abutment for Xign® 3.4 headØ 2.5mm; 4 mm Hight	TB Xi 3.4 4	13-422384	E
Ball abutment for Xign® 3.8 headØ 2.5mm; 2 mm Hight	TB Xi 3.8 2	13-422385	E
Ball abutment for Xign® 3.8 headØ 2.5mm; 4 mm Hight	TB Xi 3.8 4	13-422386	E
Ball abutment for Xign® 4.5 headØ 2.5mm; 2 mm Hight	TB Xi 4.5 2	13-422387	E
Ball abutment for Xign® 4.5 headØ 2.5mm; 4 mm Hight	TB Xi 4.5 4	13-422388	E

Insertion with HAS2 or Tool E (REF 13-462377)

ACCESSORIES FOR BALL ABUTMENTS

Description			Code	REF	Price cat.
Nylon cap transj	parent (Pull-off force ca. 1200g) for B Xi (EXTERNAL PRODUCT)	2 pieces / pack	NC	13-465028	A1
, , , , , , , , , , , , , , , , , , , ,	Pull-off force ca. 800g) for B Xi (EXTERNAL PRODUCT)	2 pieces / pack	NC 1	13-465029	A1
, , ,	w (Pull-off force ca. 500g) for B Xi (EXTERNAL PRODUCT)	2 pieces / pack	NC 2	13-465030	A1
Green, strong		2 pieces / pack	R-NC	13-465034	A1
Pink, medium	Nylon cap R-NC With increased friction strength	2 pieces / pack	R-NC1	13-465033	A1
Orange, soft	(for use with worn out Localicer®) EXTERNAL PRODUCT	2 pieces / pack	R-NC2	13-465032	A1
Sleeve for all NC	(EXTERNAL PRODUCT)		н	13-465031	A1

MULTI-UNIT ABUTMENTS



Description	Code	REF	Price cat.
MU2 17 Xi 3.8, angled, Ti6Al4V ELI. Incl. SF 44	MU2 17 Xi 3.8	13-422500	L
MU2 35 Xi 3.8, angled, Ti6Al4V ELI. Incl. SF 44	MU2 35 Xi 3.8	13-422501	L
MU2S 0.5 Xi 3.8, straight, Ti6Al4V ELI	MU2S 0.5 Xi 3.8	13-422502	G
MU2S 1.5 Xi 3.8, straight, Ti6Al4V ELI	MU2S 1.5 Xi 3.8	13-422503	G
MU2S 2.5 Xi 3.8, abutment straight, Ti6Al4V ELI	MU2S 2.5 Xi 3.8	13-422504	G
Gingivaformer, Ti6Al4V ELI. Incl. SF MU2 Hight above abutment shoulder 6 mm	GF MU2	13-418286	С
Localicer, Ti6AI4V ELI. Incl. SF MU2 Hight above abument shoulder 6.7 mm	MU2	13-418287	С

Insertion of the angled MU2-abutments with HT 1.25 Insertion of the straight MU2S-abutments with HT 1.77

ACCESSORIES FOR MULTI-UNIT ABUTMENTS



Description	Code	REF	Price cat.
Temporary base, Ti6Al4V (SF MU2 has to be ordered separately	TC MU2	13-418290	D
Transfer straight, Ti6Al4V ELI, incl. screw SFL MU2	TS MU2	13-418291	С
Castable for Multi-Unit, incl. screw	PA MU2	13-418292	А
Screw for TC MU, Ti6Al4V ELI	SE MU2	13-418293	А
Lab analogue for Multi-Unit, Ti6Al4V ELI	IA MU 2	13-418295	В
Hex-instrument, long	HT 1.25	13-425100	С
Hex-instrument, long	HT 1.77	13-425103	с

INSERTION TOOLS

	Description		Code	REF	Price cat.
	Insertion tool medium, for TempBase®	For implant	ITT2 Xi 3.0 OS	13-422408	D
	Insertion tool short, for TempBase®	Ø 3.0 mm	ITT3 Xi 3.0 OS	13-422409	D
Ta	Insertion tool long	For implant	IT1 Xi 3.4 OS	13-422410	D
	Insertion tool medium	Ø 3.4 mm	IT2 Xi 3.4 OS	13-422411	D
112 Ar (d+ 4)	Insertion tool short		IT 3 Xi 3.4 OS	13-422412	D
H	Insertion tool long, 36 mm	For implant	IT1 Xi U OS	13-422413	D
	Insertion tool medium, 26 mm	Ø 3.8 - 5.5 mm (4.5 6+2 also)	IT2 Xi U OS	13-422414	D
	Insertion tool short, 18 mm		IT3 Xi U OS	13-422415	D



TOOLS

H-----

HT 1.77 M

- 100

HOLDBOALT HAS 2 CINDER-NO. 463100

----8

1

(FE

A

De	escription	Code	REF	Price cat.
	arallel depth gage brown	PDG 3.0	13-422416	В
	rallel depth gage silver	PDG 3.4	13-422417	В
Pa	arallel depth gage yellow	PDG 3.8	13-422418	В
Pa	arallel depth gage blue	PDG 4.5	13-422419	В
	andardized X-ray measuring probe. Scaling 1 mm, ray measuring for cylindrical implant types	PDG	13-425400	А
He	ex-instrument for 1.22 - 1.25 mm screws	HT 1.25	13-425100	С
He	ex-instrument for contra-angle	HT 1.25 M	13-425112	В
	ex-instrument, extralong, length 45 mm, 1.25 mm	HTX 1.25	13-425102	С
ΗT	TW 1.25	HTW 1.25	13-425111	В
	ex-instrument for 1.77 mm rews + Localicer®	HT 1.77	13-425103	С
He	ex-instrument for contra-angle	HT 1.77 M	13-425113	В
He	ex-instrument for contra-angle, extralong	HTX 1.77	13-425104	С
He	ex-instrument for 0.9 mm screws	HT 0.9	13-422428	С
He	ex-instrument for contra-angle	HT 0.9 M	13-425114	С
He	ex-instrument for contra-angle	TT 1.25 M	13-425115	С
	rill extension contra-angle (DX2), ctends by 19 mm	DX 2	13-500704	D
Fla	at spanner (for TB Xi ballhead abutments)	HAS2	13-463109	Н
То	ool E, 20 mm long	Tool E	13-462377	D
То	orque wrench, 10 - 70 Ncm	TW2	13-425402	S

Torque wrench, 10 - 70 Ncm

We recommend to have our torque wrenches calibrated annually.

GUIDE JACKET



Description	Amount	Material	REF	Price cat.
BFH 2.0 guide jacket for pilot drill 2.0mmd	Pack of 5	Ti6Al4V	13-425410	A
BFH 2.5 guide jacket for pilot drill 2.5mmd	Pack of 5	Ti6Al4V	13-425411	A
BFH 3.0 guide jacket for pilot drill 3.0mmd	Pack of 5	Ti6Al4V	13-425412	A
BFH 3.2 guide jacket for pilot drill 3.2mmd	Pack of 5	Ti6Al4V	13-425413	A
BFH 3.5 guide jacket for pilot drill 3.5mmd	Pack of 5	Ti6Al4V	13-425414	A



a) length	5 mm
b) hight of step	0.7 mm
c) max. Ø top	3.7 - 5 mm
d) nominal Ø	3 - 4.4 mm
e) Ø of drilling in the drill template	2.05 - 3.55 mm



Model with residual teeth for the fabrication of a drill guide for creating cavities for fixating the later drill guide for implant cavities.



Drill guide for creating cavities for later fixation of the surgical dill guide.



Surgical drill guide for safe BECES® placement. The drill sleeves are designed for 2.0 mm Twist drills.



	I EATLESS ® DRILI	_ DFN / DFLN	FOR IMPLANTS V	VITH CONICAL C	ORE		
Heat		Description	Ø working space	max. drill depth	Length	REF	Price cat.
	7 11 15 19 23	Pilot drill DS 2	2.0 mm	15 mm	32 mm	13-425001	D
	9 13 17 21	DS 2.8	2.8 mm	15 mm	32 mm	13-425005	D
		DSL 2	2.0 mm	23 mm	41 mm	13-425002	D
		Form drill					
		DFN 3.0	2.7 mm	15 mm	36 mm	13-425030	Е
		DFN 3.4	3.0 mm	15 mm	36 mm	13-425031	E
		DFN 3.7	3.4 mm	15 mm	36 mm	13-425032	Е
		DFN 4.2 - 4.5	4.05 mm	15 mm	36 mm	13-425033	Е
		Form drill long					
(<u></u>	-115	DFLN 3.0	2.7 mm	18 mm	39 mm	13-425035	E
and the second second	AND	DFLN 3.4	3.0 mm	18 mm	39 mm	13-425036	E
UTLN 2.1	A CONTRACT	DFLN 3.7	3.4 mm	18 mm	39 mm	13-425037	E
and the second second	110	DFLN 4.2 - 4.5	4.05 mm	18 mm	39 mm	13-425038	E
la sur		Corticalis drills (1 C Drill Xi 3.0	use alternatively to bon 3.0 mm	e tap) 5 mm	27 mm	13-425041	D
hime	Tar and	C Drill 3.4	3.2 mm	5 mm	27 mm	13-425042	D
No.	5 mm	C Drill 3.7	max. 3.8 mm	max. 5 mm	27 mm	13-425043	D
	12245	C Drill 4.2 - 4.5	max. 4.6 mm	max. 5 mm	27 mm	13-425044	D

The actual drill diameter is smaller than the implant diameter. * Namely for the implant systems Hexacone[®] and Xign[®].

IT HAS BEEN SCIENTIFICALLY PROVEN

Heatless[®] drills by Dr. Ihde Dental generate 55% less heat than traditional bone drills from other manufacturers. This makes it possible to use higher rotational speeds: between 3,000 and 5,000 rpm are recommended with good external cooling and intermittent drill technique.

CORTICAL TAP FOR 1ST CORTICALIS, WITHOUT HEX

	Description Cortical tap		Code KSG Xi 3.0 OS	REF 13-422424	Price cat. F
	Cortical tap	For the cap use in	KSG Xi 3.4 OS	13-422425	F
10.2 mm Cortic	Cortical tap	hard cortical bone only	KSG Xi 3.8 OS	13-422426	F
	Cortical tap		KSG Xi 4.5 OS	13-422427	F

DRILLSTOPP TRAY

Not suitable for dry heat sterilizers



Description	Code	REF	Price cat.
Drill stop A		13-500881	
Drill stop C		13-500883	
Drill stop D		13-500884	
Drill stop E		13-500885	
Drill stop G		13-500887	
Drill stop I		13-500889	
Drill stop J		13-500890	
Drill stop K		13-500891	
Drill stop L		13-500892	
Form drill	DFN 3.0	13-425030	
Form drill	DFN 3.4	13-425031	
Form drill	DFN 3.7	13-425032	
Form drill	DFN 4.2 - 4.5	13-425033	
Form drill	DFN 5.5	13-425034	
Form drill	DFLN 3.0	13-425035	
Form drill	DFLN 3.4	13-425036	
Form drill	DFLN 3.7	13-425037	
Form drill	DFLN 4.2 - 4.5	13-425038	
Tray with content		13-60031-K	739.00

STARTER TRAY

Autoclavable up to 134° C, not suitable for dry heat sterilizers This surgical kit contains all drills and tools for first works with the Xign[®] system. Material: autoclavable plastic



Description	Code	REF	Price cat.
Insertion tool medium	ITT2 Xi 3.0	13-422400	
Insertion tool medium	IT2 Xi 3.4	13-422403	
Insertion tool	IT2 Xi U	13-422406	
Hex-instrument	HT 0.9	13-422428	
Hex-instrument, long	HT 1.25	13-425100	
SpiralDrill DS 2	DS 2.0	13-425001	
Form drill DFN 3.0	DFN 3.0	13-425030	
Form drill DFN 3.4	DFN 3.4	13-425031	
Form drill DFN 3.7	DFN 3.7	13-425032	
Form drill DFN 4.2-4.5	DFN 4.2-4.5	13-425033	
Standardized X-ray measuring probe	PDG	13-425400	
Torque wrench	TW2	13-425402	
Tray with content*		13-S60044-K	upon request
*empty tray upon request			

Please read our detailed instructions for cleaning and re-sterilisation of surgical instruments on

 $https://implant.com/images/Instructions/GA/instrumente/Web_996-Anleitung_Resterilisation_instrumenten-03_V008.pdf$



IMPLANT SYSTEM FOR ENDOSSEOUS DENTAL IMPLANTIONS



Secure anti-rotation through high precision internal hexagon Apical expanded bone thread Excellent stability in all bone qualities: double condensation Universal application for fixed and removable prosthodontics Abutment alignment and 100% tightness through the taper

HC2 implants have an apical thread and feature an internal 6-edge, an internal marginal taper and a US standard internal thread.

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).


Due to technical reasons the product dimensions shown in this brochure might deviate from reality. **HC2** is a registered trademark. **HC2** implants are patent-protected.

In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.

PREPARATORY WORK FOR TEMPLATE APPLICATION

- 1. Ask your laboratory to prepare a drill template with the determined drill holes for the pilot drills. To be on the safe side, you can ask the laboratory to insert guide sleeves (**REF** BFH) into the drillholes, which specify the exact drill direction. Please use a 2.0 / 2.2 mm Ø drill for the pilot drilling.
- 2. For the following drill sequences you can use drill stops, which can be attached and tightened to the drill according to the length of drilling channel. Gingival thickness and template height are taken into account as needed. Thanks to the extremely high cutting efficiency of our drills, no ascending drilling sequences will usually be required.

Recommended RPM: 2000-5000. Apply sufficient cooling and allow the cooling to reach the working blades of the drills.




General note: HC2 implants are used as compression screws. In order to achieve a good bone condensation and implant stability, the drilling should be carried out thinner than the core diameter of the implant. The minimal diameter of the drill depends on the bone density. It is therefore not possible to advise drill-sequences which fit all bone-qualities. Typically in the soft maxillary bone only small drill-diameters are used (e.g. the usage of DOS1 only, for HC2 implants with 3.3 - 5.5 mm diameter), whereas in the highly mineralized lower jaw a specific drill sequence with respect to the mineralization of the bone is necessary. For insertion under pressure use the Handgrip. Due to technical reasons HC2 2.9 mmd is not available with expanded apical thread. HC2 implants with diameters 2.9 and 3.3 mm as well as 3.7 mm are not for use as single tooth restauration.



SURGERY

1. Recommended drill sequence



2. Implant packaging



Original packaging




Open the blister using the flap. Remove the label and stick it into the patients record.



The blister (secondary packaging) contains the implant in a sterile tube (primary package).

3. Remove the implant from its packaging

- 1. Open the lid. The implant is connected to the lid through a breakable section.
- 2. Remove the implant without touching the inner walls of the tube.



4. Handling

Attach the insertion tool to the implant by holding the top, to which the implant is secured, with your other hand.


Alternative: Firmly attach the assembled contra-angle handpiece instrument IT 2.5 M to the implant. For ratchets ITL 2.5 can be used as well.


After you have attached the insertion tool, firmly hold the lid in your hand and break the implant off the top along the break joint. Then insert the implant into the drill hole as much as possible.

Insertion tool IT 1 or IT 2 Make sure the hexagon is in the correct position and that the tool is fully inserted.

HC2 Implant

Break joint Lid with implant mount







5. Insertion


Using the ratchet screw the implant clockwise into the cavity.


The endosseous part of the implant must be completely covered by the bone.

After insertion the implant can be turned by a ¼ rotation backwards in order to relieve the bone and allow blood access to the implant site.

6. Remove insertion tool from implant

Remove the insertion tool from the implant.






7. Result

Result: A correctly inserted implant




8. Post-operative treatment

Close the implant with the suitable surgical cover screw.



After healing: Remove the surgical cover screw.






409. Handgrip



Break off the holder.

Insert the implant with axial pressure while turning.





Max. insertion torque for diameter

2.9 mm	30 Ncm
3.3 mm	40 Ncm
3.7 mm	50 Ncm
4.5/5.5 mm	60 Ncm



10. Pick Up Impressions


Impression taking with an A-silicone such as Safeprint® by Dr. Ihde Dental. The use of open and closed impression tray is possible.

10.1 Pick-up-procedure with an individual impression tray.

Hex tool HT 1.25

Tightening of the impression post HLT

HC2 Implant




Impression tray

Loosen screw with HT 1.25

Window in impression tray

Jaw



10.2 Prior to the impression

For pick up impressions the tray is inserted over the impression post until the screw peaks out on the other side and becomes accessible for the HEX-tool.

The impression post HLT must not necessarily be unscrewed from the implant in order to remove the impression tray. It can be repositioned later as well

10.3 Taking the impression

Disengage HLT from the implant: HLT remains in the impression

After the impression is taken, the implant is closed with a healing abutment (Gingiva former - straight or anatomic) and the impression is sent to the laboratory.

10.4 Preparation of the impression tray for model fabrication

Screw analog against the impression post.

Fasten the laboratory analog in the impression using HT 1.25

HLT

IA

HLT

HC2 Implant

11. Closed tray impressions

11.1 Impression taking with a closed impression tray


To take impression use an adequately large impression tray.

Impression posts TS/TSL are mounted with the help of the knurled screw

Tighten the impression post with the knurled screw

TS/TSL HC

HC2 Implant




11.2 Inserting impression

The filled impression tray is positioned sufficiently deep over the impression post to also allow an impression of the mucosa.

Impression post TS/TSL HC

HC2 Implant




11.3 Removing impression

When the closed-tray method is applied, the impression post TS/TSL HC remains on the implant after the impression tray is removed. After removal of the impression tray the impression post will be unscrewed and repositioned in the impression.

After the impression is taken, the implant is closed with a healing abutment (Gingiva former - straight or anatomic) and the impression is sent to the laboratory. TS/TSL HC

HC2 Implant



11.4 Mounting the lab analog

to casting.

Screw analog IA or IA HC M to the transfer post TS HC.

A
Tighten the impression post onto the
Afterwards the impression post is repositioned in the
impression.
B
The impression can now be casted. In IA HC M
block the lower access to the lock screw out prior
TS HC

IA or IA HC M

12. Laboratory procedures

12.1

The impression is poured. Then the impression posts (HLT or TS/TSL HC) are unscrewed from the laboratory analog.

12.2

The laboratory analog is now in the proper position and orientation in the Gypsum.

12.3

Positioning of the screwable abutments TLA15 HC, thereby the optimal position and adequate angulation must be determined.

NOTE The hexagon must be completely inserted into the analog.

HT 1.25

IA or IA HC M

Insert screw

TLA 15 Take care to position the hexagon correctly

IA or IA HC M

Arrow: View from above. Turn clockwise.

Gypsum

Bone

12.4

Ensure proper position of the abutment when transferring into the mouth.

Tightening torque of the screw during fastening on the implant: 20 Ncm

12.5


If more than one angled abutment is used, your laboratory will prepare a detachable synthetic bar (e.g. from Pattern Resin) in order to facilitate the correct positioning in the mouth.

TLA 15 HC

TLA 15 HC

Bone

Pattern Resin®



HC2 IMPLANTS

HC2 implants have a roughened endosseous surface and a machined apical thread. They feature an internal hex, an internal marginal taper and a US standard internal thread.



HC2 WITH AGGRESSIVE APICAL THREAD: HC2

As a result of many years of clinical observation of products, the design of the famous HC2 implant has been revised: the broadened apical thread is fully self-cutting. Thanks to the new apical thread portion, the implant is more stable even in weak bone and higher insertion torque can be reached.

If the implant is anchored in the 2nd cortical, it may be used in immediate load protocols. Especially in the upper jaw the usage of the new handgrip (REF 13-311431, with adapter IT HC REF 13-418196,) for inserting the implant is mandatory. This tool allows to apply vertical insertion forces and will enhance the anchorage. The drill sequence remains unchanged compared to the former design of the HC2 implant. And of course all abutments and tools remain the same.

Should the first cortical be unusually firm, the insertion can be achieved using the handgrip REF 13-311431 with the adapter IT HC REF13-418196.

Dimensions HC2 4.5 + 13	
a) basal thread Ø	4.3 mm
b) height of the apical thread height of polished collar	3.2 mm 2.0 mm
c) Implant neck height	2.5 mm
d) nominal Ø	3.7 mm

Application limitations

HC2/HC implant with a diameter of 2.9 mm - 3.3 mm implants may not be placed in a loaded area, especially not in the molar or premolar area. Likewise these implants may not be used where diagonal loading (off-axis loading) occurs, i.e. not for upper anteriors. Under no circumstances may Hexacone 2.9 mm implants be used for work that involves unsupported occlusal surfaces (consoles). If used in immediate load protocols, the prosthetic construction must be safely inserted on the 2nd post-operative day, and it should not be removed within the first 6 months.

In general we recommend to use implants up to (and including) the diameter 3.7 mm with care and not to use them for single tooth replacements, unless strict force control/force distribution is guaranteed.



HC2 IMPLANTS WIDE APICAL THREAD

Maximum insertion torque see page 70

5
1
2
T
D





Description	endosseous Ø	endosseous length	REF	Price cat.
HC2 3.3 8	3.3 mm	8 mm	13-412220	G
HC2 3.3 10	3.3 mm	10 mm	13-412221	G
HC2 3.3 11.5	3.3 mm	11.5 mm	13-412222	G
HC2 3.3 13	3.3 mm	13 mm	13-412223	G
HC2 3.3 15	3.3 mm	15 mm	13-412224	G
HC2 3.7 8	3.7 mm	8 mm	13-412202	G
HC2 3.7 10	3.7 mm	10 mm	13-412203	G
HC2 3.7 11.5	3.7 mm	11.5 mm	13-412210	G
HC2 3.7 13	3.7 mm	13 mm	13-412204	G
HC2 3.7 15	3.7 mm	15 mm	13-412205	G
HC2 4.5 8	4.5 mm	8 mm	13-412206	G
HC2 4.5 10	4.5 mm	10 mm	13-412207	G
HC2 4.5 11.5	4.5 mm	11.5 mm	13-412208	G
HC2 4.5 13	4.5 mm	13 mm	13-412209	G
HC2 5.5 8	5.5 mm	8 mm	13-412211	G
HC2 5.5 10	5.5 mm	10 mm	13-412212	G
HC2 5.5 11.5	5.5 mm	11.5 mm	13-412213	G
HC2 5.5 13	5.5 mm	13 mm	13-412214	G



HC2 6+2 IMPLANTS WITH AGGRESSIVE APICAL THREAD

HC2 6+2 was especially developed for the area of the 1st and 2nd molars in the upper and lower jaw. It is possible and recommendable to use it as a compression screw implant in the upper jaw. Endosseous length 6-8 mm. The upper edge of the polished 75° reverse cone can end at bone level or slightly above it. HC2 6+2 implants have a laser-generated surface structure (no-itis® laser) in the enossal area.

The conical polished implant head (c) should be submerged into the bone.

4.5 - 5.5 mm

6 - 8 mm

2 mm

,	75° conus
د]	
	L.
	b
a	

a) endosseous Ø:

c) reverse conus:

b) endosseous length:

Description	endosseous Ø	endosseous length	REF	Price cat.
HC2 4.5 6+2	4.5 mm	6-8 mm	13-412217	G
HC2 5.5 6+2	5.5 mm	6-8 mm	13-412218	G

Delivery includes surgical screw CSTI, REF 13-418101.



TRADITIONAL HC 6+2 IMPLANTS

HC 6+2 was especially developed for the area of the 1st and 2nd molars in the upper and lower jaw. It is possible and recommendable to use it as a compression screw implant in the upper jaw. Endosseous length 6-8 mm (8 mm incl. reverse cone). The upper edge of the polished 75° reverse cone can end at bone level or slightly above it. HC 6+2 implants have a laser-generated surface structure (No-Itis[®] laser) in the enossal area.


The conical polished implant head (a) should be submerged into the bone.

^{75°} conus		Description HC 4.5 6+2 HC 5.5 6+2	endosseous Ø 4.5 mm 5.5 mm	endosseous length 6-8 mm 6-8 mm	REF 13-413217 13-413218	Price cat. G
a	b					
a) endosseous Ø:	4.5 - 5.5 mm					
b) endosseous length:	6 - 8 mm					
c) reverse conus:	2 mm					



TRADITIONAL **HC** IMPLANTS

Maximum insertion torque see page 7



Dimensions HC Implants	
a) endosseous nominal Ø	2.9 - 5.5 mm
b) endosseous length	8 - 15 mm
c) length micro thread / polished neck	2.1 mm

	Description	endoss.Ø	endoss. length	REF	Price cat.
	HC 3.38	3.3 mm	8 mm	13-413220	G
	HC 3.3 10	3.3 mm	10 mm	13-413221	G
	HC 3.3 11.5	3.3 mm	11.5 mm	13-413222	G
	HC 3.3 13	3.3 mm	13 mm	13-413223	G
	HC 3.3 15	3.3 mm	15 mm	13-413224	G
	HC 3.78	3.7 mm	8 mm	13-413202	G
	HC 3.7 10	3.7 mm	10 mm	13-413203	G
	HC 3.7 11.5	3.7 mm	11.5 mm	13-413210	G
1E	HC 3.7 13	3.7 mm	13 mm	13-413204	G
	HC 3.7 15	3.7 mm	15 mm	13-413205	G
100	HC 4.1 8	4.1 mm	8 mm	13-413300	G
	HC 4.1 10	4.1 mm	10 mm	13-413301	G
	HC 4.1 11.5	4.1 mm	11.5 mm	13-413302	G
1	HC 4.1 13	4.1 mm	13 mm	13-413303	G
	HC 4.5 8	4.5 mm	8 mm	13-413206	G
	HC 4.5 10	4.5 mm	10 mm	13-413207	G
	HC 4.5 11.5	4.5 mm	11.5 mm	13-413208	G
EF	HC 4.5 13	4.5 mm	13 mm	13-413209	G
	HC 5.5 8	5.5 mm	8 mm	13-413211	G
	HC 5.5 10	5.5 mm	10 mm	13-413212	G
	HC 5.5 11.5	5.5 mm	11.5 mm	13-413213	G
EE	HC 5.5 13	5.5 mm	13 mm	13-413214	G



Delivery includes surgical screw CSTI, REF 13-418101.

HEXACONE® PLUS MU 0° IMPLANTS

g d e		Description	max. nominal Ø / without apical thread	max. nominal Ø / with apical thread	endosseous length	REF	Price cat.
cT T	Maximum insertion	HC Plus MU 3.3 13 0°	3.3 mm	4 mm	13 mm	13-412250	G
	torque: <mark>50 Ncm.</mark> Material Ti6Al4V	HC Plus MU 3.3 15 0°	3.3 mm	4 mm	15 mm	13-412251	G
		HC Plus MU 3.3 17 0°	3.3 mm	4 mm	17 mm	13-412252	G
		HC Plus MU 3.3 19 0°	3.3 mm	4 mm	19 mm	13-412253	G
		HC Plus MU 3.3 21 0°	3.3 mm	4 mm	21 mm	13-412254	G
E f	Dimensions	HC Plus MU 3.3 23 0°	3.3 mm	4 mm	23 mm	13-412255	G
	HC Plus MU 0° 4.1 17	HC Plus MU 4.1 10 0°	4.1 mm	4.7 mm	10 mm	13-412259	G
ŭ		HC Plus MU 4.1 13 0°	4.1 mm	4.7 mm	13 mm	13-412260	G
a) endosseous maximal Ø	3.3 / 4.1 mm	HC Plus MU 4.1 15 0°	4.1 mm	4.7 mm	15 mm	13-412261	G
b) endosseous length	11.5 - 21.5 mm	HC Plus MU 4.1 17 0°	4.1 mm	4.7 mm	17 mm	13-412262	G
c) length micro thread	1.5 mm	HC Plus MU 4.1 19 0°	4.1 mm	4.7 mm	19 mm	13-412263	G
d) platform Ø	4.8 mm	HC Plus MU 4.1 21 0°	4.1 mm	4.7 mm	21 mm	13-412264	G
e) height head	2.6 mm	HC Plus MU 4.1 23 0°	4.1 mm	4.7 mm	23 mm	13-412265	G
f) height of the apical threa	d 3.2						



Description	Code	REF	Price cat.
Insertion tool incl. screw REF 418316. For Hexacone Plus MU $0^\circ\!.$	IT HCMU	13-418315	F


HEXACONE® PLUS MU 15° IMPLANTS

d		Description	max. nominal Ø / without apical thread	max. nominal Ø / with apical thread	endosseous length	REF	Price cat.
g 🖉 🛛 e	Maximum insertion	HC Plus MU 3.3 13 15°	3.3 mm	4 mm	13 mm	13-412225	G
c T Material Ti6Al4V	HC Plus MU 3.3 15 15°	3.3 mm	4 mm	15 mm	13-412226	G	
- Le le		HC Plus MU 3.3 17 15°	3.3 mm	4 mm	17 mm	13-412227	G
		HC Plus MU 3.3 19 15°	3.3 mm	4 mm	19 mm	13-412228	G
		HC Plus MU 3.3 21 15°	3.3 mm	4 mm	21 mm	13-412229	G
		HC Plus MU 3.3 23 15°	3.3 mm	4 mm	23 mm	13-412230	G
E	Dimensions	HC Plus MU 4.1 10 15°	4.1 mm	4.7 mm	10 mm	13-412235	G
	HC Plus MU 15° 4.1 17	HC Plus MU 4.1 13 15°	4.1 mm	4.7 mm	13 mm	13-412236	G
		HC Plus MU 4.1 15 15°	4.1 mm	4.7 mm	15 mm	13-412237	G
a) endosseous maximal Ø	3.3 / 4.1 mm	HC Plus MU 4.1 17 15°	4.1 mm	4.7 mm	17 mm	13-412238	G
b) endosseous length	11.5 - 21.5 mm	HC Plus MU 4.1 19 15°	4.1 mm	4.7 mm	19 mm	13-412239	G
c) length micro thread	1.5 mm	HC Plus MU 4.1 21 15°	4.1 mm	4.7 mm	21 mm	13-412240	G
d) platform Ø	4.8 mm	HC Plus MU 4.1 23 15°	4.1 mm	4.7 mm	23 mm	13-412241	G
e) height head	3.9 mm						
f) height of the apical threa	d 3.2						

g) connecting part 2 mm

Description

11



Insertion tool for KOC MU, BECES MU & Hexacone Plus MU 15°. Use with IT2 BCS, IT2 S BCS, AH MU, handgrip. CodeREFPrice cat.ITX MU1513-418203G



ACCESSORIES FOR HEXACONE® PLUS MU

7

Description	Code	REF	Price cat.
Ratchet for all hex instruments and insertion tools.	RAT2	13-425051	К
Torque wrench 10 - 70 Ncm. It is recommended to have the torque ratchets	TW2	13-425402	S
recalibrated by us once a year.			

ACCESSORIES

	Description Insertion tool for KOS MU & BCS MU Use with IT2 BCS, IT2 S BCS, AH MU	Code ITX MU15	REF 13-418203	Price cat. G
	Titanium base Use with SF K MU (REF 418164)	T-Base MU	13-418188	В
	Castable abutment for use with T-Base and SF K MU	PA2 MU	13-418189	В
	Lab analogue for MU-implants	IA K MU	13-418159	A
	Prosthetic screw for KOC [®] MU and BECES [®] MU	SF K MU	13-418164	В
	Long screw for prosthetic use or as pick-up screw for use with TS MU (Tool: HT 1.25). Material Ti6Al4V	SFL MU	13-418168	В
	Castable abutment UCLA for direct use on MU-implants. SF K MU sold separately	PA MU	13-418119	В
	Transfer Coping (Temporary base) SF K MU must be ordered separately	TC MU	13-418161	D
	Transfer for pick-up, straight Delivery incl. SFL MU	HLT MU	13-418162	С
0	Adapter for handgrip Fits ITX MU15 (REF 13-418203)	AH MU	13-900041	F
	Scan abutment for MU implants, incl. screw SSA MU. Sterilisable, two-part Material Ti6Al4V	SAB MU	13-418205	D

SCANBODIES

	Description	System	Material	Amount	Code	REF	Price cat.
	Scanbody-MU cylyndrical	BECES® MU KOC® MU Hexacone® MU	POM	Pack of 5	Scanbody-MU	13-462056	В
Top view							
	Flag-Scanbody SCB MU incl. screw SFK MU (REF 13-418164) For intra-oral scan	BECES® MU KOC® MU Hexacone® MU	POM	Pack of 1	SCB MU	13-462073	В
Top view							

Please go to http://simpladent-implant.com/en/stl to download the corresponding STL files.

SURGICAL ACCESSORIES

Application limitations Hexacone[®] 2.9 mm implants may not be placed in a loaded area, especially not in the molar or premolar area. Likewise these implants may not be used where diagonal loading (off-axis loading) occurs, i.e. not for upper anteriors. Under no circumstances may Hexacone 2.9 mm implants be used for work that involves unsupported occlusal surfaces (consoles). If used in immediate load protocols, the prosthetic construction must be safely inserted on the 2nd postoperative day, and it should not be removed within the first 6 months.

In general we recommend to use implants up to (and including) the diameter 3.7 mm with care and avoid using them for single tooth replacements, unless strict force control is guaranteed.

			Description		Code	REF	Price cat.
				for 3 mm gingival height	HSI 3	13-418111	В
			Gingivaformer	for 5 mm gingival height	HSI 5	13-418112	В
				for 3 mm gingival height	HSIW 3	13-418191	В
Y		Wide gingivaformer	for 5 mm gingival height	HSIW 5	13-418192	В	
				3 mm height, 4.5 mm width	HSI 3-4.5	13-418268	В
Y	Y	Ų	Anatomic	3 mm height, 5.5 mm width	HSI 3-5.5	13-418269	В
		gingivaformer	5 mm height, 6.7 mm width	HSI 5-6,7	13-418270	В	
			Gingivaformer	3 mm height, 3.3 mm width	HSIS 3-3.3	13-418277	В

Screwable abutments for cemented bridges, without anti-rotation protection. Trimming and grinding is possible. Tighten with HT 1.25. Recommended insertion torque 20 Ncm.



Description	Code	REF	Price cat.
Height above implantat 8.5 mm The impression is made directly on the TCA, with tool TZ HC	TCA	13-418129	В
The impression is made directly on the TCA	TCA W	13-418173	В

Superstructure with hex and screw. Straight, for cemented bridges, without anti-rotation protection. Trimming and grinding is possible. Tighten with HT 1.25. Delivery inclusive screw SF 20. Recommended insertion torque 20 Ncm.



Description



Code



Description	Code	REF	Price cat.
Abutment	TLA W	13-418193	D
Anatomical abutment	ANAB	13-418276	Е



Description	C
15°angled, 1 mm gingival height	TI
15°angled, 2 mm gingival height	TI
15°angled, 3 mm gingival height	TI
25°angled, 1 mm gingival height	TI
25°angled, 2 mm gingival height	TI
25°angled, 3 mm gingival height	TI

Code	REF	Price cat.
TLA15 HC1	13-418135	F
TLA15 HC2	13-418136	F
TLA15 HC3	13-418137	F
TLA25 HC1	13-418139	F
TLA25 HC2	13-418140	F
TLA25 HC3	13-418141	F

REF

Price cat.



Delivery inclusive screw SF 20



IMPRESSION TAKING AND LABORATORY ACCESSORIES

And a

Į

0

Description Impression post Click-on No screw is needed	Code HLTC	REF 13-418107	Price cat. C
Short Impression post For TLA, TLA 15 and TLA 25 For Pick-up, with screw	HLT	13-418108	С
Pickup-Screw For HLT REF 418108	SF HLT long	13-418185	В
Impression post for HC Height 10.6 mm	TS HC	13-418109	С
Impression post for HC Height 15.5 mm	TSL HC	13-418110	С
Long impression post With screw	HLTS	13-418118	С
Lab analogue for Hex	ІА НС	13-418113	В

DIGITAL IMPRESSION TAKING



Description	Material	Unit	Code	REF	Price cat.
Scanbody for digital impression taking Screw SF 20 is optional and must be	POM	Pack of 5	Scanbody HC	13-418288	В
ordered separately					



Screwable spacer abutment for bridges and bars. Screw in with HT 1.77. Recommended insertion torque 25 Ncm.

Ψ	Description for gingival height 3 mm for gingival height 4 mm for gingival height 5 mm for gingival height 6 mm		Code TSA 3 TSA 4 TSA 5 TSA 6	REF 13-418143 13-418144 13-418145 13-418146	Price category B B B B
Description	TSA Analog	Castable abutment 10.5 mm high Pack of 5	Prosthetic screw for PSS on BTS/TSA		
Code	BTS	PSS (white)	SF		
REF	13-418152		13-418151		
Price cat.	В	В	В		

Screwable mesostructure for bridges and bars. Screw in with HT 1.77 hex key. Recommended insertion torque 25 Ncm. The position of the TCT hex is assigned with this approach.

	Description			Code	REF	Price category
	For gingival	height 0.5 mm		TCT HC 0	13-418130	В
	For gingival	height 1.5 mm		TCT HC 1	.5 13-418131	В
	For gingival	height 2.5 mm		TCT HC 2	2.5 13-418132	В
				A		
Description	Transfer post	Ung Screw	TCT analog	Castable abutment 12 mm high inside circular Pack of 5	Castable abutment 12 mm high inside edged Pack of 5	Fastening screw
Code	TST	SFL	BTT	PSTR (gray)	PSTA	SF
REF	13-418147	13-420428	13-418100	13-418124	13-418123	13-418151
Price cat.	В	А	В	В	В	В

This set contains all necessary components for the mesiostructure. For bridges and bars. Screwable (anti-rotation).

	Description Screw for PSTA Castable abutment, 12 m Mesiostructure for bridge		Code SF TCTL PSTA TCTL 0.5	REF 13-418165 13-418123 13-418138	Price category A B D
Hex				13-418263	F
			T		
Description	Laboratory analog for HC, HC2, with internal Hex	Long transfer post for HC, HC2, anti-rotation	Short transfer post for HC, HC2, anti-rotation	Castable al 12 mm hig Pack of 5	butment, round, h
Code	IA HC	HLTS	HLT	PSTR	
REF	13-418113	13-418118	13-418108	13-418124	
Price cat.	В	С	С	В	

ABUTMENTS

This abutment converts the internal hexagon of the HC2 implant into an external standard-hexagon. The prosthetic screw is screwed through. It tightens the prosthetic and the abutment at the same time.

	Description	Material	Code	REF	Price category
SF 275					
	Tempbase for HRA HC	PEEK	трв е	13-418274	С
	HEX reverse abutment incl. screw SF 275 REF 13-418275	Ti6Al4V	HRA HC	13-418273	D



LOCALICER°

Suitable tool: HT 1.77. We recommend a minimum of six implants per jaw and the use of a single denture as splint when using LOC abutments and KOC LOC.

Description		Hight	Code	REF	Price category
Localicer® fo	r Hexacone®	2 mm	LOC HC 2	13-418116	С
Localicer [®] fo	r Hexacone®	4 mm	LOC HC 4	13-418117	с

ACCESSORIES FOR LOCALICER®





Pull off force Yellow 600 g, Pink 1.200 g, Clear 1.800 g, Violet 2.700 g. Black has no retention and is designed for temporary solutions for up to one month.

MULTI-UNIT ABUTMENTS

Insertion of the angled MU2 abutments with HT 1.25; Insertion of the straight MU2S-Abutments with HT 1.77. Not for use on single implant constructions.

	Description	Material	Code	REF	Price cat.
h=3 mm	MU2 17 HC, angled, incl. SF 20	Ti6Al4V	MU2 17 HC	13-418281	L
1 h= 3 mm	MU2 35 HC, angled, incl. SF 20	Ti6Al4V	MU2 35 HC	13-418282	L
1.	MU2S 0.5 HC, straight	Ti6Al4V	MU2S 0.5 HC	13-418283	G
۲ 4	MU2S 1.5 HC, straight	Ti6Al4V	MU2S 1.5 HC	13-418284	G
V - 💎 -	MU2S 2.5 HC, straight	Ti6Al4V	MU2S 2.5 HC	13-418285	G
6 <u>mm</u>	GF MU2 Gingivaformer incl. SF MU2 Height above abutment shoulder 6 mm	Ti6Al4V	GF MU 2	13-418286	С
6,7 mm	MU2 Localicer incl. SF MU2 Height above abutment shoulder 6.7 mm	Ti6Al4V	MU 2	13-418287	С
	Prosthetic screw for MU2	Ti6Al4V	SF 20	13-420943	А

ACCESSORIES FOR MULTI-UNIT ABUTMENTS

Screw

rice cat.



ACCESSORIES

9	R	
-1	ŧ	

Description Ball abutment for fitting prostheses Application on TSA 3-6 abutments only Head diameter 2.5 mm	Code SB	REF 13-418153	Price cat. B
ΨΨΨΨ			



Description Hight Code REF Price cat. 0.5 mm TB 0.5 13-418126 В Frh Frh Frh Frh 0.5 mm TB 2 13-418127 В Ball abutment head - Ø 2.5 mm Screw in with HT 1.25 For use with NC - caps 0.5 mm TB 4 13-418128 В

Description	Pull-off force	Code	REF	Price cat.
	ca. 1200 g, transparent	NC	13-465028	A1
Nylon cap NC 2 pieces / pack	ca. 800 g, pink	NC 1	13-465029	A1
(EXTERNAL PRODUCT)	ca. 500 g, yellow	NC 2	13-465030	A1
Nylon cap R-NC	grün, strong	R-NC	13-465034	A1
with increased friction strength (for use with worn out Localicer®)	rosa, medium	R-NC 1	13-465033	A1
2 pieces / pack (EXTERNAL PRODUCT)	orange, soft	R-NC 2	13-465032	A1
. ,				_
Sleeve for all NC (EXTERNAL PRODUCT)		Н	13-465031	В

88

INSERTION TOOLS



Description

ITL 2.5

ITM 2.5

IT 2.5 M

ITWH 2.5 M

Туре	Code	REF	Price cat.
8 mm, click-on, hexagon (six edges)	IT 2.5	13-418174	В
22 mm, click-on, hexagon (six edges)	ITL 2.5	13-418175	В
20 mm, click-on, hexagon (six edges)	ITM 2.5	13-418176	В
Insertion tool for contra-angle handpiece	IT 2.5 M	13-418150	В
Insertion tool contra-angle for HC	ITWH 2.5 M	13-418184	с

TOOLS

	Description	Туре	Code	REF	Price cat.
	Hex tool 1.25	long, 21 mm	HT 1.25	13-425100	С
	Torx tool 1.25	long, for all screws, 21 mm	TT 1.25	13-425105	С
	Hex tool 1.25	for contra-angle handpiece, 45 mm	HTW 1.25	13-425111	В
	Hex tool 1.25	short, 14 mm	HTS 1.25	13-425101	С
	Hex tool 1.77	for all superstructures, 19 mm	HT 1.77	13-425103	С
M	Hex tool 1.25 M	for contra-angle handpiece, long, 26.1 mm	HT 1.25 M	13-425112	В
	Hex tool 1.77 M	for contra-angle handpiece, long, 28.6 mm	HT 1.77 M	13-425113	В
	Hex tool	extra long, 45 mm	HTX 1.25	13-425102	С
	Hex tool	for contra-angle, extra long, 45 mm	HTX 1.77	13-425104	С
	Punch	for contra-angle handpiece, 4.9 mm Ø	PUW1	13-425404	С
	Punch	manual, 5.2 mm Ø	PU	13-425406	С
	Standardized probe	Scale 1 mm for X-ray measurements 22 mm	PDG	13-425400	A
8	Drill extension contra-angle	extends by 19 mm	DX2	13-500704	D
_	Guide sleeve	for pilot drill, Titanium, 10 mm, 2.2 mm Ø Pack of 5	BFH	13-425401	A
660	X-ray measuring sphere	Surgical steel, 0.5 mm Ø Pack of 5	RM	13-425403	A




GUIDE JACKET





Description	Amount	Material	REF	Price cat.
BFH 2.0 guide jacket for pilot drill 2.0mmd	Pack of 5	Ti6Al4V	13-425410	A
BFH 2.5 guide jacket for pilot drill 2.5mmd	Pack of 5	Ti6Al4V	13-425411	A
BFH 3.0 guide jacket for pilot drill 3.0mmd	Pack of 5	Ti6Al4V	13-425412	A
BFH 3.2 guide jacket for pilot drill 3.2mmd	Pack of 5	Ti6Al4V	13-425413	A
BFH 3.5 guide jacket for pilot drill 3.5mmd	Pack of 5	Ti6Al4V	13-425414	A



a) length	5 mm
b) hight of step	0.7 mm
c) max. Ø top	3.7 - 5 mm
d) nominal Ø	3 - 4.4 mm
e) Ø of drilling in the drill template	2.05 - 3.55 mm



Model with residual teeth for the fabrication of a drill guide for creating cavities for fixating the later drill guide for implant cavities.



Drill guide for creating cavities for later fixation of the surgical dill guide.



Surgical drill guide for safe BECES® placement. The drill sleeves are designed for 2.0 mm Twist drills.



Description	Туре	Code	REF	Price cat.
Titan base incl. screw	Abutment basefor zirkonium. Anti-rotation (anti-rotation). Material Ti6Al4V	MB HC	13-418267	D

HANDGRIP-TRAY





THE THE MANAGE

Description	Туре	Code	REF	Price cat.
Adapter	for all contra angle handpiece instruments, compatible with Handgrip	Adapter Wst	13-310530	С
Ratchet	for all hex tools and insertion tools	RAT2	13-425051	K
Torque ratchet*	heavy duty for all Insertion-, hex- and torx tools 10 - 70 Ncm	TW 2	13-425402	S
Adapter for Handgrip		Adapter IT HC	13-418196	С
Handgrip **	selflocking	Handgriff A	13-311431	V

* We recommend having the torque ratchet recalibrated by us once a year.

** To clean this tool a heatable ultrasonic bath and a thermo disinfector (i.e. Miele TD-Serie) are required. If these devices are not available in the dental office the handle with REF 13-311431 should be purchased instead.



HANDGRIP-TRAY empty Size of closed tray: W 90 mm / L 195 / H 45 REF 13-60043 Price category G

For safe storage and sterilization of handgrips (max. 3 pieces) and adapters (max. 8 pieces). Plastic, autoclaveable up to 134° C, not suitable for dry heat sterilizers.



HEATLESS® DRILLS FOR IMPLANTS WITH CONICAL CORE

Surgical steel, color-coded, depth-coded and autoclaveable. The drill is marked with laser depth markings. Use between 3,000 and 5,000 rpm with good cooling and intermittent drill technique. Due to the extremely high cutting performance, you can work without pressure. For the implant systems HC2 and Xign[®]. Drill types DFN 3.0 - DFN 4.2-4.5.

	đ. li	and the standard	to to the set		C. I.	0.55	Distant
	Ø working range 0.1 - 1.5 mm	max. working depth	total length 31.7 mm	color code yellow	Code BCD 1	REF 13-900240	Price cat. C
n 00	0.1 1.5 mm	13 1111	51.7 mm	yenow		15 500240	C
	0.1 - 1.5 mm	15 mm	42 mm	yellow	BCDX 1	13-900243	C
100 0 0 10	2.0 / 3.6 mm	13 mm	30 mm		DFN 2.9 13	13-418102	E
	2.0 / 3.6 mm	15 mm	32 mm		DFN 2.9 15	13-418103	E
C7285	2.0 mm	17 mm	36.5 mm		DS 2	13-425001	D
	2.8 mm	17 mm	36.5 mm		DS 2.8	13-425005	D
H	4.5 mm	25 mm	44.5 mm		DSL+ 2.8	13-425015	E
	2.7 mm	18 mm	36 mm		DFN 3.0	13-425030	E
	3.0 mm	18 mm	36 mm		DFN 3.4	13-425031	Е
	3.4 mm	18 mm	36 mm		DFN 3.7	13-425032	Е
DFN 4.1	3.5 mm	18 mm	36 mm		DFN 4.1	13-425049	E
DFN 4,2-4,5	4.05 mm	18 mm	36 mm		DFN 4.2 - 4.5	13-425033	Е
Hanna	4.4 mm	18 mm	36 mm		DFN 5.5	13-425034	E
	2.7 mm	18 mm	39 mm		DFLN 3.0	13-425035	E
	3.0 mm	18 mm	39 mm		DFLN 3.4	13-425036	E
	3.4 mm	18 mm	39 mm		DFLN 3.7	13-425037	E
	4.05 mm	18 mm	39 mm		DFLN 4.2 - 4.5	13-425038	E
DFLN+SA	3 mm	25 mm	43.5 mm		DFLN+ 3.4	13-425029	E
No. of Concession, Name	3.4 mm	11.5 mm	30 mm		DFSN 3.7	13-425039	D
	3.9 mm	11.5 mm	30 mm		DFSN 4.2 - 4.5	13-425040	D
	max. 3.8 mm	max. 5 mm	27 mm		C Drill 3.7	13-425043	D
C Drill 4.1	max. 4.1 mm	2.5 mm	27 mm		C Drill 4.1	13-425050	D
CDnii424.5	max. 4.6 mm	max. 5 mm	27 mm		C Drill 4.2 - 4.5	13-425044	D
C DATE S. S. S. S.	max. 5.5 mm	2.5 mm	27 mm		C Drill 5.5	13-425045	D

IT HAS BEEN SCIENTIFICALLY PROVEN

that Heatless® Drills generate 55 % less heat compared to traditional bone drills by other manufacturers. This enables higher rotational speeds: We recommend between 3.000 and 5.000 RPM with good external cooling and intermittent drill technique.

TRAY

Autoclaveable up to 134° C. Not suitable for dry heat sterilizers. Size of closed tray: W 175 mm T 145 mm H 65 mm Please read our detailed instructions for cleaning and re-sterilization of surgical instruments on www.implant.com/en/downloads



Tray with content: Tray empty:

REF 13-S60017-K REF 13-60017-K

Description	Code	REF	Description	Code	REF
Twist drill	BCD 1	13-900240	Insertion tool short	IT 2.5	13-418174
Twist drill	DS 2	13-425001	Insertion tool medium	IT 2.5 M	13-418150
Twist drill	DS 2.8	13-425005	Universal adapter	UAW	13-425107
Form drill	DFN 2.9 13	13-418102	Hex tool 1.25 long	HT 1.25	13-425100
Form drill	DFN 2.9 15	13-418103	Hex tool 1.25 short	HTS 1.25	13-425101
Form drill	DFN 3.0	13-425030	Hex tool 1.77	HT 1.77	13-425103
Form drill	DFN 3.7	13-425032	Punch	PUW 1	13-425404
Form drill	DFN 4.2 - 4.5	13-425033	Drill extension	DX 2	13-500704
Form drill	DFN 5.5	13-425034	Standardized probe	PDG	13-425400
Form drill	DFSN 3.7	13-425039	Standardized probe	PDG	13-425400
Form drill	DFSN 4.2 - 4.5	13-425040	Standardized probe	PDG	13-425400
Cortical drill	C Drill 3.7	13-425043	Twist drill	DFLN 3.0	13-425035
Cortical drill	C Drill 4.2 - 4.5	13-425044	Twist drill	DFLN 3.7	13-425037
Cortical drill	C Drill 5.5	13-425045	Twist drill	DFLN 4.2 - 4.5	13-425038
Insertion tool long	ITL 2.5	13-418175	Torque ratchet	TW 2	13-425402



DRILLSTOP - TRAY

Not suitable for dry heat sterilizers.



Description	Code	REF	Price €
Drillstop A		13-500881	
Drillstop C		13-500883	
Drillstop D		13-500884	
Drillstop E		13-500885	
Drillstop G		13-500887	
Drillstop I		13-500889	
Drillstop J		13-500890	
Drillstop K		13-500891	
Drillstop L		13-500892	
Formdrill	DFN 3.0	13-425030	
Formdrill	DFN 3.4	13-425031	
Formdrill	DFN 3.7	13-425032	
Formdrill	DFN 4.1	13-425049	
Formdrill	DFN 4.2 - 4.5	13-425033	
Formdrill	DFN 5.5	13-425034	
Formdrill	DFLN 3.0	13-425035	
Formdrill	DFLN 3.4	13-425036	
Formdrill	DFLN 3.7	13-425037	
Formdrill	DFLN 4.2 - 4.5	13-425038	
Drillstop Tray with content		13-60031-K	739.00

STARTER TRAY

This surgical kit contains all drills and tools for first works with the HC2 system. Material: Plastic autoclaveable up to 134° C Not suitable for dry heat sterilizers.



Description	Code	REF	Price €
Insertion tool	IT 2.5	13-418174	
Insertion tool	ITL 2.5	13-418175	
Insertion tool	ITM 2.5	13-418176	
Hex tool long	HT 1.25	13-425100	
Twist drill	DS 2.0	13-425001	
Twist drill	DS 2.8	13-425005	
Formdrill	DFN 3.0	13-425030	
Formdrill	DFN 3.4	13-425031	
Formdrill	DFN 3.7	13-425032	
Formdrill	DFN 4.1	13-425049	
Formdrill	DFN 4.2-4.5	13-425033	
Corticalis drill 3.7	C-Drill 3.7	13-425043	
Corticalis drill 4.1	C-Drill 4.1	13-425050	
Corticalis drill 4.2 - 4.5	C-Drill 4.2 - 4.5	13-425044	
Torque ratchet	TW2	13-425402	

Starter Tray for HC2 with content
Starter Tray for HC2 empty

13-S60021-K 13-60021-K upon request upon request

(The products of this catalogue are CE marked (class I) and CE 1936 marked (class IIa and IIb) according to 93/42/EC Directive).

Commercial products that are not monitored by our notified body are declared as third-party products.

Basal implants may only be used and operated by qualified persons with valid authorisation (para. 2 MedProdAnw Verordnung).

We are certified according to DIN EN ISO 13485 and Annex II of Directive 93/42 EEC.

The product dimensions shown in this brochure may differ from reality for technical reasons.

If implants are reprocessed, there is a risk of the development of infections, because no validated method for processing exists. Implants therefore may not be reprocessed.

Compilation and explanation of symbols on the packaging:



Batch No.



Sterilized by



radiation



Non-sterile

Intended for use by dentists or surgeons only

Rx ONLY

Single use product

i



Expiry date

Instruction for use



REF

Store in a dry place

Store tightly keep closed

2

Do not use if packing is damaged

Do not

resterilize



Manufacturer

Production date

Catalogue number





Simpladent GmbH Dorfplatz 11 CH - 8737 Gommiswald / SG SWITZERLAND Phone +41 (0)55 293 23 70 implants@simpladent.ch www.simpladent-implants.com

PLEASE ASK YOUR LOCAL DISTRIBUTOR FOR THE VALID PRICE LIST

SIMPLADENT®



Dr. Ihde Dental GmbH Erfurter Str. 19 D-85386 Eching / München Tel +49 (0)89 319 761 0 info@ihde-dental.de