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DEUTSCH

DE

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ENGLISH

EN

Instructions for Use - Please read carefully

## Single part implant systems

### I. General instructions for dental implants

Familiarity with the relevant scientific and technical literature is a prerequisite for any implantological treatment. Any dentist undertaking such treatment must be familiar with the surgical procedure and experienced in evaluating normal and pathological treatment outcomes. Prosthodontic treatment planning must precede implantation – even though the surgical situation and outcome may necessitate subsequent adaptations to the treatment plan. The descriptions given in these Instructions for Use by themselves are insufficient to allow immediate use of the dental implant system. Have an experienced implantological surgeon give you an introduction to the system you selected. Seminars and trainings related to the implant system are offered by SIMPLADENT GmbH on a regular basis.

### II. Description

Single piece implant systems are systems for enossal dental implantation. They consist of the implants themselves, boneexpanding screws, impression abutments, a ratchet, a torque wrench and drills. Single implant systems are designed for a single-step dental implant approach.

### III. Main components

All implants are made of titanium alloy (Ti6Al4V, ASTM F 136-13, ISO 5832-3) or pure titanium. Analogues: Ti6Al4V, Aluminium  
Impression utilities: Ti6Al4V, Grade 4, ASTM F 136-13, ISO 5832-3  
Ratchet, torque wrench: Stainless surgical steel

### IV. Indications

Endosseous jaw-implant system.

### V. Contraindications

#### Absolute

- Patients whose systemic disease does not permit implantation (for a detailed understanding of this extensive contraindication, reference is made to studies and the relevant literature).
- Allergies to implant system components
- Single tooth implantations with implants that (at least partially) have a conical thread that is designed for lateral compression of bone along the enossal implant axis, provided that the nominal diameter of the implant is less than 3.8 mm.

#### Relative

- Diabetes, pregnancy
- Discomfort in the temporomandibular joint (although numerous temporomandibular joint complaints can be treated with implant therapy as the lower jaw position is stabilised with a fixed restoration)
- Root remains, centres of inflammation remaining in the jaw
- In children and young people who do not yet have a fully developed jaw

Little bone availability: If bone augmentation is necessary, and/or if there is a threat of immediate proximity to sensitive structures (nerves), it should be checked whether, for example, a corticobasal implant is less risky to use than a two-part implant. To avoid the occurrence of perimplantitis, an implant with a polished/machined implant surface or a one-piece implant should be given preference, if this possibility exists.

### VI. Diagnosis and patient information

Diagnosis consists of the patient's medical history, clinical examination and radiological examination using conventional dental radiographs and orthopantomographs and, if necessary, CBCT scans. It is recommended to have an experienced general practitioner perform a complete medical check-up. Implant procedures are elective procedures requiring comprehensive patient education. The patient must be familiarized with treatment alternatives and informed of possible consequences of the implant treatment. In particular, the patient must be informed of the importance of meticulous oral hygiene and of the cost of post-treatment. Please consult the pertinent but constantly changing laws, regulations and court decisions regarding the manner and extent of patient education required.

### VII. Adverse effects

#### Transient complaints

- Pain, swellings, difficulties while speaking, gingival inflammation.

#### More persistent complaints

- As with all other dental implants, a loss of bone substance in the horizontal and vertical dimension may occur even with a completely appropriate surgical and prosthodontic procedure. The nature and extent of the bone loss can not be predicted.
- If implants end up being located in the vicinity of certain special structures (nerves, maxillary sinus, adjacent teeth etc.), this may result in reversible or irreversible adverse effects on these structures.
- To prevent infection, antibiotic preventive and/or local disinfection measures are recommended.
- Chronic pain has been reported in connection with dental implants.
- In the event that osseointegration is unsuccessful, exfoliation (loss) of the implant may occur.

### VIII. Application

See system application brochure. Basal implants may only be handled and used by dentists/surgeons with valid authorisation issued by the manufacturer. This restriction includes all advices for patients, creating treatment plans, placement of implants, prosthetic treatments on implants, postoperative care (throughout the following years).

### IX. Interactions

There have been reports in the literature of interactions between titanium surfaces and oral care products containing fluoride. The patient must be informed of this. Because crestal implants, basal implants and natural teeth have different degrees of rigidity, restorations supported by any combination of several of these may have very different effects on the static balance of the system. Extensive implant/restoration systems with crestal implants may reduce the relative flexion and mobility of the cranial bones. There have been reports of pain in this context.

### X. Notes

- Following implantation, enter the type of implant used and its batch number in the patient documentation.
- Do not use implants beyond the expiration date given on the package.
- Store all implants in their closed packages in a dry place. Do not open the package until immediately prior to inserting the implant.
- Avoid all contact with extraneous substances before and during the insertion of the implant. Never touch the enossal part of the implant with your hand.
- If implants, abutments or accessories are swallowed accidentally, locate the affected object by appropriate means (e.g. a radiograph) and immediately take the requisite medical measures.
- The healing of the implant must be followed at regular implants, e.g. by radiological means. To avoid bone loss, non-osseointegrated or infected implants must be removed in good time. The time of implant removal will be determined by the dentist.
- The polished abutment heads must be roughened prior to the cementation of the prosthetic workpieces.
- In patients whose history shows periodontal diseases, smooth implants should be evaluated as an alternative to implants. A strict antibiotic regimen should also be evaluated.

### XI. Storage / Shelf life

The shelf life is 5 years from the date of sterilization. Store in a cool and dry place. Only implants in their original closed package and labelled "sterile" can be considered sterile. The manufacturer will not be liable for implants re-sterilized by the implantologist or third parties, regardless of the method used for sterilization. Sterile products are marked with the STERILE symbol. The expiration date is marked with the hourglass symbol. The LOT symbol indicates the batch or lot number. Implants are strictly intended for single use only.

SIMPLADENT GmbH reserves the right to change the design of its products and components or their packaging, to edit instructions or application notes or to re-negotiate prices or terms of delivery. Liability is limited to the replacement of the defective product. Any other remedy beyond replacement is excluded. Applications of our products may be described orally, in writing, through electronic media or by way of practical demonstrations. The user has the duty to personally determine whether or not any product is suitable for a particular purpose, indication or procedure.

The application of this product is beyond the control of SIMPLADENT GmbH and is therefore subject to the user's own responsibility. We warrant the quality of our implant systems within the framework of our Terms and Conditions for Sales and Deliveries.

996-1TEILIG-03

Revision issued 2021-07

In case that implants would be reprocessed (cleaned, restilizated) infections could occur, because no validated procedures for reprocessing are available.

When processing re-sterilizable medical devices, the rules of the RKI-Guideline and the European standard EN ISO 17664 are relevant and must be obeyed. Guidelines regarding the processing of our instruments are available in the brochure accompanying our instruments. At the same time these guidelines are displayed on the website [www.simpladent-implants.com](http://www.simpladent-implants.com).

Keep this and all medicines out of the reach of children. To be used by dentists or surgeons only.

*It is necessary that you carefully read the system application prospectus (SAP) for the implant system which you are planning to use. This prospectus contains important information about the application of the systems components.*

**BECES\* implants are only for sale or use to the authorized treatment provider.**

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Dr. Ihde Dental GmbH
Simpladent Implant Solutions Pvt Ltd