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Vertrieb

Simpladent Implant Solutions Pvt Ltd
Unit no 417, 4th Floor, KM Trade Tower,
H-3, Sector 14, kaushambi, Ghaziabad
201010 Uttar Pradesh
India

EC REP
Dr. Ihde Dental GmbH
Erfurter Str. 19
D - 85386 Eching / München
Tel +49 (0)89 319 761 0
info@ihde-dental.de

Two-part implant systems

General instructions for dental implants

Knowledge of the extensive specialized literature is prerequisite for the performance of implantological treatments. Surgical experience in the evaluation of the normal and pathological treatment is mandatory. The prosthetic planning must be performed prior to implantation - also if the surgical situation may cause a later alteration in the planning. These instructions alone are not sufficient for the immediate use of the Implant System.

Description

The Implant System is a system for enossal dental implantations. It consists of implants, surgical and prosthetic instruments as well as abutments. The Implant System is suitable for single and two step implantation methods.

Constituents

The materials constituent the implant can be found on the respective package label or in the brochure concerning the use of the system.

Indications

Endosseous jaw-implant system.

Contraindications

Absolute

- Patients whose systemic disease does not permit implantation (for a detailed understanding of this extensive contraindication, reference is made to studies and the relevant literature).
- Allergies to implant system components
- Single tooth implantations with implants that (at least partially) have a conical thread that is designed for lateral compression of bone along the enossal implant axis, provided that the nominal diameter of the implant is less than 3.8 mm.

Relative

- Diabetes, pregnancy
- Discomfort in the temporomandibular joint (although numerous temporomandibular joint complaints can be treated with implant therapy as the lower jaw position is stabilised with a fixed restoration)
- Root remains, centres of inflammation remaining in the jaw
- In children and young people who do not yet have a fully developed jaw

Little bone availability: If bone augmentation is necessary, and/or if there is a threat of immediate proximity to sensitive structures (nerves), it should be checked whether, for example, a corticobasal implant is less risky to use than a two-part implant. To avoid the occurrence of periimplantitis, an implant with a polished/machined implant surface or a one-piece implant should be given preference, if this possibility exists.

Use

Diagnosis / Explanation

The diagnosis consists of the anamnesis, the clinical examination, the radiological examination with small format x-ray and orthopantomogram as well as CT examinations. It is recommended to have a medical check-up performed by an experienced medical doctor.

Implantations are optional interventions, which require extensive explanation to the patient. Treatment alternatives and results of the implantations as well as results of a possible failure must be explained. In particular, the importance of oral hygiene and the costs of a following treatment must be explained. For type and extent of the explanation, reference is made to the respective case law.

Approach

Under local anesthesia, preferably under flap formation, an implant bed is established with rotating or displacing instruments. For preparation of an adequate bone bed, the use of the drill and the compliance with a drilling sequence is recommended:

A corresponding additional distance to endangered structures must be complied with. For a detailed scheme concerning this aspect, please refer to our elaborate product brochures.

During the preparation of the bone cavity, sufficient cooling of the implant bed is to be ensured (i.e. with cooled (sterile) physiological saline solution).

Never drill above 500-800 rpm. Always use a sharp drill (max. 10 applications). Use intermittent reaming and drilling technique.

The implant is removed from the sterile packaging immediately prior to introduction and stably inserted in the bone bed. It is mandatory that a good primary stability is achieved. Screw implants are screwed in with the help of a ratchet (RAT2), the torque wrench (TW2) or a screwdriver.

To minimize frictional warmth, it must be screwed in slowly and without great pressure. Depending on the available bone and the intended planning, the implants are inserted with the neck portion (polished shaft) extending from the bone or completely covered by the bone. If a removable dental prosthesis will be worn over the implants covered with mucous membrane during the healing phase, the implant must be inserted in the bone in order to avoid exposure dur-

ing the healing phase as well to prevent chewing pressure on the implant. In the selection of the implant, the height and width as well as the quality of the present bone must be taken into consideration. It is recommended not to remove too much soft bone during the drilling, but rather to compress the bone to the side with displacing osteotomes.

Afterwards, the implant is either closed with the surgical screw on the mouth side (two step method), or a healing screw in adequate length (mucosa height) is screwed in (single step method). If indicated, an adequate wound closure and radiological control follow. During the healing phase (3-6 months are recommended depending on bone quality), the surgical situs is to be controlled in regular periods. In incomplete healing of the mucous membrane and / or premature mucous membrane penetration, it is recommended to replace the surgical screw with a healing cap. As long as the static of the total construction as well as the situation during the insertion allows, implants can also be loaded directly after the insertion. This is decided by the treating dentist. In case of doubt, a healing period of 3-5 months without loading should be complied with - depending on the quality of the bone. After completion of the healing phase, a radiological and clinical control is performed.

Abutments from the Implant System are available for the prosthetic treatment. The table below provides information concerning the fitting parts. For further information please refer to the respective brochure concerning the use of the system. Basal implants may only be handled and used by dentists / surgeons with valid authorisation issued by the manufacturer. This restriction includes all advices for patients, creating treatment plans, placement of implants, prosthetic treatments on implants, postoperative care (throughout the following years).

Side effects

Temporary discomforts

Pain, feeling of tension, swelling, difficulties speaking and inflammation of the gingiva.

Long term disorders

- Oroantral and oronasal fistula
- Localized or systemic infections
- Crater shaped bone degeneration over the duration of use
- A horizontal and vertical bone degeneration is to be expected, also in proper surgical and prosthetic proceedings - as in all other dental implants. Type and extent of the bone loss is not foreseeable. Functionally caused and inflammatory bone degeneration can accumulate. If implants are placed in the vicinity of certain structures (nerves, maxillary sinus, neighboring teeth etc), a reversible or irreversible impairment (damage) of these structures can occur
- To prevent infection, antibiotic preventive measures and/or local disinfection are recommended

Interactions

The interactions between titanium surfaces and oral hygiene products containing fluoride has been described in the literature. The patient is to be informed thereof. Due to the different rigidity of crestal implants and natural teeth, a greatly differing influence of the static is to be expected in combined works.

Extensive implant prosthetic systems on crestal implants can impair the flexion and mobility of the cranial bones against one another. Pain conditions have been described in this connection.

Note

- After the implantation, the type of the used implant and the charge number must be noted in the patient file in writing.
- Implants may only be used in the period prior to expiration.
- Implants must be stored dry and sealed. The seal is to be opened immediately prior to introduction. All contact with foreign substances prior to insertion of the implant is to be avoided. The enossal part of the implant may not be touched.
- Upon accidental swallowing of implants, abutments or accessories, the remainder of the object is to be secured (x-ray or similar) and the necessary medical measures taken.
- The healing of the implant must be followed at regular intervals, e.g. by radiological means. To avoid bone loss, non-osseointegrated or infected implants must be removed in good time. The time of implant removal will be determined by the dentist.
- The polished abutment heads must be roughened prior to the cementation of the prosthetic workpieces.
- In patients whose history shows periodontal diseases, smooth implants should be evaluated as an alternative to implants. A strict antibiotic regimen should also be evaluated.

Shelf life / Storage of the implants

The shelf life is 5 years after sterilization; see expiration date on the package. Store tightly closed and dry.

Comment

Only originally sealed implants are sterile during the shelf life period. If implants are re-sterilized by the final consumer, no responsibility is assumed - independent of the method of sterilization. Sterile products are marked with STERILE. The expiration date is noted with the hour glass symbol. The LOT states

Legend

Read instructions

Expiration date


(Gamma) -sterilized

Only use once

Do not re-sterilize

This product is only for sale to dentists, oral surgeons and maxillo-facial surgeons or on their behalf.

Charge number

Store in dry place

Store tightly keep closed

Do not use if packing is damaged


Manufacturer

Catalogue number

-5° - 55°

Authorised representative in the European Community

CE 1936

the charge number. Implants may only be used once.

SIMPLADENT GmbH reserves the right to change the design of the product and components or the packaging thereof, to adapt instructions of use as well as to newly agree on prices or delivery conditions. The liability is limited to the replacement of the defective product. Further claims of all types are excluded.

996-S-2TEILIG-03

Revision issued 2021-07

In case that implants would be reprocessed (cleaned, re-sterilized) infections could occur, because no validated procedures for reprocessing are available.

When processing re-sterilizable medical devices, the rules of the RKI-Guideline and the European standard EN ISO 17664 are relevant and must be obeyed. Guidelines regarding the processing of our instruments are available in the brochure accompanying our instruments. At the same time these guidelines are displayed on the website www.simpladent-implant.com.

Keep this and all medicines out of the reach of children. To be used by dentists or surgeons only.

It is necessary that you carefully read the system application prospectus (SAP) for the implant system which you are planning to use. This prospectus contains important information about the application of the systems components.

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Dorfplatz 11
CH-8737 Gommiswald
SWITZERLAND
Tel. +41 (0)55 293 23 70
Fax +41 (0)55 293 23 00
simpladent@implant.com
www.simpladent-implant.com

Distributed by

Simpladent Implant Solutions Pvt Ltd

Unit no 417, 4th Floor, KM Trade Tower,
H-3, Sector 14, Kaushambi, Ghaziabad
201010 Uttar Pradesh
India

Dr. Ihde Dental GmbH
Erfurter Str. 19
D - 85386 Eching / München
Tel +49 (0)89 319 761 0
info@ihde-dental.de

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